

# Iron Metabolism Hepcidin 25 ELISA




**DRG**

# New Hepcidin 25 bioactive ELISA

## The gold standard in Hepcidin measurement

### Benefits

- Easy and straight forward assay procedure (no extraction or centrifugation)
- Total assay time < 2 hours
- All reagents ready to use
- Two controls included in the kit
- High sensitivity
- Good correlation to SELDI-TOF-MS



### Clinical relevance

Hepcidin deficiency can result in hereditary hemochromatosis (body iron overload) which can progress to liver fibrosis and cirrhosis. In contrast, defects in the HFE gene reduce hepcidin synthesis resulting in extremely low transferrin saturation and low to normal ferritin concentrations. Furthermore, infections and inflammation lead to increased hepcidin concentrations, causing iron deficiency that ultimately prevents growth of extracellular growing microorganisms.

### Assay characteristics

- Assay Principle: Competitive ELISA
- Dynamic Range: 0.354 - 80 ng/mL of Hepcidin
- Total Assay Time: approx. 2 hours (60/30/20 min.)
- Sample Volume: 20 µl of Serum or Plasma (EDTA, Citrate, Heparin)
- Mean Intra Assay Precision: 5.12 %
- Mean Inter Assay Precision: 10.33 %

### Ordering informations

Description	Code	Size
Hepcidin 25 bioactive ELISA	EIA-5258	96 Wells

### Intended use


The **DRG Hepcidin 25 bioactive ELISA** is an enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of Hepcidin in serum and plasma.

### Background

Hepcidin is a 25-amino acid, cysteine-rich peptide hormone, produced by the liver. Hepcidin controls plasma iron levels by acting on the Fe-transporter ferroportin. This reduces the absorption of iron from the intestine and the release of iron in the macrophage and hepatocyte. Hepcidin is secreted in response to iron overload and inflammation, while its concentration decreases during iron depletion.

### Method comparison

DRG Hepcidin 25 bioactive ELISA EIA-5258 showed good correlation to SELDI-TOF-MS (Surface-enhanced laser desorption/ionization -time-of-flight- mass spectrometry) ( $r=0.89$ )




# Hepcidin 25 bioactive ELISA

## Principle of the test

The DRG Hepcidin bioactive ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding. The microtiter wells are coated with a polyclonal antibody directed towards an antigenic site of the Hepcidin molecule. Endogenous Hepcidin of a patient sample competes with a biotinylated Hepcidin conjugate for binding to the coated antibody. After incubation, the unbound conjugate is washed off. Bound biotinylated Hepcidin is detected by streptavidin-horseradish peroxidase complex. After addition of the substrate solution, the intensity of colour developed is inversely proportional to the concentration of Hepcidin in the patient sample.

## Example of a typical standard curve

Standard	Optical Units (450 nm)
Standard 0 (0 ng/mL)	2.13
Standard 1 (2 ng/mL)	1.84
Standard 2 (6.5 ng/mL)	1.42
Standard 3 (25 ng/mL)	0.65
Standard 4 (45 ng/mL)	0.41
Standard 5 (80 ng/mL)	0.25



## Reproducibility

Sample	n	Mean (ng/mL)	CV (%)
1	20	66.4	3.3
2	20	63.5	9.9
3	20	10.1	2.1

## Inter-Assay

The between assay variability is shown below:

Sample	n	Mean (ng/mL)	CV (%)
1	40	4.9	11.5
2	40	23.8	12.0
3	40	63.1	7.5

## Linearity

	Sample 1	Sample 2	Sample 3
Concentration (ng/mL)	18.0	67.3	75.4
Average Recovery	99.1	96.3	96.8
Range of Recovery (%)	from	93.3	93.3
	to	111.1	99.9

## Recovery

	Sample 1	Sample 2	Sample 3
Concentration	11.0	68.0	81.0
Average Recovery	103.2	98.9	100.6
Range of Recovery (%)	from	101.6	93.2
	to	104.8	105.5

## Sensitivity

The analytical sensitivity of the DRG Hepcidin 25 bioactive ELISA is 0.35 ng/mL.

## Dynamic range

The dynamic range of the DRG Hepcidin 25 bioactive ELISA is between 0.35 - 80 ng/mL.

## Specificity

Hepcidin 100 %  
Pro-Hepcidin 0.04 %

# DRG ELISAS

## Oncology

CYFRA 21-1  
CA 72-4  
CA 15-3  
CA 125  
CA 19-9  
CEA  
TPS  
TPA  
PSA  
free PSA  
NSE  
Chromogranin

## Gyn. Endocrinology

Estradiol  
Progesterone  
17a-OH Progesterone  
DHEA-S  
Testosterone  
DHEA  
Estrone  
Androstendione  
DHT  
SHBG  
DHEA  
LH, FSH, PRL

## Prenatal Supervision

PAPP-A  
Free  $\beta$  HCG  
AFP  
Free Estriol  
HCG  
HPL  
PLGF

## Saliva Diagnostics


Cortisol  
Estradiol  
Testosterone  
DHEA  
Progesterone  
17a-OH Progesterone

## Diabetes/Obesity

Insulin  
C-Peptid  
Proinsulin  
Leptin

## Iron Metabolism

Hepcidin  
Pro-Hepcidin



## Bone Metabolism

25-OH Vitamin D Total

## Hypertension

Renin  
Aldosterone

## ELISAS that perform

DRG develops and manufactures ELISAS for use in clinical and research laboratories.

The experience of our production and management team guarantees to provide high quality products, competitive prices and excellent customer service.

DRG works to DIN EN ISO 9001:2000, ISO 13485:2003 and ISO 13485:2003 under CMDCAS standard, certified by TÜV Rheinland Product Safety GmbH, an indication of our commitment to customer service, quality control and improved health care.

## DRG Diagnostics

DRG Instruments GmbH, founded in 1973 by Dr. Geacintov, subsidiary of DRG Intl. Inc., USA, is a diagnostics manufacturer of ELISAS.

The DRG Group operates through a network of DRG subsidiaries in Germany, Poland, Russia and the Czech Republic and through distributors worldwide.



## DRG Subsidiaries

**Czech Republic**  
DRG spol. s.r.o.,  
Brno

**China**  
DRG Intl. China,  
Beijing

**Poland**  
DRG MedTek South Poland,  
Gliwice

DRG MedTek Sp. zo.o.,  
Warszawa

**Russia**  
DRG TechSystems,  
Moscow

DRG Biomed 0.0.0.,  
St. Petersburg



DRG Instruments GmbH, Germany  
Frauenbergstraße 18  
D-35039 Marburg  
Tel. +49 (0) 64 21/17 00 0,  
Fax +49 (0) 64 21/17 00 50  
Internet: www.drg-diagnostics.de  
E-mail: drg@drg-diagnostics.de

Distributed by



DRG International Inc. USA  
1167 U.S. Highway 22 East  
Mountainside, N.J. 07092 USA  
Phone: +1 (908) 233-2079  
Fax +1 (908) 233-0758  
Internet: www.drg-international.com  
E-mail: corp@drg-international.com

New

New