User's Manual

Free PSA ELISA

Enzyme Immunoassay for the quantitative determination of free Prostate Specific Antigen (PSA) in human serum or plasma

EIA-4189

96

DRG 

DRG Instruments GmbH, Germany

Frauenbergstr. 18, D-35039 Marburg

Telefon: +49 (0)6421-1700 0, Fax: +49-(0)6421-1700 50

Internet: www.drg-diagnostics.de

E-Mail: drg@drg-diagnostics.de

Version 9.0
Effective, November 2012

**Please use only the valid version of the package insert provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Arbeitsanleitung.
Si prega di usare la versione valida dell'inserto del pacco a disposizione con il kit.
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos

1	INTENDED USE	2
2	INTRODUCTION.....	2
3	PRINCIPLE OF THE TEST.....	2
4	MATERIALS PROVIDED WITH THE KIT.....	3
5	MATERIAL REQUIRED BUT NOT PROVIDED	3
6	STORAGE AND STABILITY.....	3
7	PRECAUTIONS	4
8	GUIDELINE FOR SAMPLE COLLECTION; PREPARATION AND STORAGE.....	4
9	ASSAY PROCEDURE	5
10	VALIDITY OF THE ASSAY	6
11	QUALITY CONTROL	6
12	EXPECTED VALUES LIMITATIONS OF THE PROCEDURE	7
13	PERFORMANCE CHARACTERISTICS.....	7
14	LEGAL ASPECTS.....	10
15	SUGGESTED READING	10
1	VERWENDUNGSZWECK	11
2	EINFÜHRUNG	11
3	TESTPRINZIP	11
4	IM TEST KIT ENTHALTENE MATERIALIEN	12
5	BENÖTIGTE NICHT IM TESTUMFANG ENTHALTENE MATERIALIEN	12
6	LAGERUNG UND STABILITÄT	12
7	VORSICHTSMASSNAHMEN	13
8	ORIENTIERUNGSHILFE FÜR PROBENNAHME; AUFBEREITUNG UND LAGERUNG	14
9	TESTDURCHFÜHRUNG.....	15
10	GÜLTIGKEIT DES ASSAYS.....	16
11	QUALITÄTSKONTROLLE	16
12	ERWARTETE WERTE UND GRENZEN DER METHODE	17
13	LEISTUNGSMERKMALE	17
14	RECHTLICHE GRUNDLAGEN.....	20
15	EMPFOHLENE LITERATUR	20
	SYMBOLS USED WITH DRG ASSAYS.....	21

1 INTENDED USE

The free PSA ELISA is used for the quantitative determination of free prostate specific antigen (f-PSA) in human serum or plasma samples. The determination of f-PSA levels is generally used in conjunction with a total PSA (t-PSA) measurement to determine the ratio between f-PSA and t-PSA. This ratio helps to estimate the risk for prostate cancer and to discriminate between elevated t-PSA levels caused by cancerous or non-cancerous conditions. F-PSA determinations are especially recommended for men with elevated t-PSA levels and negative results with digital rectal examination (DRE) in order to decide if a second prostate biopsy is indicated.

2 INTRODUCTION

Prostate cancer is the most frequent type of cancer found in man and is the second cause of death due to cancer in males. Until recently, digital rectal examination (DRE) was frequently used as only diagnostic modality for the detection of early stages of prostate cancer. In the recent years the determination of serum PSA levels has become the most accepted method to improve the diagnostic specificity of DRE. Although PSA is a tissue specific protein and is not solely tumor specific, it has become the most important marker for prostate carcinoma, showing a better specificity than other biochemical markers used in this context (PAP, total alkaline phosphatase, carcinoembryonic antigen, etc.)

In 1979, Wang et al. isolated a specific antigen for normal prostate tissue and called this protein PSA. PSA is a 33 kDa serine proteinase. Immunohistological studies have shown that PSA is localized in the cytoplasm of prostate acinar cells, ductal epithelium and in the secretion on the ductal lumina, present in normal, benign hyperplastic and malignant prostate tissues as well metastatic prostate cancer and in seminal plasma. If the structural integrity of the prostate is disturbed and/or the gland size is increased, the amount of PSA in the blood plasma may become elevated. In the blood plasma, most of the PSA forms complexes with various proteinase inhibitors. Only a small fraction of PSA circulates as free inactive PSA. Basically three major forms of PSA can be distinguished, only two of which are immunoreactive. The predominant form of PSA is a complex with α 1-antichymotrypsin (ACT-PSA). Inactive free PSA (f-PSA) represents around 10-40% of the immunologically detectable PSA. The total amount of immunoreactive PSA is known as total PSA (t-PSA). PSA complexed with α -2-macroglobulin cannot be detected by immunological assays and is therefore frequently called occult PSA (o-PSA).

Current methods of screening men for prostate cancer utilize the detection of t-PSA. Levels of 4.0 ng/ml or higher are strong indicators of the possibility of prostatic cancer and are an indication for follow-up examinations of the patient. However, elevated serum PSA levels are frequently also attributed to benign prostatic hyperplasia, leading to a high percentage of false positive screening results. A potential solution to this problem involves the determination of free PSA levels. Studies have suggested that the percentage of free PSA is lower in patients with prostate cancer than those with benign prostatic hyperplasia. Thus, the measurement of free serum PSA in conjunction with total PSA, can improve specificity of prostate cancer screening in selected men with elevated total serum PSA levels, which would subsequently reduce unnecessary prostate biopsies with minimal effects on cancer detection rates.

3 PRINCIPLE OF THE TEST

This f-PSA ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA). The microtiter wells are coated with an anti-f-PSA monoclonal antibody, directed towards an epitope of an antigen molecule. An aliquot of patient serum is incubated in the coated well with enzyme conjugated second antibody (E-Ab), directed towards a different region of the antigen molecule. After incubation the unbound E-Ab is washed off. The amount of bound E-Ab is proportional to the concentration of antigen in the sample. After adding the substrate solution, the intensity of colour developed is proportional to the antigen concentration in the sample. The measured ODs of the standards are used to construct a calibration curve against which the unknown samples are calculated.

4 MATERIALS PROVIDED WITH THE KIT

Each kit contains reagents sufficient for 96 determinations.

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells.
Wells coated with anti-free PSA antibody (monoclonal).
2. **Zero Standard/Sample Diluent**, 1 vial, 10.0 mL, ready to use.
Contains non-mercury preservative.
3. **Standard (Standard 1-5)**, 5 vials, 0.5 mL, ready to use;
Concentrations: 0.75 – 1.5 – 3.0 – 6.0 – 12.0 ng/mL
The standards are calibrated against NIBSC (WHO) Standard 96/668
Contain non-mercury preservative.
4. **Control Low & High**, 2 vial, 0,5 mL each, ready to use;
For control values and ranges please refer to vial label or QC-Datasheet.
Contain non-mercury preservative.
5. **Assay Reagent**, 1 vial, 6 mL, ready to use;
Contains non-mercury preservative.
6. **Enzyme Conjugate**, 1 vial, 6 mL, ready to use,
Anti-PSA antibody conjugated to horseradish peroxidase;
Contains non-mercury preservative.
7. **Substrate Solution**, 1 vial, 12 mL, ready to use,
Tetramethylbenzidine (TMB).
8. **Stop Solution**, 1 vial, 14 mL, ready to use,
contains 0.5 M H₂SO₄,
Avoid contact with the stop solution. It may cause skin irritations and burns.
9. **Wash Solution**, 1 vial, 25 mL (40X concentrated),
Dilute with distilled water before use.
see „Reagent Preparation“.

5 MATERIAL REQUIRED BUT NOT PROVIDED

- Precision micropipettes (volume: 25 µL and 100 µL) with disposable tips
- Distilled water
- ELISA photometer with 450 nm- and 630 nm-filters
- Timer with 60 min. range or higher
- Microplate washer (optional)
- Vortex or similar mixing tools
- Container for the proper handling of waste and samples after use

6 STORAGE AND STABILITY

- Store the kit and components at 2 °C to 8 °C
- Bring to room temperature (18 °C - 25 °C) at least 30 minutes before use. After use put back into the refrigerator. Avoid long time storage at room temperature.
- Do not use the kit or components after the expiry date. For expiry date of the original packed kit see kit label.
- Close the bottles immediately after use.
- Store the plate incl. desiccant in the provided zip-lock pouch. Modules that are not used should always be stored under this condition.
- Ensure that kit components do not freeze.
- Opened kits retain activity for 8 weeks if stored as described above.

7 PRECAUTIONS

- ELISA kits are only for in vitro diagnostic use by professionals.
- Serum and plasma samples should be treated as potentially infectious materials. Wear gloves and proper laboratory attire when handling sample materials. Do not eat, drink or smoke in areas where specimen or kit reagents are handled. Do not pipette with the mouth. In case of skin contact, wash with a germicidal soap and copious amounts of water. Seek medical advice if indicated.
- The PSA standards and controls are of human origin. They have been tested and confirmed negative for HIV, HBsAg and HCV. However, all standards should be treated as potential biohazards.
- Due to the potentially infectious character of samples and kit components all materials that have come in contact with these materials should be sterilized and disposed of according to local regulations. This also includes the liquid waste.
- The assay reagents contain preservatives, TMB, H₂O₂ or sulphuric acid and may be harmful if ingested. A direct skin or mucosa contact should be avoided. In case of skin contact, wash thoroughly with water and seek medical attention if required.
- The stop solution contains H₂SO₄. Since the H₂SO₄ used to terminate the color reaction is corrosive, the instrumentation employed to dispense it should be thoroughly cleaned after use.
- Do not interchange reagents from different LOT# or different suppliers.
- Avoid reagent or sample carry-over by using fresh tips for solutions and samples.
- Do not use test kit if zip lock pouch or bottles have been damaged.

8 GUIDELINE FOR SAMPLE COLLECTION; PREPARATION AND STORAGE

8.1 Sample collection

Blood samples are collected by veinpuncture. As different factors could influence the PSA level in blood, doctors should ensure that the patient has avoided the following conditions before taking the blood sample

The following conditions may lead to an increase of PSA levels

- biking
- sexual intercourse (ejaculation)
- Manipulation of the prostate during medical examinations like DRE, transrectal prostatic ultrasound etc. Please note that especially DRE might lead to an increase of f-PSA only, so that PSA carcinomas might be overlooked
- Prostatitis
- Liver dysfunction

The following conditions may lead to a decrease of PSA levels

- Intake of 5-alpha-reductaseinhibitors, antiandrogens, or GnRH analoga

8.2 Sample preparation

The preparation of serum or plasma samples is performed according to standard techniques. Serum or plasma should be prepared as soon as possible to avoid hemolysis and to improve the stability of PSA.

8.3 Storage of samples

For the assay either fresh serum or plasma samples can be used. F-PSA is not as stable as t-PSA.

Literature (Thomas, 2008) recommends to measure samples within 24 hours. In case of longer storage, freeze at -20 °C. A repeated freezing and thawing of samples should be avoided.

Note

- Highly lipemic or hemolytic samples can give incorrect analytical results.
- Samples must be free of microbial contaminations.
- Samples containing high titers of rheumatoid factor and human anti-mouse antibodies (HAMA) could give erroneous results.

9 ASSAY PROCEDURE

9.1 Reagent preparation

Before starting the assay the wash buffer must be diluted to the right concentration.

Per well approximately 2 mL of diluted buffer are needed. Calculate the volume of buffer required for the testing.

Take 1/40 of the volume wash buffer concentrate and dilute with 39/40 of the volume distilled water.

9.2 Test procedure

Note: It is highly recommended to perform all measurements as duplicates. An independent standard curve should be made for each series of measurements. For best results it is important that the solutions are always added to the wells in the same order to minimize variations in incubation times.

1. Prior to use bring all reagents, standards, controls, and samples to room temperature (18-25°C).
2. Check that all components are not expired and take care that bottles and plate (inclusive pouch) are not damaged.
3. Format the required microplate wells. Keep in mind that all measurements should be performed as duplicate. Document position of wells and respective samples, standards and controls to ensure later identification. Put any unused microwell modules back into the zip lock bag with the desiccant, seal bag and store at 2-8 °C.
4. Pipette 25 µL of standards, controls or samples into each well. Samples with an expected f-PSA value higher than 12 ng/ml should be diluted with the sample diluent.
5. Add 50 µL of Assay Reagent into each well and mix by moving plate on the table (10 sec.)
6. Incubate 1 h at room temperature (18-25°C)
7. Add 50 µL of conjugate into each well and mix by moving plate on the table (10 sec.)
8. Incubate 1 h at room temperature (18-25°C)
9. Remove solution from the wells by aspirating the liquid or by decanting it. If decanting, tap plate on adsorbent paper to remove residual liquid.
10. For washing fill plate with diluted wash buffer and wait 10sec before removing the buffer; repeat wash 5 times (for a total of 6 times)
We recommend the following procedure: wash wells 6-times with 250-300 µL/well diluted wash buffer. Preferably use an automated washing procedure. If washing manually, take care that the wash buffer remains in each well for the same time. This is necessary to receive lowest possible CV-values.
11. Pipette 100 µL TMB-substrate solution into each well
12. Incubate 15 min at room temperature (18-25°C)
13. Add 100 µL/well stopping solution (same order as substrate solution)
14. Read absorbencies (OD) at 450 nm (blinking 630nm)

9.3 Results


1. Calculate the mean absorbance for each duplicate.
2. Subtract the absorbance value of the zero standard from the mean absorbance values of standards, control and samples.
3. Draw the standard curve on lin-lin or log-log graph paper by plotting absorbance values of standards against appropriate PSA concentrations.
4. Read off the f-PSA concentrations for the control and the samples.

10 VALIDITY OF THE ASSAY

1. The OD 450 nm of the blanking well is lower than 0.150. Higher values indicate a chromogen/substrate contamination. In such a case, repeat the assay carefully checking the reagent.
2. The OD 450 nm of the highest standard (12 ng/mL) must be higher than 0.9. Lower values indicate kit or control decay. In such a case, check the expiry date of the kit before repeating the assay.
3. The control provided should not differ by more than 15% from the concentration stated on the label of the vial if run at least in duplicate.
4. Worksheet and standard curve of typical assay: Not to be used for calculation of actual test results.

Example

Wells	Identity	A 450 nm		Conc. ng/mL
1-2	St 0.00ng/mL	0.011	0,009	
3-4	St 0.75ng/mL	0.109	0.096	
5-6	St 1.5ng/mL	0.185	0.179	
7-8	St 3.0ng/mL	0.393	0.392	
9-10	St 6.0ng/mL	0.871	0.837	
11-12	St 12.0ng/mL	1.901	1.844	
13-14	control	0.218	0.231	2.05

Note that the absolute OD values for the standards might vary due to temperature influences or age of the conjugate. As long as the OD values form a standard curve and remain within the specifications and the control shows the expected value, results for unknown PSA samples are valid.

11 QUALITY CONTROL

- It is recommended that internal controls are used in every assay.
- Control results should be within established ranges and should represent preferably low, medium and high concentrations.

12 EXPECTED VALUES LIMITATIONS OF THE PROCEDURE

When total PSA values are between 4-10 ng/mL the ratio of f-PSA to t-PSA might be employed to increase the diagnostic specificity of PSA testing. Generally the likeliness of a cancerous condition increases the lower the f-PSA portion is.

A ratio of f-PSA to t-PSA of 25 % usually indicates a high probability of BPH (benign prostate hyperplasia). Low f-PSA likely signals prostate cancer. Most men with prostate cancer have an f- PSA value below 15 %. If free PSA is below 7% prostate cancer is most likely. According to the American Cancer Society and National Cancer Institute men with f-PSA below 7 % should undergo biopsy.

Please note that an isolated f-PSA concentration is of no diagnostic value.

For a detailed discussion of probabilities of prostate cancer and percent free PSA also see Gion et al. (1998).

The above values are only for user's guidance. If possible, it is recommended for each laboratory to establish its own specific values that take into consideration a population indigenous to the area where the laboratory is located.

Note: PSA values and the f-PSA to t-PSA can only be used to estimate the cancer risk. They should always be interpreted in conjunction with other clinical findings and should not be used as a sole basis for prostate cancer diagnosis.

13 PERFORMANCE CHARACTERISTICS

13.1 Detection limit

The limit of detection for this kit is 0.1 ng/mL

13.2 Precision

Intra- and inter-assay precisions were established by analysing three patient sera of different PSA concentrations.

The results are summarized in Tables 1 and 2.

Table 1) *Intra-assay precision*

Patients	Number of replicates	Mean ng/mL	SD ng/mL	CV %
1	10	0.31	0.024	7.7
2	10	1.75	0.107	6.1
3	10	9.63	0.623	6.5

Table 2) *Inter-assay precision*

Patients	Number of replicates	Mean ng/mL	SD ng/mL	CV %
1	10	2.86	0.17	6.1
2	10	0.75	0.06	7.8
3	10	2.1	0.16	7.8

13.3 Recovery

Spiked serum samples were prepared by adding aliquots of samples with highly elevated free PSA to normal male serum samples. The recovery of the antigen was in between 91.8 – 116% and in average 102%.

Table 3) *Recovery*

sample	Expected value (ng/mL)	Observed value (ng/mL)	Recovery %
1	5.48	5.47	99,8
	2.74	3.19	116
	5.48	5.03	91.8
2	2.78	2.81	101
3	5.57	5.68	102
	2.74	2.82	103

13.4 Hook effect

No hook effect has been noticed with samples up to 5000 ng/mL.

13.5 Correlation

The f-PSA ELISA (EIA-4189) was compared to a commercially available and CE marked two step Free PSA Elisa:

$$Y = 0.8922x + 0.11; R^2 = 0.8977$$

In a second study the f-PSA ELISA (EIA-4189) was compared to a second reference f-PSA ELISA:

$$Y = 1.0708 x; R^2 = 0.7969$$

13.6 Calibration

The free PSA ELISA (EIA-4189) is calibrated against WHO Standard 96/668.

13.7 Specificity

The assay is highly specific for free PSA, with a relatively low cross-reactivity to other proteins and polypeptides, lipids or chemotherapeutic agents that might be present in patient samples.

Table 4) Specificity

Antigens	Amount added	Cross reaction
Proteins		
AFP	10 µg/mL	No
CEA	10 µg/mL	No
HCG	10 µg/mL	No
Lactalbumin	10 µg/mL	No
PAP	1 µg/mL	No
Interfering substances		
Bilirubin	0.2 mg/mL	No
Triglyceride	15 mg/mL	No
Hemoglobin*	0.1 mg/ml	No
Chemotherapeutic Agents		No
Cyclophosphamid	800 µg/mL	No
Doxorubicin * HCl	20 µg/mL	No
Diethylstilbestrol	2 µg/mL	No
Flutamide	10 µg/mL	No
Methotrexate	50 µg/mL	No

* at higher concentration hemoglobin results in too high OD values, hemolytic samples should thus be avoided.

14 LEGAL ASPECTS

14.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

14.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 14.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

14.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 14.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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