User's Manual

Amyloid β 1-40 (human) ELISA

IVD

REF

EIA-523 I

96

Legal Manufacturer:

DRG

DRG Instruments GmbH, Germany

Division of DRG International, Inc

Frauenbergstr. 18, D-35039 Marburg

Telefon: +49 (0)6421-17000 Fax: +49-(0)6421-1700 50

Internet: www.drg-diagnostics.de

E-mail: drg@drg-diagnostics.de

Distributed by:

Contents / Inhaltsverzeichnis

1	INTRODUCTION.....	2
2	PRINCIPLE OF THE TEST	2
3	WARNINGS AND PRECAUTIONS	3
4	REAGENTS	4
5	SPECIMEN COLLECTION AND PREPARATION.....	6
6	ASSAY PROCEDURE	6
7	EXPECTED NORMAL VALUES.....	8
8	QUALITY CONTROL	8
9	PERFORMANCE CHARACTERISTICS.....	9
10	LIMITATIONS OF USE	10
11	LEGAL ASPECTS	10
12	REFERENCES / LITERATURE.....	11

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1 INTRODUCTION

1.1 Intended Use

The **DRG Amyloid β 1-40 (human) ELISA** (Abeta-40 ELISA) is an enzyme immunoassay for the quantitative in vitro measurement of human Amyloid β peptide 1-40 in serum, plasma and CSF.

1.2 Summary and Explanation

The pathological hallmarks of Alzheimer's disease (AD) are amyloid plaques, neurofibrillary tangles, synaptic degeneration and neuronal loss (1). Amyloid plaques are mainly composed of amyloid-beta ($A\beta$) 40 and 42 peptides derived from the proteolytic cleavage of amyloid precursor protein (APP) by β -site APP cleavage enzyme 1 (BACE1) (2,3) and the γ -secretase (4). The endosome and the endocytic pathway have been proposed as possible sites for the β and γ cleavage sites of APP (5), and the resulting $A\beta$ peptides are secreted by both neuronal and non-neuronal cells (6). Recently, soluble forms of $A\beta$ have been implicated in neurotoxicity (7), and may correlate better with cognition than amyloid plaque burden (8). As $A\beta$ is considered to play an early and pivotal role in AD pathogenesis, it may be a useful tool in diagnosing AD in the preclinical/early stages, as well as for monitoring potential $A\beta$ modifying therapies (9). While human CSF $A\beta$ levels have mostly shown reduction with disease progression (10), much of the data on human plasma $A\beta$ remain unclear. Studies have shown that subjects who exhibit overproduction of $A\beta$ such as individuals with familial forms of AD and Down's syndrome have higher plasma $A\beta$ levels compared to controls (11). In addition to variability found with age and disease severity (12), $A\beta$ plasma levels are further modulated by production in large peripheral organs such as skeletal muscle, liver, kidney, skin, and lung (13,14), movement from brain interstitial fluid into blood stream (15), and peripheral uptake and clearance of $A\beta$ in the liver by presence of low density lipoprotein receptor-related protein-1 (16) and receptors for advanced glycation end products (RAGE) (17). This may explain why much of the plasma $A\beta$ data in cross-sectional studies do not show significant differences between sporadic AD and controls. Longitudinal studies seem more promising, showing initially elevated $A\beta$ 42 levels in those that eventually develop AD in the future (18). Furthermore, subjects who manifest sporadic AD within the time frame of the study start with higher plasma $A\beta$ 42 levels, and the $A\beta$ 42/ $A\beta$ 40 ratios also appear to be significantly different from those who remain asymptomatic (19).

2 PRINCIPLE OF THE TEST

The DRG Abeta-40 ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the sandwich principle.

The microtiter wells are coated with a monoclonal [mouse] antibody directed towards a unique antigenic site on the Abeta-40 peptide. An aliquot of patient sample containing endogenous Abeta-40 peptide is incubated in the coated well with enzyme conjugate, which is a biotinylated anti-Abeta-40 antibody. After incubation the unbound conjugate is washed off. Finally, Enzyme Complex, which is streptavidin conjugated with horseradish peroxidase is added, and after incubation, unbound enzyme complex is washed off.

The amount of bound peroxidase is proportional to the concentration of Abeta-40 peptide in the sample.

Having added the substrate solution, the intensity of colour developed is proportional to the concentration of Abeta-40 in the patient sample.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21 °C to 26°C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with *Stop Solution* containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
18. Some reagents may contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Material Safety Data Sheets. Material Safety Data Sheets for this product are available upon request directly from DRG.

4 REAGENTS

4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;
Wells coated with anti-Abeta-40 antibody (mouse monoclonal).
2. **Standard (Standard 0-6)**, 7 vials, (lyophilized), 1 mL each,
Concentration: 0 - 20 - 50 - 125 - 250 - 500 - 750 pg/mL
see „Reagent Preparation“
Contain non-mercury preservative.
3. **Control Low & High**, 2 vials, (lyophilized) 1 mL each,
see „Reagent Preparation“
For control values and ranges please refer to vial label or QC-Datasheet.
Contain non-mercury preservative.
4. **Assay Buffer**, 1 vial, 14 mL, ready to use,
Contains non-mercury preservative.
5. **Enzyme Conjugate**, 1 vial, 7 mL, ready to use,
Anti-Abeta-40 antibody conjugated to horseradish peroxidase;
Contains non-mercury preservative.
6. **Enzyme Complex**, 1 vial, 14 mL, ready to use,
Streptavidin conjugated to horseradish peroxidase;
Contains non-mercury preservative.
7. **Substrate Solution**, 1 vial, 14 mL, ready to use,
Tetramethylbenzidine (TMB).
8. **Stop Solution**, 1 vial, 14 mL, ready to use,
contains 0.5 M H₂SO₄,
Avoid contact with the stop solution. It may cause skin irritations and burns.
9. **Wash Solution**, 1 vial, 30 mL (40X concentrated),
see „ Reagent Preparation “.

Note: Additional *Assay Buffer* for sample dilution is available upon request.

4.2 Materials required but not provided

- A microtiter plate calibrated reader (450 ± 10 nm) (e.g. the DRG Instruments Microtiter Plate Reader).
- Calibrated variable precision micropipettes.
- Absorbent paper.
- Distilled or deionized water
- Timer
- Semi logarithmic graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for two months if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature prior to use.

Standard / Control

Reconstitute the lyophilized content with 1 mL distilled water and let stand for 10 minutes in minimum. Mix the standards and controls several times before use.

Note: *The reconstituted standards and controls can be stored for 1 week at 2 °C to 8 °C.
For longer storage, reconstituted controls should be apportioned and stored at -20 °C.*

Wash Solution

Add distilled water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated Wash Solution with 1170 mL distilled water to a final volume of 1200 mL.
The diluted Wash Solution is stable for 2 weeks at room temperature.

4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Material Safety Data Sheets (see chapter 13).

4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum, plasma (EDTA-, heparin- or citrate plasma) and spinal cord fluid (CSF) can be used in this assay. Do not use haemolytic, icteric or lipaemic specimens.

Please note: Samples containing sodium azide should not be used in the assay.

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette # 02.1388.001), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma:

Whole blood should be collected into centrifuge tubes containing anti coagulant and centrifuged immediately after collection.

(E.g. for EDTA plasma Sarstedt Monovette – red cap - # 02.166.001;
for Heparin plasma Sarstedt Monovette – orange cap - # 02.165.001;
for Citrate plasma Sarstedt Monovette – green cap - # 02.167.001.)

CSF:

Collect CSF by lumbar puncture.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 1 day at 2 °C to 8 °C prior to assaying. Specimens held for a longer time (up to 12 months) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

CSF:

CSF should be diluted **at least 50-fold** with *Assay Buffer* **before** analysis.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

Dilution 1:50: 5 μ L CSF + 245 μ L *Assay Buffer*

Serum and plasma:

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Assay Buffer* and reassayed as described in *Assay Procedure*.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10: 25 μ L sample + 225 μ L *Assay Buffer* (mix thoroughly)
b) dilution 1:100: 25 μ L dilution a) 1:10 + 225 μ L *Assay Buffer* (mix thoroughly).

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Test Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the holder.
2. Dispense **50 μ L Assay Buffer** into each well.
3. Dispense **100 μ L** of each **Standard, Control** and **sample** with new disposable tips into appropriate wells.
4. Dispense **50 μ L Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
5. Incubate for **180 minutes** at room temperature on a plate shaker with \sim 700 rpm.
6. Briskly shake out the contents of the wells.
Rinse the wells **4 times** with diluted *Wash Solution* (400 μ L per well). Strike the wells sharply on absorbent paper to remove residual droplets.
Important note:
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
7. Add **100 μ L of Enzyme Complex** to each well.
8. Incubate for **30 minutes** at room temperature.
9. Rinse the wells **4 times** with diluted *Wash Solution* (400 μ L per well). Strike the wells sharply on absorbent paper to remove residual droplets.
10. Add **100 μ L of Substrate Solution** to each well.
11. Incubate for **20 minutes** at room temperature.
12. Stop the enzymatic reaction by adding **100 μ L of Stop Solution** to each well.
13. Determine the absorbance (OD) of each well at **450 \pm 10 nm** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Using linear graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred method. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as $>$ 750 pg/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Units (450 nm)
Standard 0 (0 pg/mL)	0.09
Standard 1 (20 pg/mL)	0.14
Standard 2 (50 pg/mL)	0.22
Standard 3 (125 pg/mL)	0.47
Standard 4 (250 pg/mL)	0.77
Standard 5 (500 pg/mL)	1.51
Standard 6 (750 pg/mL)	2.27

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently normal healthy adults, using the DRG Abeta-40 ELISA the following values are observed:

Population	Valid N	Age (years)	Median (pg/mL)	Mean (pg/mL)	5% - 95% Percentile
apparently healthy males and females	40	25-52	144.0	135.1	72.8 – 203.1

The results are in good agreement with published normal values (80 - 300 pg/mL).

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 1.66 pg/mL – 750 pg/mL.

9.2 Specificity of Antibodies (Cross Reactivity)

No cross reactivity was observed with related peptides Abeta 1-42 and Abeta 12-28.

9.3 Sensitivity

The analytical sensitivity was calculated by adding 2 standard deviations to the mean of 20 replicate analyses of *Standard 0* and was found to be 1.66 pg/mL.

9.4 Reproducibility

9.4.1 Intra Assay

The within assay variability is shown below:

Sample	n	Mean (pg/mL)	CV (%)
1	20	84.0	9.3
2	20	161.7	7.5
3	20	280.9	9.3

9.4.2 Inter Assay

The between assay variability is shown below:

Sample	n	Mean (pg/mL)	CV (%)
1	40	100.0	11.8
2	40	218.5	11.2
3	40	273.0	9.5

9.5 Recovery

Samples have been spiked by adding Abeta-40 solutions with known concentrations in a 1:1 ratio. The expected values were calculated by addition of half of the values determined for the undiluted samples and half of the values of the known solutions. The % Recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100.

		Sample 1 (Serum)	Sample 2 (RPMI + 10% FCS)	Sample 3 (EDTA plasma)
Concentration [pg/mL]		222.5	171.4	342.2
Average Recovery		91.8	104.8	93.7
Range of Recovery [%]	from	90.3	96.4	89.8
	to	94.2	110.5	97.4

9.6 Linearity

		Sample 1 (Serum)	Sample 2 (RPMI + 10% FCS)	Sample 3 (EDTA plasma)
Concentration [pg/mL]		442.0	671.9	533.9
Average Recovery		107.9	89.9	102.3
Range of Recovery [%]	from	103.6	85.8	97.7
	to	112.1	93.5	112.4

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 30 mg/mL) have no influence on the assay results.

The assay contains reagents to minimize interference of HAMA and heterophilic antibodies. However, extremely high titers of HAMA or heterophilic antibodies may interfere with the test results.

10.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of Abeta-40 in a sample.

10.3 High-Dose-Hook Effect

No hook effect was observed in this test up to 96,000 pg/mL of Abeta-40.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test. The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12 REFERENCES / LITERATURE

1. Price D.L., Sisodia S.S.: Mutant genes in familial Alzheimer's disease and transgenic models. *Annu. Rev. Neurosci.*; 1998; 21:479–505.
2. Sinha S. et al.: Purification and cloning of amyloid precursor protein beta-secretase from human brain. *Nature*; 1999; 402: 537–540.
3. Vassar R. et al.: Beta-secretase cleavage of Alzheimer's amyloid precursor protein by the transmembrane aspartic protease BACE. *Science* 1999; 286:735–741.
4. De Strooper B.: Aph-1, Pen-2, and Nicastrin with Presenilin generate an active gamma-Secretase complex. *Neuron*; 2003; 38:9–12.
5. Small S.A., Gandy S.: Sorting through the cell biology of Alzheimer's disease: intracellular pathways to pathogenesis. *Neuron* 2006; 52:15–31.
6. Selkoe D.J.: Clearing the Brain's Amyloid Cobwebs. *Neuron*; 2001; 32:177–180.
7. Lesne S. et al.: A specific amyloidbeta protein assembly in brain impairs memory. *Nature*; 2006; 440:352-357.
8. McLean C.A. et al.: Soluble pool of Abeta amyloid as a determinant of severity of neurodegeneration in Alzheimer's disease. *Ann. Neurol.*; 1999; 46:860–866.
9. Galasko D.: Biomarkers for Alzheimer's disease-clinical needs and application. *J Alzheimers Dis* 2005; 8:339–346.
10. Jensen M. et al.: Cerebrospinal fluid A beta42 is increased early in sporadic Alzheimer's disease and declines with disease progression. *Ann. Neurol.*; 1999; 45:504–511.
11. Iwatsubo T.: Amyloid beta protein in plasma as a diagnostic marker for Alzheimer's disease. *Neurobiol. Aging*; 1998; 19:161–163.
12. Giedraitis V. et al.: The normal equilibrium between CSF and plasma amyloid beta levels is disrupted in Alzheimer's disease. *Neurosci. Lett.*; 2007; 427:127–131.
13. Vassar R. et al.: Beta-secretase cleavage of Alzheimer's amyloid precursor protein by the transmembrane aspartic protease BACE. *Science*; 1999; 286:735–741.
14. Sinha S. et al.: Purification and cloning of amyloid precursor protein beta-secretase from human brain. *Nature*; 1999; 402:537–540.
15. Shibata M. et al.: Clearance of Alzheimer's amyloid-ss(1–40) peptide from brain by LDL receptor-related protein-1 at the blood-brain barrier. *J. Clin. Invest.* 2000; 106:1489–1499.
16. Deane R. et al.: LRP/amyloid beta-peptide interaction mediates differential brain efflux of Abeta isoforms. *Neuron*; 2004; 43:333–344.
17. Deane R. et al.: RAGE mediates amyloid-beta peptide transport across the blood-brain barrier and accumulation in brain. *Nat. Med.*; 2003; 9:907–913.
18. Mayeux R. et al.: Plasma A[beta]40 and A[beta] 42 and Alzheimer's disease: relation to age, mortality, and risk. *Neurology*; 2003; 61:1185–1190.
19. Graff-Radford N.R. et al.: Association of low plasma Abeta42/Abeta40 ratios with increased imminent risk for mild cognitive impairment and Alzheimer disease. *Arch.Neurol.*; 2007; 64:354–362.

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