## Instructions for Use

# Hepcidin 25 (bioactive) HS ELISA

IVD

REF EIA-5782

$\Sigma$  96

**DRG**

DRG Instruments GmbH, Germany  
Frauenbergstraße. 18, D-35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

Distributed by:

**DRG**

DRG International, Inc., USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Contents / Inhaltsverzeichnis / Contenuti / Contenido / Contenu**

1	INTRODUCTION .....	2	1	INTRODUCCIÓN .....	27	
2	PRINCIPLE OF THE TEST .....	2	2	FUNDAMENTO DEL ENSAYO .....	27	
3	WARNINGS AND PRECAUTIONS .....	3	3	PRECAUCIONES .....	27	
4	REAGENTS.....	4	4	COMPONENTES DEL KIT.....	28	
5	SPECIMEN COLLECTION AND PREPARATION.....	6	5	MUESTRAS .....	30	
6	ASSAY PROCEDURE.....	6	6	PROCEDIMIENTO DE ENSAYO .....	30	
7	EXPECTED NORMAL VALUES.....	8	7	VALORES ESPERADOS.....	32	
8	QUALITY CONTROL.....	8	8	CONTROL DE CALIDAD .....	32	
9	PERFORMANCE CHARACTERISTICS.....	9	9	CARACTERÍSTICAS DEL ENSAYO.....	32	
10	LIMITATIONS OF USE.....	10	10	LIMITACIONES DE USO .....	33	
11	LEGAL ASPECTS .....	11	11	ASPECTOS LEGALES .....	33	
12	REFERENCES / LITERATURE .....	11	12	REFERENCIAS / BIBLIOGRAFÍA.....	33	
1	EINLEITUNG .....	12	1	INTRODUCTION.....	34	
2	TESTPRINZIP .....	12	2	PRINCIPE DU TEST .....	34	
3	VORSICHTSMAßNAHMEN.....	13	3	PRECAUTIONS D'UTILISATION.....	34	
4	BESTANDTEILE DES KITS .....	14	4	COMPOSITION DU KIT .....	35	
5	PROBENVORBEREITUNG.....	16	5	ECHANTILLON .....	37	
6	TESTDURCHFÜHRUNG .....	16	6	RÉALISATION DU TEST .....	37	
7	ERWARTETE WERTE .....	18	7	VALEURS ATTENDUES.....	39	
8	QUALITÄTS-KONTROLLE.....	18	8	CONTROLE DE QUALITE .....	39	
9	ASSAY CHARACTERISTIKA .....	18	9	CARACTERISTIQUES DU TEST .....	39	
10	GRENZEN DES TESTS .....	19	10	LIMITES D'UTILISATION.....	40	
11	RECHTLICHE GRUNDLAGEN .....	19	11	ASPECTS LEGAUX.....	40	
12	REFERENZEN / LITERATUR .....	19	12	REFERENCES / BIBLIOGRAPHIE .....	40	
1	INTRODUZIONE .....	20	SYMBOLS USED.....			41
2	PRINCIPIO DEL TEST .....	20				
3	PRECAUZIONI .....	20				
4	COMPONENTI DEL KIT.....	21				
5	CAMPIONI.....	23				
6	ATTUAZIONE DEL TEST.....	23				
7	VALORI NORMALI .....	25				
8	CONTROLLO QUALITÀ.....	25				
9	CARATTERISTICHE DEL TEST .....	25				
10	LIMITAZIONE DEL TEST .....	26				
11	ASPETTI LEGALI .....	26				
12	BIBLIOGRAFIA.....	26				

## 1 INTRODUCTION

### 1.1 Intended Use

The **DRG Hepcidin 25 (bioactive) HS ELISA** is a **high sensitive** enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of Hepcidin-25 in serum or plasma (EDTA-, heparin- or citrate plasma).

### 1.2 Summary and Explanation

Hepcidin is an iron homeostasis regulator peptide. The bioactive peptide Hepcidin-25 is generated predominantly in the liver by proteolytic cleavage of the C-terminal 25 amino acids of prohepcidin (1). Subsequent N-terminal processing of Hepcidin-25 results in smaller peptides of 20-24 amino acids that show greatly reduced activity and accumulate in the urine (2).

Although originally identified as antimicrobial peptide (3), Hepcidin-25 is now established as a major regulator of dietary iron absorption and cellular iron release (4). Hepcidin exerts its regulatory function by counteracting the function of ferroportin, the major cellular iron exporter in the membrane of macrophages, hepatocytes, and the basolateral site of enterocytes. Hepcidin-25 induces the internalization and degradation of ferroportin, resulting in increased intracellular iron stores, decreased dietary iron absorption, and decreased circulating iron concentrations (5).

Hepatocellular hepcidin synthesis decreases under conditions of increased demand for circulating iron like iron deficiency, hypoxia, anemia, and erythropoiesis. In contrast, hepcidin synthesis is induced by inflammation and infection (6).

Serum Hepcidin-25 has been shown to add value to identify and differentiate specific disease conditions. Hepcidin deficiency causes hereditary hemochromatosis, characterized by body iron overload that may progress to liver cirrhosis (7). In addition, low Hepcidin-25 concentration can be induced by iron loading anemias (8) and chronic hepatitis C (9). In contrast, high Hepcidin-25 levels have been found in iron-refractory iron-deficiency anemia (10), during infection, chronic kidney disease (11), and after intensive exercise, explaining the high iron deficiency among athletes (12).

#### Product Patent Protected:

Pat. US 7,320,894 B2; US 8,017,737 B2, EP 2109624, EP 1578254,

Japan 4638350, Russia 2 359 268 C2, China 200380108964.8, Hong Kong 1114419, Canada 2,506,668

## 2 PRINCIPLE OF THE TEST

The DRG Hepcidin 25 (bioactive) HS ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the principle of competitive binding.

The microtiter wells are coated with a monoclonal (mouse) antibody directed towards an antigenic site of the Hepcidin-25 molecule. Endogenous Hepcidin-25 of a patient sample competes with a Hepcidin-25-biotin conjugate (Enzyme Conjugate) for binding to the coated antibody.

After incubation the microtiter plate is washed to stop the competition reaction. In the following incubation the bound biotin molecules are detected with streptavidin peroxidase (Enzyme Complex).

After incubation the plate is washed the second time.

After addition of the substrate solution, the intensity of colour developed is inversely proportional to the concentration of Hepcidin-25 in the patient sample.

### 3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21 °C - 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with *Stop Solution* containing 0.5 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

## 4 REAGENTS

### 4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;  
Wells coated with anti-Hepcidin-25 antibody (monoclonal).
2. **Standard (Standard 0 - 5)**, 6 vials, lyophilized, 0.5 mL;  
Concentrations: 0 – 1 – 3 – 9 – 27 – 81 ng/mL  
Conversion: 1 ng/mL = 0.358 nmol/L  
See "Reagent Preparation".  
Contain non-mercury preservative.
3. **Control Low & High**, 2 vials, lyophilized, 0.5 mL;  
For control values and ranges please refer to vial label or QC-Datasheet.  
See "Reagent Preparation".  
Contain non-mercury preservative.
4. **Sample Diluent**, 1 vial, 3 mL, ready to use;  
Contains non-mercury preservative
5. **Enzyme Conjugate**, 1 vial, 7 mL, ready to use,  
Hepcidin-25 conjugated to biotin;  
Contains non-mercury preservative.
6. **Enzyme Complex**, 1 vial, 14 mL, ready to use,  
Streptavidin conjugated to horseradish peroxidase  
Contains non-mercury preservative.
7. **Substrate Solution**, 1 vial, 14 mL, ready to use,  
Tetramethylbenzidine (TMB).
8. **Stop Solution**, 1 vial, 14 mL, ready to use,  
contains 0.5 M H<sub>2</sub>SO<sub>4</sub>,  
Avoid contact with the stop solution. It may cause skin irritations and burns.
9. **Wash Solution**, 1 vial, 30 mL (40X concentrated),  
See "Reagent Preparation".

**Note:** Additional *Sample Diluent* for sample dilution is available upon request.

### 4.2 Materials required but not provided

- A microtiter plate calibrated reader (450 ± 10 nm) (e.g. the DRG Instruments Microtiter Plate Reader).
- Calibrated variable precision micropipettes.
- Absorbent paper.
- Distilled or deionized water
- Timer
- Semi logarithmic graph paper or software for data reduction

### 4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for two months if stored as described above.

#### 4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature prior to use.

##### **Standards**

Reconstitute the lyophilized contents of each standard vial with 0.5 mL deionized water and let stand for 10 minutes in minimum. Mix several times before use.

**Note:** *The reconstituted standards are stable for 2 months at 2 °C - 8 °C.  
For longer storage freeze at -20 °C.*

##### **Control**

Reconstitute the lyophilized content with 0.5 mL deionized water and let stand for 10 minutes in minimum. Mix the control several times before use.

**Note:** *The reconstituted controls are stable for 2 months at 2 °C - 8 °C.  
For longer storage freeze at -20 °C.*

##### **Wash Solution**

Add deionized water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL deionized water to a final volume of 1200 mL.

*The diluted Wash Solution is stable for 2 weeks at room temperature.*

#### 4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Material Safety Data Sheet.

#### 4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

## 5 SPECIMEN COLLECTION AND PREPARATION

Serum or plasma (EDTA-, heparin- or citrate plasma) can be used in this assay.

Do not use haemolytic, icteric or lipaemic specimens.

*Please note:* Samples containing sodium azide should not be used in the assay.

### 5.1 Specimen Collection

#### Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

#### Plasma:

Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

### 5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 4 days at 2 °C to 8 °C prior to assaying.

Specimens held for a longer time (up to 12 months) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

### 5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Sample Diluent* and reassayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

#### Example:

- a) dilution 1:10:           10 µL sample + 90 µL *Sample Diluent* (mix thoroughly)
- b) dilution 1:100:        10 µL dilution a) 1:10 + 90 µL *Sample Diluent* (mix thoroughly).

## 6 ASSAY PROCEDURE

### 6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

## 6.2 Test Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the frame holder.
2. Dispense **20 µL** of each **Standard, Control** and **samples with new disposable tips** into appropriate wells.
3. Dispense **50 µL Enzyme Conjugate** into each well.  
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Incubate for **60 minutes** at room temperature.
5. Briskly shake out the contents of the wells. Rinse the wells  
**4 x with 400 µL diluted Wash Solution** per well (if a plate washer is used) - or  
**4 x with 300 µL diluted Wash Solution** per well for manual washing.  
Strike the wells sharply on absorbent paper to remove residual droplets.  
**Important note:**  
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
6. Dispense **100 µL of Enzyme Complex** into appropriate wells.
7. Incubate for **30 minutes** at room temperature.
8. Briskly shake out the contents of the wells. Rinse the wells  
**4 x with 400 µL diluted Wash Solution** per well (if a plate washer is used) - or  
**4 x with 300 µL diluted Wash Solution** per well for manual washing.  
Strike the wells sharply on absorbent paper to remove residual droplets.
9. Add **100 µL of Substrate Solution** to each well.
10. Incubate for **20 minutes** at room temperature.
11. Stop the enzymatic reaction by adding **100 µL of Stop Solution** to each well.
12. Determine the absorbance (OD) of each well at **450 ± 10 nm** with a microtiter plate reader.  
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

## 6.3 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Using semi-logarithmic graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 81 ng/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

### 6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Units (450 nm)
Standard 0 (0 ng/mL)	2.03
Standard 1 (1 ng/mL)	1.70
Standard 2 (3 ng/mL)	1.29
Standard 3 (9 ng/mL)	0.77
Standard 4 (27 ng/mL)	0.37
Standard 5 (81 ng/mL)	0.17

## 7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently normal healthy adults, using the DRG Hepcidin 25 (bioactive) HS ELISA the following values are observed:

Population	Valid N	Range (ng/mL)	Mean (ng/mL)	2.5 <sup>th</sup> - 97.5 <sup>th</sup> Percentile (ng/mL)	Median (ng/mL)
Males and Females	75	0.25 - 47.66	16.45	1.49 - 41.46	13.47

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

## 8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

## 9 PERFORMANCE CHARACTERISTICS

### 9.1 Assay Dynamic Range

The range of the assay is between 0.153 ng/mL - 81 ng/mL.

### 9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross reactivity of the assay:

Analyte	% Cross-Reactivity
Prohepcidin	< 0.001
Insulin	< 0.001
Hepcidin-22	24.2
Hepcidin-20	87.7

### 9.3 Sensitivity

The analytical sensitivity of the DRG ELISA was calculated by subtracting 2 standard deviations from the mean of 20 replicate analyses of the Zero Standard and was found to be 0.153 ng/mL.

### 9.4 Reproducibility

#### 9.4.1 Intra Assay

The within assay variability is shown below:

Sample	n	Mean (ng/mL)	CV (%)
1	20	5.78	5.31
2	20	9.48	5.68
3	20	20.20	5.50

#### 9.4.2 Inter Assay

The between assay variability is shown below:

Sample	n	Mean (ng/mL)	CV (%)
1	40	3.20	14.35
2	40	21.93	9.50
3	40	59.63	13.62

#### 9.4.3 Inter-Lot

The inter-assay (between-lots) variation was determined by repeated measurements of 3 samples in 3 different kit lots.

Sample	n	Mean (ng/mL)	CV (%)
1	18	0.95	14.99
2	18	22.67	7.47
3	18	30.25	0.62

**9.5 Recovery**

Samples have been spiked by adding Hepcidin 25 solutions with known concentrations in a 1:1 ratio. The % recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100 (expected value = (endogenous Hepcidin 25 + added Hepcidin 25) / 2; because of a 1:2 dilution of serum with spike material).

	Sample 1	Sample 2	Sample 3
Concentration [ng/mL]	20.8	27.3	59.2
Average Recovery [%]	99.6	94.4	98.0
Range of Recovery [%]	from	88.2	88.4
	to	111.8	105.8

**9.6 Linearity**

	Sample 1	Sample 2	Sample 3
Concentration [ng/mL]	14.7	27.3	59.0
Average Recovery [%]	103.1	92.1	98.4
Range of Recovery [%]	from	92.5	87.9
	to	114.6	95.2

**10 LIMITATIONS OF USE**

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

**10.1 Interfering Substances**

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

**10.2 Drug Interferences**

Until today no substances (drugs) are known to us, which have an influence to the measurement of Hepcidin 25 in a sample.

**10.3 High-Dose-Hook Effect**

No hook effect was observed in this test.

## 11 LEGAL ASPECTS

### 11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

### 11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient. Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

## 12 REFERENCES / LITERATURE

1. Scamuffa N, et al. Regulation of prohepcidin processing and activity by the subtilisin-like proprotein convertases Furin, PC5, PACE4 and PC7. *Gut*. 2008; 57(11):1573-82.
2. Kemna EH, Tjalsma H, Podust VN, Swinkels DW. Mass spectrometry-based hepcidin measurements in serum and urine: analytical aspects and clinical implications. *Clin. Chem*. 2007; 53(4):620-8.
3. Krause A et al. LEAP-1, a novel highly disulfide-bonded human peptide, exhibits antimicrobial activity. *FEBS Lett*. 2000; 480(2-3):147-50.
4. Kroot J et al. Hepcidin in Human Iron Disorders: Diagnostic Implications. *Clin. Chem*. 2011; 57:121650–1669.
5. Nemeth E, Ganz T. The role of hepcidin in iron metabolism. *Acta Haematol*. 2009; 122(2-3):78-86.
6. Hentze MW, Muckenthaler MU, Galy B, Camaschella C. Two to tango: regulation of Mammalian iron metabolism. *Cell*. 2010; 142:24–38.
7. Swinkels DW, Jorna AT, Raymakers RA. Synopsis of the Dutch multidisciplinary guideline for the diagnosis and treatment of hereditary haemochromatosis. *Neth. J. Med*. 2007; 65(11):452-5.
8. Camaschella C, Silvestri L. Molecular mechanisms regulating hepcidin revealed by hepcidin disorders. *Scientific World Journal*. 2011; 11:1357-66.
9. Girelli D et al. Reduced serum hepcidin levels in patients with chronic hepatitis C. *J. Hepatol*. 2009; 51(5):845-52.
10. Gardenghi S et al. Ineffective erythropoiesis in beta-thalassemia is characterized by increased iron absorption mediated by downregulation of hepcidin and up-regulation of ferroportin. *Blood*. 2007; 109:5027–35.
11. Uehata T et al. Serum hepcidin-25 levels and anemia in non-dialysis chronic kidney disease patients: a cross-sectional study. *Nephrol Dial. Transplant*. 2012; 27(3):1076-83.
12. Peeling P, Dawson B, Goodman C, Landers G, Trinder D. Athletic induced iron deficiency: new insights into the role of inflammation, cytokines and hormones. *Eur. J. Appl. Physiol*. 2008; 103: 381–91.

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