DRG:HYBRID•XL[®]

Instructions for Use

HbA1c

IVD

REF

HYC-5325

80 tests

DRG Instruments GmbH, Germany
Frauenbergstraße. 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

ASSAY PROTOCOL BARCODE (APB)

(Version 2.60 Software or later / ab Softwareversion 2.60 oder höher / Versione Software 2.60 o superiore /
Versión de Software 2.60 o posterior)

HYC-5325- v2.60.3

The barcode must be used to install the assay protocol into the DRG:HYBRiD-XL software via the SCAN NEW LOT page.

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Please refer to section **3: Routine Procedures: "Installing a new assay product"** of the User Manual v2.60 or later.

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Contents / Inhaltsverzeichnis / Contenuti / Contenido / Contenu

1	INTRODUCTION	4
2	PRINCIPLE OF THE TEST	4
3	WARNINGS AND PRECAUTIONS	4
4	REAGENTS	4
5	SPECIMEN COLLECTION AND PREPARATION	5
6	ASSAY PROCEDURE	5
7	QUALITY CONTROL	6
8	EXPECTED NORMAL VALUES	6
9	LIMITATIONS OF USE	6
10	PERFORMANCE CHARACTERISTICS	6
11	METHOD COMPARISON	6
12	LEGAL ASPECTS	7

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1 INTRODUCTION

1.1 Intended Use

The **DRG:HYBRID-XL HbA1c** is an immunoturbidimetric assay for the quantitative *in vitro diagnostic* measurement of Hemoglobin A1c (HbA1c) in **whole blood**.

Only for use with the DRG:HYBRID-XL Analyzer.

1.2 Summary and Explanation

HbA1c is a glycosylated form of haemoglobin formed by the attachment of glucose residues in the blood to the haemoglobin molecules. In the diabetic population where blood glucose levels are abnormally elevated the level of HbA1c also increases.

The level of HbA1c is proportional to the level of glucose in the blood and has been widely accepted as an indicator of the mean blood glucose concentration in the preceding 6-8 weeks. It is therefore a long-term indicator of diabetic control. For routine use HbA1c levels should be monitored every 3-4 months. However in gestational diabetes and after a change in therapy it may be useful to measure HbA1c more frequently at 2-4 week intervals.

2 PRINCIPLE OF THE TEST

The DRG:HYBRID-XL HbA1c is an **immunoturbidimetric assay**.

The sample (haemolysis sample) is added to the unsensitised latex particles, and the surface of the latex adsorbs total haemoglobin in the sample.

Anti-human HbA1c mouse monoclonal antibodies interact with the latex-bound HbA1c. Reaction with anti-mouse IgG goat antibody leads to an amplification of this reaction and the forming of agglutination complexes.

The amount of agglutination caused depends on the amount of HbA1c that adsorbs the surface of the latex, this agglutination is measured as a turbidity.

The concentration of HbA1c (mmol/mol) in the sample is determined by referring to the calibration curve.

3 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only.
For professional use only.
- This kit can only be used in combination with the DRG:HYBRID-XL Analyzer.
- Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is clear and understood.
- Do not remove, exchange, discard or damage any of the barcode labels provided with each kit and its components. All barcodes build an integral system for the kit lot.**
- Respect the general safety measures for use of laboratory reagents.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear appropriate disposable gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may cause false results.
- Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Unused reagent cartridges must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant provided.
- Optimal test results are only obtained when using calibrated pipettes.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange reagent cartridges of different kits even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the wells in the reagent cartridges may differ slightly.
- Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

- For information on hazardous substances included in the kit please refer to Safety Data Sheets.
For professional users the Safety Data Sheet for this product is available upon request directly from DRG.

4 REAGENTS

4.1 Reagents provided

4.1.1 Reagent Cartridges

80 pieces containing the following:

- Reagent 1**, 260 µL
Suspension of latex particles;
Contains non-mercury preservative.
- Reagent 2**, 130 µL
Buffer solution containing anti-human HbA1c antibody (monoclonal, mouse) and anti-mouse IgG goat antibody
Contains non-mercury preservative.

4.1.2 Re-Calibrator 1 & 2

2 vials, 0.5 mL each, lyophilized;

For re-calibration of the quantitative DRG:HYBRID-XL HbA1c test.

Concentrations are lot-specific.

Re-Calibrators are standardised against the following reference material: IFCC reference material
See "Reagent Preparation".

Contain non-mercury preservative.

4.1.3 Control 1 & 2

2 vials, 0.5 mL each, lyophilized;

For control values and ranges please refer to the bar code on vial label or to the QC-Datasheet.

See "Reagent Preparation".

Contain non-mercury preservative.

4.2 Materials required but not provided

- General needed laboratory equipment
- Ultra-pure water
DRG recommends to use Clinical Laboratory Reagent Water (CLRW) according to CLSI guideline 3C-A4 with the following specifications:
Resistivity at 25 °C [MΩ·cm]: > 10
Conductivity at 25 °C [µS·cm⁻¹]: < 0.1
Total Organic Carbon/p.p.b.[µg/L]: < 50
Colloids [µg/mL]: <0.05
- [REF]** HYB-5252 DRG:HYBRID-XL Analyzer
- [REF]** HYI-5392: *System Solution 5L*, 5000 mL;
(Instrument Feed Water according to CLSI guideline 3C-A4 with the following specification can also be used:
Resistivity at 25 °C [MΩ·cm]: > 1
Conductivity at 25 °C [µS·cm⁻¹]: < 1
Total Organic Carbon/p.p.b.[µg/L]: < 200
Colloids [µg/mL]: <0.1)
- [REF]** HYI-5394: *Wash Buffer*, 40x concentrate, 25 mL
- [REF]** HYI-5395: *Needle Cleaning Solution*, 30 mL. Cleaning solution for the pipetting needle (daily and weekly maintenance, see also user manual)
- [REF]** HYI-5387: *Cuvettes*, (2 x 360 pieces)

For use of the *Secondary Sample Holder* for secondary tubes the following tubes are required:

- [REF]** HYI-5391: *Sample Tubes (Secondary)*, 2500 pcs.

4.3 Storage Conditions

All kit components should be stored at 2 °C to 8 °C to ensure product performance until the defined expiry date.

When stored at 2 °C to 8 °C, **unopened kits** will retain reactivity until expiration date. Do not use reagents beyond this date.

- Cartridges (stored at 2 °C to 8 °C) in the supplied and unopened zip/foil bags will retain reactivity until expiration date.
- Unopened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity until expiration date.

Opened reagents and the reagent cartridges must be stored at 2 °C to 8 °C.

Once the plastic bag has been opened, care should be taken to tightly close it again along with the supplied desiccant bag.

Immediately after end of each run the Re-Calibrator and Control vials have to be removed from the instrument, tightly capped and stored at 2 °C to 8 °C.

- Unused cartridges in opened zip/foil bags (stored at 2 °C to 8 °C) will retain reactivity until expiration date, if stored as described above.
- Pierced or open cartridges must be disposed of immediately.

Opened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity for 8 weeks.

4.3.1 On-board Stability

For Re-Calibrators and Controls the on-board stability has been evaluated under controlled laboratory conditions at room temperature (20 °C to 25 °C).

Due to the differences in laboratory environmental conditions and reagent volumes, the on-board stability may deviate from the declared value.

On-board stability:	48 h
---------------------	------

4.4 Reagent Preparation

Bring all reagents, such as controls and re-calibrators, to room temperature (20 °C to 25 °C) prior to use. Reagent Cartridges can be used directly after storage in the refrigerator.

Re-Calibrator 1 & 2

Reconstitute the lyophilized content of each vial with 0.5 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

Note: *The reconstituted re-calibrators are stable for 2 months at 2 °C to 8 °C.*

Control 1 & 2

Reconstitute the lyophilized content of each vial with 0.5 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

Note: *The reconstituted controls are stable for 2 months at 2 °C to 8 °C.*

Wash Buffer (not included in the kit)

For Wash Buffer (1x) dilute 25 mL of Wash Buffer (40x) with 975 mL ultra-pure water to a final volume of 1000 mL.

The diluted Wash Buffer (1x) is stable for 2 weeks at room temperature.

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Whole blood (EDTA and heparin) can be used in this assay.

A minimum of 68 µL of pretreated sample is needed for one determination. This includes 8 µL sample and 60 µL dead volume. (See chapter 5.2.1 "Sample preparation".)

Attention:

- This test was not verified with blood collection tubes of all available manufacturers.
- Sample Collection Systems of some manufacturers may contain different materials which in isolated cases could affect the test results.
- If primary tubes for sample collection are used, please follow the instructions of the manufacturer.
- Do not use heat inactivated samples.
- Do not use standards or external controls stabilized with azide.

5.1 Specimen Collection

Whole blood should be collected into tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate anti-coagulant).

5.2 Specimen Storage and Preparation

Whole blood should be capped and may be stored for up to 7 days at 2 °C to 8 °C prior to performing the assay.

Whole blood stored for a longer time (up to two months) should be frozen only once at -20 °C prior to the assay. Thawed samples should be inverted several times prior to testing.

5.2.1 Sample Preparation

To determine HbA1c a hemolysate must be prepared for each sample.

1. Dispense **1000 µL** ultra-pure water into a tube.
2. Add **20 µL** of **well-mixed** whole blood and mix.
3. Allow to stand for 5 minutes or until complete lysis is evident.

Hemolysed specimens may be stored for up to 24 hours at 2 °C to 8 °C prior to performing the assay.

Note: Reconstituted Re-Calibrators and Controls need no further preparation.

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents, such as controls and re-calibrators and specimens, must be allowed to come to room temperature (20 °C to 25 °C) before use.
All reagents and samples must be mixed without foaming. Reagent Cartridges can be used directly after storage in the refrigerator.
- Samples, controls and re-calibrators should be measured within 2 hours in order to avoid possible evaporation effects.
- The *Secondary Sample Holder* (HYI-5437) for secondary tubes has the capacity for a maximum of 20 samples including controls and re-calibrators. They all have to be pipetted into the secondary tubes, and the respective barcodes of control/re-calibrator vials and, if available, the sample barcodes have to be read with the external barcode scanner.

6.2 Test Procedure

- **The total assay time for DRG:HYBRiD-XL HbA1c is 10 minutes.**
- To ensure proper operation of the test, the instructions in the user manual for the DRG:HYBRiD-XL should strictly be followed.
- All test specific information required for the correct operation is included in the respective barcodes of the reagents.
Take care not to damage these bar codes!
- **It is recommended to tap the bottom of the Cartridge Segments containing the reagent cartridges once on the bench before placing them on the rotor. This is to avoid foam and adhering of the liquid on the sealing of the reagent cartridge.**
- Place reagent cartridges on the rotor of the unit. The heating to 37 °C incubation temperature is performed automatically in the unit.

6.3 Calibration

Traceability:

This method was standardized against the following commercially available test material: *"Human Blood Hemolysate BCR-405, Institute for Reference Materials and Measurements (IRMM) Geel Belgium"*

Each DRG:HYBRiD-XL reagent contains a barcode with the specific information for recalibration of the reagent lot. The Master Curve is printed as a 2-D barcode on the outer label of the kit package and on the QC-Datasheet and has to be scanned with the external barcode scanner prior to the first use of the respective kit lot.

Recalibration is recommended:

- if a new kit lot is used. Each new lot should be verified by running the kit internal re-calibrators and controls before routine use.
- if one or both assay controls are found outside the specified range.
- after 4 weeks of use of the same reagent kit on the unit.

6.4 Calculation of Results

The analyte concentrations are calculated automatically by the DRG:HYBRiD-XL's system software.

7 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

It is also recommended to participate in national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analysing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following: expiration dates and storage conditions of reagents, operational reliability of the analyser. In addition, it is indicated to perform a Recalibration.

In case of further questions please contact your local distributor or DRG directly.

7.1 Internal Controls

For Quality Control it is necessary to use the two internal controls provided with each kit.

Acceptance ranges for both internal controls (*Control 1 & 2*) were established by the manufacturer and are summarized in the QC-Datasheet added to the kit. Note that the expected values and acceptance ranges stated in the QC-Datasheet always refer to the current kit lot.

Internal controls should be run in single determination:

- on a routine basis (e.g. once per 24 h)
- if re-calibration is required (if one or both internal controls are out of range)
- if a new kit lot is used (in order to avoid any negative impact on the kit performance by improper transport or to detect improper storage during transport).

7.2 External Controls

Use controls at both normal and pathological levels.

The control intervals and control ranges for external controls should be adapted to the individual requirements of each laboratory. All results must be within the defined limits.

Each laboratory should establish corrective measures to be taken if values of external controls are not found in the acceptance range.

8 EXPECTED NORMAL VALUES

Following values are recommended by the American Diabetes Association (ADA) (Diabetes Care 2012)

A non-diabetic person	< 39 mmol/mol
Diabetes	> 47 mmol/mol
Increased risk of developing diabetes in the future	39 - 46 mmol/mol

It is strongly recommended that each laboratory should determine its own normal and pathological values.

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated with other clinical observations and diagnostic tests.

9 LIMITATIONS OF USE

Reliable and reproducible results will be obtained, when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10 PERFORMANCE CHARACTERISTICS

10.1 Assay Dynamic Range

The dynamic range of the assay is defined by the minimum and the maximum value of the Master Curve.

Values found below the measuring range are indicated as "< 15 mmol/mol".

Values found above the measuring range are indicated as "> 150 mmol/mol".

The measuring range of the assay is between 15 mmol/mol – 150 mmol/mol.

10.2 Sensitivity

The sensitivity study was designed according to CLSI guideline EP17-A2.

The Limit of Blank (LoB) is 4.1 mmol/mol.

10.3 Precision Performance

The precision study was designed according to CLSI guideline EP5-A2.

10.3.1 Within-Run (Intra Assay)

Samples were analysed 1 times daily in duplicates for 20 days on 3 different instruments. CV (%) is given as mean of 40 duplicate CVs. The within-run variability is shown below:

Sample	n	Mean (mmol/mol)	CV (%)
1	40	28.3	2.9
2	40	31.3	2.2
3	40	42.6	1.8
4	40	108.7	1.6

10.3.2 Total Precision (Inter Assay)

Samples were analysed 2 times daily in duplicates for 20 days (n=80) on 3 different instruments. The total precision is shown below:

Sample	n	Mean (mmol/mol)	CV (%)
1	80	28.3	6.9
2	80	31.3	5.1
3	80	53.7	5.0
4	80	110.3	3.9

10.3.3 Inter-Lot

The between-lots variation was determined by 6 measurements of 4 samples with 3 different kit lots.

Sample	n	Mean (mmol/mol)	CV (%)
1	18	37.0	1.3
2	18	47.5	6.6
3	18	59.7	5.2
4	18	119.2	2.3

10.4 Interfering Substances

Haemoglobin (up to 5 mg/mL), Bilirubin (up to 0.3 mg/mL) and Triglyceride (up to 7 mg/mL) have no influence on the assay results. Until today no substances (drugs) are known to us, which have an influence to the measurement of HbA1c in a sample.


10.5 High-Dose-Hook Effect

Hook effect was not observed in this test up to a concentration of 166 mmol/mol of HbA1c.

11 METHOD COMPARISON

A comparison of DRG:HYBRID-XL HbA1c Test HYC-5325 (y) and a commercially available reference method (SEBIA, Capillary Electrophoresis) (x) using clinical samples gave the following correlation:

$$\begin{aligned} n &= 50 \\ r &= 0.946 \\ y &= 1.3978x - 28.67 \end{aligned}$$

12 LEGAL ASPECTS

Only for countries where the declaration of European Conformity (CE mark) is applicable.

12.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different kit lots could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

For further information please refer to the User Manual of the DRG:HYBRiD-XL, analyser-specific application sheets, product information and package inserts of all necessary components.

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

13 REFERENCES / LITERATURE / LITERATUR / BIBLIOGRAFIA / LITERATURA

1. Clinical Laboratory Diagnostics/ Lothar Thomas. 8th ed. Frankfurt: TH-Books Verlagsgesellschaft; 2012. 236 ff.
2. Tietz Fundamentals of Clinical Chemistry/Carl A. Burtis, Edward R. Ashwood - 5th ed. Philadelphia: W.B. Saunders Company; 2001. 427 ff.
3. Henry's Clinical Diagnosis and Management. 22nd ed. Philadelphia: W.B. Saunders Company; 2011. 510 ff.
4. Mazurek et al.; Prevalence of HbA1c greater than 6.5% and 7.0% among hospitalized patient without known diagnosis of diabetes in an urban inner city hospital. J Clin Endocrinol Metab 2010; 95, 1344-8.
5. American Diabetes Association. Standards of medical care in diabetes. Diabetes Care 2012, 35, 11-63.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitätskennzeichnung	Conformità europea	Conformidad europea	Conforme aux normes européennes
	Consult instructions for use*	Gebrauchsanweisung beachten*	Consultare le istruzioni per l'uso	Consulte las Instrucciones	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device*	<i>In-vitro</i> -Diagnostikum*	Per uso Diagnostica in vitro	Para uso Diagnóstico in vitro	Usage Diagnostic in vitro
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
	Catalogue number*	Artikelnummer*	No. di Cat.	Número de catálogo	Référence
	Batch code*	Chargencode*	Lotto no	Número de lote	No. de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit*	Temperaturbegrenzung*	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis*	Data di scadenza	Fecha de caducidad	Date limite d'utilisation
	Manufacturer*	Hersteller*	Fabbricante	Fabricante	Fabricant
	Caution*	Achtung			
<i>Distributed by</i>	Distributor	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Numéro
<i>Reagent Cartridge</i>	Reagent Cartridge	Reagenzien-Cartridge	Cartucce di reagente	Cartucho de reactivo	Cartouche de réactif
<i>Re-Calibrator</i>	Re-Calibrator	Re-Kalibrator	Re-Calibratore	Re-Calibrador	Re-calibrateurs
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Tracciante enzimatico	Conjugado enzimático	Conjugué enzymatique
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complesso enzimatico	Complex enzimático	
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungsmedium	Diluyente dei campioni	Solución para dilución de la muestra	Diluant d'échantillon
<i>Reagent 1</i>	Reagent 1	Reagenz 1	Reagente 1	Reactivo 1	Réactif 1
<i>Reagent 2</i>	Reagent 2	Reagenz 2	Reagente 2	Reactivo 2	Réactif 2
<i>Wash Buffer</i>	Wash Buffer	Waschpuffer	Tampone di lavaggio	Tampón de lavado	Tampon de lavage
<i>System Solution</i>	System Solution	Systemlösung		Solución de sistema	Solution système
<i>Needle Cleaning Solution</i>	Needle Cleaning Solution	Nadel-Reinigungslösung	Soluzione lavaggio ago	Solución de lavado de la aguja	Solution de nettoyage des aiguilles
<i>Denaturation Buffer</i>	Denaturation Buffer	Denaturierungspuffer	Tampone die denaturazione	Tampón de denaturalización	Tampon de dénaturation
<i>Neutralization Buffer</i>	Neutralization Buffer	Neutralisierungspuffer	Tampone di neutralizzazione	Tampón de neutralización	Tampon de neutralisation
<i>Secondary Sample Holder</i>	Secondary Sample Holder	Sekundärprobenhalter	Sostegno secondario dei campioni	Soporte de muestras para tubos secundarios	Support d'échantillons secondaires
<i>Sample Tubes (Secondary)</i>	Sample Tubes (Secondary)	Sekundärprobenröhrchen	Tubetti di campioni secondari	Tubos de muestra secundarios	Tubes échantillon secondaires
<i>Dilution Cartridge</i>	Dilution Cartridge				Cartouches de dilution
<i>Vial Adapter</i>	Vial Adapter	Fläschchen-Adapter	l'adattatore dei tubetti	Adaptador de tubos	l'Adaptateur de flacons