DRG:HYBRID-XL®

Instructions for Use

HbA1c

IVD

CE

REF

HYC-5325

80 tests

DRG

DRG Instruments GmbH, Germany

Frauenbergstraße 18, 35039 Marburg

Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50

Website: www.drg-diagnostics.de

E-mail: drg@drg-diagnostics.de

DRG

DRG International, Inc., USA

841 Mountain Ave., Springfield, NJ 07081

Phone: (973) 564-7555, Fax: (973) 564-7556

Website: www.drg-international.com

E-mail: corp@drg-international.com

ASSAY PROTOCOL BARCODE (APP)

ASSAYPROTOKOLL-BARCODE / CODICE A BARRE DEL PROTOCOLLO DEL DOSAGGIO / CÓDIGO DE BARRAS DEL PROTOCOLO DEL ENSAYO / CODE-BARRES DU PROTOCOLE DE TEST

(Software Version 2.61 or later / ab Softwareversion 2.61 oder höher / Versione software 2.61 o successiva / Versión 2.61 del software o posterior / Version de logiciel 2.61 ou ultérieure / Wersja oprogramowania 2.61 lub późniejsza)

HYC-5325 – v2.61.1

The barcode must be used to install the assay protocol into the DRG:HYBRiD-XL software via the menu “SCAN NEW LOT”.

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Please refer to section 6.2.6 “Assay – Scanning a new kit” of the User Manual v2.61 or later.

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Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas / Modifications apportées / Modificações introduzidas / Wprowadzone modyfikacje											
The following changes have been made in comparison to the previous version: Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen: Rispetto alla versione precedente, sono state apportate le seguenti modifiche: Se han introducido los siguientes cambios en comparación con la versión anterior: Les modifications suivantes ont été apportées par rapport à la version précédente : Foram efetuadas as seguintes alterações em comparação com a versão anterior: W porównaniu z poprzednią wersją wprowadzono następujące zmiany:											
In general, the IFU is updated according to the requirements of the IVD Regulation (EU) 2017/746. Generell wurde die IFU gemäß den Anforderungen der IVD-Verordnung (EU) 2017/746 aktualisiert. In generale, le Istruzioni per l'uso sono state aggiornate secondo i requisiti del Regolamento IVD (UE) 2017/746. En términos generales, las IFU se han actualizado conforme a los requisitos del reglamento sobre productos sanitarios para diagnóstico in vitro (UE) 2017/746. De façon générale, mise à jour des instructions d'utilisation conformément aux exigences du Règlement (UE) 2017/746 relatif aux DIV. De um modo geral, as Instruções de Utilização (IFU) estão atualizadas de acordo com os requisitos do regulamento de IVD (UE) 2017/746.											
<table border="1"> <tbody> <tr> <td>1 INTENDED USE</td> <td>New wording</td> </tr> <tr> <td>2 PRINCIPLE OF THE TEST</td> <td>New wording</td> </tr> <tr> <td>4.1 Materials provided with the Kit</td> <td>Re-Calibrator 1 & 2: new ref. material</td> </tr> <tr> <td>5.2 Sample preparation</td> <td>New volumes (now 500µL water + 10 µL whole blood [old: 1000+20µL])</td> </tr> <tr> <td>11 METHOD COMPARISON</td> <td>New data</td> </tr> </tbody> </table>		1 INTENDED USE	New wording	2 PRINCIPLE OF THE TEST	New wording	4.1 Materials provided with the Kit	Re-Calibrator 1 & 2: new ref. material	5.2 Sample preparation	New volumes (now 500µL water + 10 µL whole blood [old: 1000+20µL])	11 METHOD COMPARISON	New data
1 INTENDED USE	New wording										
2 PRINCIPLE OF THE TEST	New wording										
4.1 Materials provided with the Kit	Re-Calibrator 1 & 2: new ref. material										
5.2 Sample preparation	New volumes (now 500µL water + 10 µL whole blood [old: 1000+20µL])										
11 METHOD COMPARISON	New data										

ENGLISH.....	4
DEUTSCH	9
ITALIANO	14
ESPAÑOL.....	19
FRANÇAIS.....	24
REFERENCES / LITERATURE	29
SYMBOLS USED	30

1 INTENDED USE

The DRG:HYBRID-XL HbA1c is an **automated** immunoturbidimetric assay for the **quantitative measurement** of Hemoglobin A1c (HbA1c) in human whole blood. The reported HbA1c value is indicated as the ratio of the glycosylated hemoglobin (Hemoglobin A1c), against the total hemoglobin concentration in the sample.

For *in vitro* diagnostic use.

For laboratory professional use, only with the DRG:HYBRID-XL Analyzer.

The device is **intended to be used** for the monitoring of long-term blood glucose control in patients previously diagnosed with diabetes mellitus.

The device is **not intended for** screening or diagnosis of diabetes or neonatal use.

1.1 Scientific Validity Report

HbA1c is a glycated form of hemoglobin formed by the attachment of glucose residues in the blood to the hemoglobin molecules. In the diabetic population where blood glucose levels are abnormally elevated the level of HbA1c also increases.

The level of HbA1c is proportional to the level of glucose in the blood and has been widely accepted as an indicator of the mean blood glucose concentration in the preceding 6-8 weeks. It is therefore a long-term indicator of diabetic control. For routine use HbA1c levels should be monitored every 3-4 months. However, in gestational diabetes and after a change in therapy it may be useful to measure HbA1c more frequently at 2-4 week intervals.

2 PRINCIPLE OF THE TEST

The DRG:HYBRID-XL HbA1c Kit is an **immunoturbidimetric assay**.

The sample (hemolysis sample) is added to the unsensitised latex particles, and the surface of the latex adsorbs total hemoglobin in the sample.

Monoclonal (mouse) anti-human HbA1c antibodies interact with the latex bound HbA1c. Reaction with (goat) anti-mouse IgG antibody leads to an amplification of this reaction and the forming of agglutination complexes.

The amount of agglutination caused depends on the amount of HbA1c that adsorbs the surface of the latex, this agglutination results in an increase in turbidity and the absorbance is measured at 645 nm. The change in absorbance signal is proportional to the amount of HbA1c in the sample. HbA1c concentration (mmol/mol) in the sample is determined by comparison with a standard curve.

3 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only.
For laboratory professional use only.
- This kit can only be used in combination with the DRG:HYBRID-XL Analyzer.
- Before starting the assay, read the instructions for use completely and carefully. Use the valid version of the instructions for use provided with the kit. Be sure that everything is understood.
- Do not remove, exchange, discard, or damage any of the barcode labels provided with each kit and its components.
All barcodes build an integral system for the kit lot.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange reagent cartridges of different kits even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the wells in the reagent cartridges may differ slightly.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Reagents of other manufacturers must not be used together with the reagents of this test kit.
- All reagents in this kit are clear liquids, substrate solution is clear and colorless. Changes in its appearance may affect the performance of the test. In that case, contact DRG.
- Optimal test results are only obtained when using calibrated pipettes.
- Microbial contamination of reagents or samples may give false results.

General precautions

- Follow laboratory quality assurance and laboratory safety guidelines.

- Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
- Do not smoke, eat, drink, or apply cosmetics in areas where samples or kit reagents are handled.
- Wear lab coats and disposable latex gloves when handling samples and reagents and where necessary safety glasses.

Biohazard information

- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. However, no known test method can offer total assurance that no infectious agent is present.
- The device contains material of animal origin, which is certified apparently free of infectious or contagious diseases and injurious parasites.
- Bovine components originate from countries where BSE (Bovine spongiform encephalopathy) has not been reported.
- All materials and samples of human or animal origin must be handled as if capable of transmitting infectious diseases.
- Handling must be done in accordance with the procedures defined by appropriate national biohazard and safety guideline or regulation. Waste must be discarded according to local rules and regulations.

Information to chemical hazards and hazard classification

- Some reagents contain preservatives in non-declarable concentrations. Nevertheless, in case of contact with eyes or skin, flush immediately with water.
- Substrate Solution contains an ingredient in non-declarable concentrations which causes serious eye irritation. In case of possible contact with eyes, rinse immediately carefully and thoroughly with eye wash or water. After contact with skin, wash with plenty of water. Take-off contaminated clothing and wash it before reuse
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national safety guidelines or regulations.
- This product does not contain substances which have carcinogenic, mutagenic, or toxic for reproduction (CMR) properties.

All reagents of this test kit do NOT contain hazardous substances in concentrations to be declared, a classification and labelling is not required.

4 MATERIALS

4.1 Materials provided with the kit

Symbol	Quantity/ Preparation	Description
Reagent Cartridges	4 x 20 pieces in a resealable plastic bag with desiccant Ready to use	Reagent Cartridge containing the following: – Reagent 1, 260 µL , Suspension of latex particles; – Reagent 2, 130 µL Buffer solution containing anti-human HbA1c antibody (monoclonal, mouse) and anti-mouse IgG (goat) antibody.
Re-Calibrator R1 & R2	2 x 0.5 mL Lyophilized See "Reagent Preparation".	Recalibrators * For recalibration of the quantitative DRG:HYBRID-XL HbA1c test. <u>Concentrations are lot-specific.</u> Conversion: NGSP (%) = DCCT (%) = $[0.09148 * \text{IFCC (mmol/mol)}] + 2.152$ Standardized against the following reference material: ERM®- AD500/IFCC
Control C1 & C2	2 x 0.5 mL Lyophilized; See "Reagent Preparation".	Controls * For control values and ranges please refer to the vial label or to the CoA.
	1 x	Instructions for Use
	1 x	Certificate of Analysis (CoA)

* Contain non-mercury preservative.

4.2 Materials required but not provided

- General needed laboratory equipment
- Ultra-pure water
DRG recommends to use Clinical Laboratory Reagent Water (CLRW) according to CLSI guideline 3C-A4 with the following specifications:
Resistivity at 25 °C [MΩ·cm]: > 10
Conductivity at 25 °C [µS·cm⁻¹]: < 0.1
Total Organic Carbon/p.p.b.[µg/L] : < 50
Colloids [µg/mL]: <0.05
- **REF** HYB-5252 DRG:HYBRID-XL Analyzer
- **REF** HYI-5392: System Solution 5L, 5000 mL;
(Instrument Feed Water according to CLSI guideline 3C-A4 with the following specification can also be used:
Resistivity at 25 °C [MΩ·cm]: > 1
Conductivity at 25 °C [µS·cm⁻¹]: < 1
Total Organic Carbon/p.p.b.[µg/L] : < 200
Colloids [µg/mL]: <0.1)
- **REF** HYI-5394: Wash Buffer, 40x concentrate, 25 mL
- **REF** HYI-5395: Needle Cleaning Solution, 30 mL.
Cleaning solution for the pipetting needle (daily and weekly maintenance, see also user manual)
- **REF** HYI-5387: Cuvettes, (2 x 360 pieces)

For use of the Secondary Sample Holder for secondary tubes the following tubes are required:

- **REF** HYI-5391: Sample Tubes (Secondary), 2500 pieces

4.3 Storage and Stability of the Kit

All kit components must be handled and stored as described below.

Unopened and unused components or kits

	Storage Temperature	Stability
– Unopened kits – Unopened reagents – Cartridges in the supplied <u>unopened</u> plastic bag – Unused reagent cartridges	2 °C to 8 °C	Until the expiration date printed on the label. Do not use reagents beyond this date!

Unused cartridges must be stored in the blister pack in the sealed foil pouch including the desiccant. Once the foil bag has been opened, care must be taken to close it tightly again.

Opened components

	Storage Temperature	Stability
Opened – Re-Calibrators – Controls	2 °C to 8 °C	(For reconstituted reagents refer to "4.4 Reagent Preparation".)

Immediately after end of each run the Re-Calibrator and Control vials have to be removed from the instrument, tightly capped, and stored at 2 °C to 8 °C.

Pierced or open cartridges must be disposed of immediately.

4.3.1 On-board Stability

For Re-Calibrators and Controls the on-board stability has been evaluated under controlled laboratory conditions at room temperature (20 °C to 26 °C).

Due to the differences in laboratory environmental conditions and reagent volumes, the on-board stability may deviate from the declared value.

On-board stability:	48 h
---------------------	------

4.4 Reagent Preparation

Bring all reagents, such as controls and re-calibrators, to room temperature (20 °C to 26 °C) prior to use. Reagent Cartridges can be used directly after storage in the refrigerator.

Re-Calibrator 1 & 2

Reconstitute the lyophilized content of each vial with 0.5 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

Stability after reconstitution:	at 2 °C to 8 °C	8 weeks
	at -20 °C (in aliquots)	8 weeks

Control 1 & 2

Reconstitute the lyophilized content of each vial with 0.5 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

Stability after reconstitution:	at 2 °C to 8 °C	8 weeks
	at -20 °C (in aliquots)	8 weeks

Wash Buffer (not included in the kit)

For Wash Buffer (1x) dilute 25 mL of Wash Buffer (40x) with 975 mL ultra-pure water to a final volume of 1000 mL.

Stability after dilution:	at 20 °C to 26 °C	1 week
---------------------------	-------------------	--------

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SAMPLE COLLECTION, STORAGE AND PREPARATION

The following sample material can be used in this test:

Whole blood (K3 EDTA or lithium heparin).

A minimum volume of 68 µL of sample is needed for one determination. This includes 8 µL sample and 60 µL dead volume.

Attention:

- This test was not verified with blood collection tubes of all available manufacturers.
- Sample collection systems of some manufacturers may contain different materials which in isolated cases could affect the test results.
- If primary tubes for sample collection are used, please follow the instructions of the manufacturer.
- In general, it should be avoided to use hemolytic, icteric, or lipemic samples. For further information refer to chapter "*Interfering Substances*".
- Samples containing precipitates have to be centrifuged prior to the test run.
- Do not use heat inactivated samples.
- Samples or external controls containing sodium azide should not be used in the assay.

5.1 Sample Collection

Whole blood: Whole blood should be collected in centrifuge tubes containing anticoagulant (K₃-EDTA or Li-Heparin) (e.g. Sarstedt Monovette with the appropriate plasma preparation) and not centrifuged after collection.

5.2 Sample Storage

Samples must be stored tightly capped prior to performing the assay. If stored frozen, freeze only once. Thawed samples must be inverted several times prior to testing.

Stability	at 2 °C to 8 °C	7 days
	at -20 °C (in aliquots)	up to 12 months

5.3 Sample Preparation

To determine HbA1c, a hemolysate must be prepared for each sample.

1. Dispense **500 µL** ultra-pure water in a tube.
2. Add **10 µL** of **well-mixed** whole blood and mix.
3. Allow to stand for 5 minutes or until complete lysis is evident.

Note: Reconstituted Re-Calibrators and Controls need no further preparation.

6 ASSAY PROCEDURE

6.1 Procedural Notes

- All reagents, such as controls and re-calibrators, and samples, must be allowed to come to room temperature (20 °C to 26 °C) before use. All reagents must be mixed without foaming. Reagent Cartridges can be used directly from the cold storage (refrigerator).
- Avoid bubbles/foam in the calibrator and control vials, especially when using them directly in the *Primary Sample Holder* (HYI-5436). Bubbles/foam can lead to wrong liquid level detection and hence false results.
- Do not interchange caps of reagent vials to avoid cross-contamination.
- To avoid possible evaporation effects of the recalibrators, controls and samples placed on the instrument, the chapter "*On-Board Stability*" must be considered.
- For primary samples tubes the *Primary Sample Holder* with a capacity for 12 samples can be used. For low volume samples the *Secondary Sample Holder* (HYI-5437) can be used. It has the capacity for a maximum of 20 samples for secondary tubes. The samples must be transferred manually into the secondary tubes (*Sample Tubes (Secondary)*).

6.2 Test Procedure

Total assay time: 10 minutes

Detailed information about operating the DRG:HYBRID-XL analyzer must be taken from the system's User Manual.

Test specific information required for automatic processing is contained in distinct barcodes of the kit:

- **Assay Protocol Barcode (APB):**
2D Data Matrix code containing the specific assay protocol/method. The APB is found on page 2 of this IFU.
- **Lot barcode:**
2D Data Matrix code containing lot specific information including the master curve (factory calibration). The lot barcode is found on the label of the kit package and CoA.
- **Re-calibrator and control vial barcodes:**
2D Data Matrix code to identify the re-calibrators and controls. This code is found on the vial labels and CoA.
- **Cartridge barcode:**
Reagent cartridges are sealed with a 1D barcode for cartridge identification. Take care not to damage this barcode when inserting the cartridges.
 - If it is the first use of this **test** on the **analyzer** the **Assay Protocol Barcode (APB)** on page 2 of this instructions for use must be loaded with the external barcode scanner and stored.
The APB includes all information about the test procedure itself (specific protocol data, including volumes, incubation times, units etc.).
Once installed it will remain until a new APB will become available by the manufacturer.
 - If the **lot** is used the first time or a new lot of this test is available, the lot barcode on the outer label of the kit package or on the CoA must be loaded with the external barcode scanner and stored.
 - Program the number of required samples and continue with the test assignment.
 - Insert the number of required cartridges and cuvettes with care in the segment. **Make sure cartridges and cuvettes are inserted correctly.**
 - Load the segments by placing them on the reagent rotor. Make sure the inserted segment is seated correctly (not tilted!).
 - Continue to start the run.

6.3 Recalibration

Detailed information how to perform the recalibration is found in the system's User Manual.

Recalibration is recommended:

- if a new kit lot is used. Each new lot should be verified by running the kit internal re-calibrators and controls before routine use.
- if one or both assay controls are found outside the specified range.
- after 4 weeks of use of the same reagent kit on the unit.

Always both recalibrators (R1 & R2) are required for recalibration.

Always two internal controls are required to validate the recalibration.

Recalibrators and controls are run together in the same run. In addition, samples can be included.

In the "*Load-Samples*" menu scan the respective barcodes for re-calibrators and controls found on the vials or CoA.

Test assignment for re-calibrators and controls is done automatically.

Continue with loading the cartridges and cuvettes required.

Start the run.

6.4 Calculation of Results

The analyte concentrations are calculated automatically by the DRG:HYBRID-XL's system software.

For the result calculation a Cubic Spline Fit is used.

7 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results.

If available, it is also recommended to participate in national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analyzing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following: expiration dates and storage conditions of reagents, operational reliability of the analyzer. In addition, it is indicated to perform a Recalibration.

After checking the above-mentioned items without finding any error contact your distributor or DRG directly.

7.1 Internal Controls

For Quality Control it is necessary to use the two internal controls provided with each kit.

Acceptance ranges for both internal controls (*Control 1 & 2*) were established by the manufacturer and are summarized in the Certificate of Analysis added to the kit. Note that the expected values and acceptance ranges stated in the Certificate of Analysis always refer to the current kit lot.

Internal controls should be run in single determination:

- on a routine basis (e.g. once per 24 h)
- if recalibration is required (if one or both internal controls are out of range)
- if a new kit lot is used (in order to avoid any negative impact on the kit performance by improper transport or to detect improper storage during transport).

7.2 External Controls

Use controls at both normal and pathological levels.

The control intervals and control ranges for external controls should be adapted to the individual requirements of each laboratory. All results must be within the defined limits.

Each laboratory should establish corrective measures to be taken if values of external controls are not found in the acceptance range.

8 REFERENCE VALUES

It is strongly recommended that each laboratory should determine its own reference values.

Following values are recommended by the American Diabetes Association (ADA) (Diabetes Care 2012):

A non-diabetic person	< 39.00 mmol/mol
Diabetes	> 47.00 mmol/mol
Increased risk of developing diabetes in the future	39.00 – 46.00 mmol/mol

Values above or below the reference range should be considered as suspicious and require additional testing.

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated with other clinical observations and diagnostic tests.

9 LIMITATIONS OF THE PROCEDURE

Reliable and reproducible results will be obtained, when the assay procedure is performed with a complete understanding of the instructions for use and in compliance with the laboratory quality assurance guidelines.

Any improper handling of samples or modification of this test might influence the results.

10 PERFORMANCE CHARACTERISTICS

10.1 Sensitivity

The sensitivity study was designed according to CLSI guideline EP17-A2. Due to the unique Test properties, it is not possible to determine the LoD and LoQ by diluting samples.

Limit of Blank (LoB)	4.140 mmol/mL
Measuring range	4.14 – 150.00 mmol/mol
Linear range	15.00 – 150.00 mmol/mol
	Indicated as:
Values below the measuring range	< 15.00 mmol/mL
Values above the measuring range	> 150.00 mmol/mL

10.2 Specificity of Antibodies (Cross-Reactivity)

No relevant substances (should be related structurally or biochemically to the Assay antigen) are known to us.

10.3 Precision Performance

The precision study was designed on the basis of to CLSI guideline EP5-A3.

10.3.1 Within-run precision

Calculated from the mean value of 40 replicates with 1 device.

Sample	n	Mean (mmol/mL)	CV (%)
1	40	28.30	2.4
2	40	42.57	1.8
3	40	53.72	2.6
4	40	110.34	2.6

10.3.2 Between-run precision

Calculated by 80 replicates each with 1 device.

Sample	n	Mean (mmol/mL)	CV (%)
1	80	28.30	6.9
2	80	42.57	7.9
3	80	53.72	5.0
4	80	110.34	3.9

10.3.3 Between-lot precision

Determined with 3 different kit lots.

Sample	n	Mean (mmol/mL)	CV (%)
1	18	37.03	1.3
2	18	59.75	5.2
3	18	47.53	6.6
4	18	119.23	2.3

10.4 Recovery

Due to the unique Test properties, it is not possible to determine the Recovery.

10.5 Linearity

Due to the unique Test properties, it is not possible to determine the Linearity.

10.6 Interfering Substances

Hemoglobin (up to 5.0 mg/mL), bilirubin (up to 0.3 mg/mL) and triglyceride (up to 7.0 mg/mL) have no influence on the assay results.

10.7 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence on the measurement of Hemoglobin A1c (HbA1c) in a sample.


10.8 High-Dose Hook Effect

Due to the unique Test properties, it is not possible to determine the High Dose Hook Effect.

11 METHOD COMPARISON

A comparison of DRG:HYBRID-XL HbA1c Test (HYC-5325) (y) and the Roche Cobas c513 - Tina-quant® Hemoglobin A1cDx Gen.3 (x) using clinical samples gave the following correlation:

$$\begin{aligned} n &= 40 \\ r &= 0.97 \\ y &= 0.9182x + 2.128 \end{aligned}$$

12 LEGAL ASPECTS

12.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the laboratory quality assurance guidelines and applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. If there is any doubt or concern regarding a result, please contact DRG.

12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different kit lots could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12.4 Reporting of Serious Incident

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

For further information please refer to the User Manual of the DRG:HYBRID-XL, analyzer-specific application sheets, product information and package inserts of all necessary components.

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12.4 Notification des incidents graves

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1. Clinical Laboratory Diagnostics/ Lothar Thomas. 8th ed. Frankfurt: TH-Books Verlagsgesellschaft; 2012. 236 ff.
2. Hempe JE and Hsia DS. Variation in the hemoglobin glycation index. Journal of Diabetes and its Complications 2022; 36(7), 108223.
3. Khor XY et al. Individualized diabetes care: Lessons from the real-world experience. 2023; 11(13), 2890-2902.
4. Vlachos G et al. Glycemic Index (GI) or Glycemic Load (GL) and Dietary Interventions for Optimizing Postprandial Hyperglycemia in Patients with T2 Diabetes: A Review Nutrients; 2020; 12, 1561.
5. Mazurek et al.; Prevalence of HbA1c greater than 6.5 % and 7.0 % among hospitalized patient without known diagnosis of diabetes in an urban inner city hospital. J Clin Endocrinol Metab 2010; 95, 1344-8.
6. American Diabetes Association. Standards of medical care in diabetes. Diabetes Care 2012, 35, 11-63.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français	Português
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes	Conformidade Europeia
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation	Consultar as instruções de uso
	In vitro diagnostic medical device *	In-vitro-Diagnostikum *	Diagnostica in vitro	Diagnóstico in vitro	Diagnostic in vitro	Dispositivo médico para diagnóstico in vitro
	Catalogue number *	Katalognummer *	No. di Cat.	No de catálogo	Référence	Número de catálogo
	Batch code *	Chargen-bezeichnung *	Lotto no	Número de lote	No. de lot	Código do lote
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests	Suficiente para <n> determinações
	Temperature limit *	Temperaturgrenzwerte *	Temperatura di conservazione	Temperatura de conservacion	Temperature de conservation	Limites de temperatura
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation	Prazo de validade
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant	Fabricante
	Distributor *	Vertriebspartner *	Distributore	Distribuidor	Distributeur	Distribuidor
	Date of manufacture *	Herstellungsdatum *	Data di produzione	Fecha de fabricación	Date de production	Data de fabricação
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques	Riscos biológicos
	Caution *	Achtung *	Attenzione	Precaución	Attention	Cuidado
	Unique device Identifier *	eindeutige Produktidentifizierung *				
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches	
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement	Conteúdo
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité	Volume / Quantidade
<i>Reagent Cartridge</i>	Reagent Cartridge	Reagenzien-Cartridge	Cartucce di reagenti	Cartuchos de reactivos	Cartouche de réactif	
<i>Re-Calibrator</i>	Re-Calibrator	Re-Kalibrator	Re-Calibratore	Re-Calibrador	Re-calibrateurs	
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle	Controlo
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungs-medium	Diluente per campioni	Diluyente de muestras	Diluant d'échantillon	Diluente de Amostra
<i>Reagent 1</i>	Reagent 1	Reagenz 1	Reagente 1	Reactivo 1	Réactif 1	
<i>Reagent 2</i>	Reagent 2	Reagenz 2	Reagente 2	Reactivo 2	Réactif 2	
<i>Dilution Cartridge</i>	Dilution cartridge	Verdünnungscartridge	Cartucce di diluizione	Cartuchos de dilución	Cartouches de dilution	
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Coniugato enzimatico	Conjugado enzimático	Conjugué enzymatique	Koniugat enzymatyczny
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complesso enzimatico	Complex enzimático	Complexe enzymatique	Kompleks Enzymu
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat	Roztwór Substratu
<i>Wash Buffer</i>	Wash Buffer	Waschpuffer	Tampone di lavaggio	Tampón de lavado	Tampon de lavage	
<i>System Solution</i>	System Solution	Systemlösung	Soluzione di sistema	Solución de sistema	Solution système	
<i>Needle Cleaning Solution</i>	Needle Cleaning Solution	Nadel-Reinigungslösung	Soluzione lavaggio ago	Solución de limpieza de la aguja	Solution de nettoyage des aiguilles	
<i>Decontamination Solution</i>	Decontamination Solution	Dekontaminierungslösung				
<i>Neutralization Buffer</i>	Neutralization Buffer	Neutralisierungspuffer	Tampone di neutralizzazione	Tampón de neutralización	Tampon de neutralisation	
<i>Denaturation Buffer</i>	Denaturation buffer	Denaturierungspuffer	Tampone di denaturazione	Tampón de denaturalización	Tampon de dénaturation	
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampone del test	Tampón de ensayo	Tampon d'essai	

Symbol	English	Deutsch	Italiano	Español	Français	Português
<i>Conjugate Buffer</i>	Conjugate Buffer	Konjugatpuffer	Tampone del coniugato	Tampón de conjugado	Tampon de conjugué	
<i>Secondary Sample Holder</i>	Secondary sample holder	Sekundärprobenhalter	Portacampioni secondario	Soporte de muestras para tubos secundarios	Support d'échantillons secondaires	
<i>Sample Tubes (Secondary)</i>	Sample Tubes (secondary)	Sekundärprobenröhrchen	Provette campione (secondarie)	Tubos de muestra (secundarios)	Tubes échantillon (secondaires)	
<i>Vial Adapter</i>	Vial Adapter	Fläschchen-Adapter	L'adattore dei tubetti	Adaptor de tubos	L'adaptateur de flacons	