DRG:HYBRID•XL[®]

Instructions for Use

CRP-HS

IVD

REF

HYC-5613

80 tests

DRG Instruments GmbH, Germany
Frauenbergstraße. 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

ASSAY PROTOCOL BARCODE (APB)

**(Version 2.60 Software or later / ab Softwareversion 2.60 oder höher / Versione Software 2.60 o superiore /
Versión de Software 2.60 o posterior / Version de logiciel 2.60 et supérieure /
Wersja oprogramowania 2.6 lub późniejsza)**

HYC-5613– v2.60.4

The barcode must be used to install the assay protocol into the DRG:HYBRiD-XL software via the SCAN NEW LOT page.
Der Barcode muss in dem Menü „NEUES LOT SCANNEN“ eingelesen werden, um das Protokoll in der DRG:HYBRiD-XL-Software zu installieren.

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Please refer to section **3: Routine Procedures: "Installing a new assay product"** of the User Manual v2.60 or later.

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Contents / Inhaltsverzeichnis / Contenuti / Contenido / Contenu / Spis treści

| | | |
|----|---|---|
| 1 | INTRODUCTION | 4 |
| 2 | PRINCIPLE OF THE TEST | 4 |
| 3 | WARNINGS AND PRECAUTIONS | 4 |
| 4 | REAGENTS | 4 |
| 5 | SPECIMEN COLLECTION AND PREPARATION | 5 |
| 6 | ASSAY PROCEDURE | 5 |
| 7 | QUALITY CONTROL | 5 |
| 8 | EXPECTED NORMAL VALUES | 6 |
| 9 | LIMITATIONS OF USE | 6 |
| 10 | PERFORMANCE CHARACTERISTICS | 6 |
| 11 | METHOD COMPARISON | 7 |
| 12 | LEGAL ASPECTS | 7 |

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1 INTRODUCTION

1.1 Intended Use

The DRG:HYBRID-XL CRP-HS is an immunoturbidimetric assay for the quantitative *in vitro diagnostic* measurement of CRP in serum and plasma (EDTA-, heparin- or citrate plasma). Only for use with the DRG:HYBRID-XL Analyzer.

1.2 Summary and Explanation

CRP (C-reactive protein) is an acute phase protein whose concentration is seen to increase as a result of the inflammatory process, most notably in response to pneumococcal (bacterial) infectious, histolytic disease and a variety of disease states. Originally discovered by Tillett et al. in 1930 inpatient sera with acute infection, CRP has now come to be used as a marker or general diagnostic indicator of infections and inflammation, in addition to serving as a monitor of patient response to therapy and surgery. Furthermore, regular measurements of CRP in infants can be a useful aid in the early diagnosis of infectious disease.

2 PRINCIPLE OF THE TEST

The DRG:HYBRID-XL CRP-HS Kit is an **immunoturbidimetric assay**.

Sample is reacted with a buffer and anti-CRP coated latex.

The formation of the antibody-antigen complex during the reaction results in an increase in turbidity and the absorbance is measured at 570 nm. The change in absorbance signal is proportional to the amount of CRP in the sample. CRP concentration in the sample is determined by comparison with a standard curve.

3 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only. For professional use only.
- This kit can only be used in combination with the DRG:HYBRID-XL Analyzer
- Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is clear and understood.
- Do not remove, exchange, discard or damage any of the barcode labels provided with each kit and its components. All barcodes build an integral system for the kit lot.**
- Respect the general safety measures for use of laboratory reagents.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear appropriate disposable gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may cause false results.
- Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Unused reagent cartridges must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant provided.
- Optimal test results are only obtained when using calibrated pipettes.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange reagent cartridges of different kits even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the wells in the reagent cartridges may differ slightly.
- Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.
- For information on hazardous substances included in the kit please refer to Safety Data Sheets. For professional users the Safety Data Sheet for this product is available upon request directly from DRG.

4 REAGENTS

4.1 Reagents provided

4.1.1 Reagent Cartridges

80 pieces containing the following:

- Reagent 1**, 220 µL
Buffer solution;
Contains non-mercury preservative.
- Reagent 2**, 220 µL
Suspension of latex particles coated with anti-human CRP antibody (monoclonal);
Contains non-mercury preservative.

4.1.2 Re-Calibrator 1 & 2

2 vials, 0.5 mL each, ready to use;

For re-calibration of the quantitative DRG:HYBRID-XL CRP-HS test.

Concentrations are lot-specific.

Re-Calibrators are standardised against the following reference material: WHO International Standard Human C-reactive Protein (Code 85/506)

Contain non-mercury preservative.

4.1.3 Control 1 & 2

2 vials, 0.5 mL each, ready to use;

For control values and ranges please refer to the bar code on vial label or to the QC-Datasheet.

Contain non-mercury preservative.

4.2 Materials required but not provided

- General needed laboratory equipment
- Ultra-pure water
DRG recommends to use Clinical Laboratory Reagent Water (CLRW) according to CLSI guideline 3C-A4 with the following specifications:
Resistivity at 25 °C [MΩ·cm]: > 10
Conductivity at 25 °C [µS·cm⁻¹]: < 0.1
Total Organic Carbon/p.p.b.[µg/L]: < 50
Colloids [µg/mL]: <0.05
- [REF] HYB-5252 DRG:HYBRID-XL Analyzer
- [REF] HYI-5392: *System Solution 5L*, 5000 mL;
(Instrument Feed Water according to CLSI guideline 3C-A4 with the following specification can also be used:
Resistivity at 25 °C [MΩ·cm]: > 1
Conductivity at 25 °C [µS·cm⁻¹]: < 1
Total Organic Carbon/p.p.b.[µg/L]: < 200
Colloids [µg/mL]: <0.1)
- [REF] HYI-5394: *Wash Buffer*, 40x concentrate, 25 mL
- [REF] HYI-5395: *Needle Cleaning Solution*, 30 mL. Cleaning solution for the pipetting needle (daily and weekly maintenance, see also user manual)
- [REF] HYI-5387: *Cuvettes*, (2 x 360 pieces)

For use of the *Secondary Sample Holder* for secondary tubes the following tubes are required:

- [REF] HYI-5391: *Sample Tubes (Secondary)*, 2500 pcs.

4.3 Storage Conditions

All kit components should be stored at 2 °C to 8 °C to ensure product performance until the defined expiry date.

When stored at 2 °C to 8 °C, **unopened kits** will retain reactivity until expiration date. Do not use reagents beyond this date.

- Cartridges (stored at 2 °C to 8 °C) in the supplied and unopened zip/foil bags will retain reactivity until expiration date.
- Unopened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity until expiration date.

Opened reagents and the reagent cartridges must be stored at 2 °C to 8 °C.

Once the plastic bag has been opened, care should be taken to tightly close it again along with the supplied desiccant bag.

Immediately after end of each run the Re-Calibrator and Control vials have to be removed from the instrument, tightly capped and stored at 2 °C to 8 °C.

- Unused cartridges in opened zip/foil bags (stored at 2 °C to 8 °C) will retain reactivity until expiration date, if stored as described above.
- Pierced or open cartridges must be disposed of immediately.

- Opened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity for 12 weeks.

4.3.1 On-board Stability

For Re-Calibrators and Controls the on-board stability has been evaluated under controlled laboratory conditions at room temperature (20 °C to 25 °C).

Due to the differences in laboratory environmental conditions and reagent volumes, the on-board stability may deviate from the declared value.

| | |
|---------------------|-----|
| On-board stability: | 8 h |
|---------------------|-----|

4.4 Reagent Preparation

Bring all reagents, such as controls and re-calibrators, to room temperature (20 °C to 25 °C) prior to use. Reagent Cartridges can be used directly after storage in the refrigerator.

Wash Buffer (not included in the kit)

For Wash Buffer (1x) dilute 25 mL of Wash Buffer (40x) with 975 mL ultra-pure water to a final volume of 1000 mL.

The diluted Wash Buffer (1x) is stable for 2 weeks at room temperature.

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum or plasma (EDTA-, heparin- or citrate plasma) can be used in this assay.

A minimum of 75 µL of sample is needed for one determination.

This includes 15 µL sample and 60 µL dead volume.

Attention:

- This test was not verified with blood collection tubes of all available manufacturers.
- Sample Collection Systems of some manufacturers may contain different materials which in isolated cases could affect the test results.
- If primary tubes for sample collection are used, please follow the instructions of the manufacturer.
- Do not use haemolytic, icteric or lipaemic specimens.
- Samples containing precipitates have to be centrifuged prior to the test run.
- Do not use heat inactivated samples.
- Do not use standards or external controls stabilized with azide.

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma:

Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 9 days at 2 °C to 8 °C prior to performing the assay.

Specimens stored for a longer time (up to six months) should be frozen only once at -20 °C prior to the assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

5.3.1 Manual Sample Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Sample Diluent* and reassayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10: 10 µL sample + 90 µL *Sample Diluent*
(mix thoroughly)
- b) dilution 1:100: 10 µL dilution a) 1:10 + 90 µL *Sample Diluent*
(mix thoroughly).

* *Sample Diluent for manual dilution is not included in this kit, but can be ordered on request (REF HYC-5613-DIL, 20 mL).*

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents, such as controls and re-calibrators and specimens, must be allowed to come to room temperature (20 °C to 25 °C) before use.
All reagents and samples must be mixed without foaming. Reagent Cartridges can be used directly after storage in the refrigerator.
- Samples, controls and re-calibrators should be measured within 2 hours in order to avoid possible evaporation effects.
- The *Secondary Sample Holder* (HYI-5437) for secondary tubes has the capacity for a maximum of 20 samples including controls and re-calibrators. They all have to be pipetted into the secondary tubes, and the respective barcodes of control/re-calibrator vials and, if available, the sample barcodes have to be read with the external barcode scanner.

6.2 Test Procedure

- **The total assay time for DRG:HYBRiD-XL CRP-HS is 10 minutes.**
- To ensure proper operation of the test, the instructions in the user manual for the DRG:HYBRiD-XL should strictly be followed.
- All test specific information required for the correct operation is included in the respective barcodes of the reagents.
Take care not to damage these bar codes!
- **It is recommended to tap the bottom of the Cartridge Segments containing the reagent cartridges once on the bench before placing them on the rotor. This is to avoid foam and adhering of the liquid on the sealing of the reagent cartridge.**
- Place reagent cartridges on the rotor of the unit. The heating to 37 °C incubation temperature is performed automatically in the unit.

6.3 Calibration

Traceability:

This method was standardized against the DRG CRP-HS ELISA (EIA-3954).

Each DRG:HYBRiD-XL reagent contains a barcode with the specific information for recalibration of the reagent lot. The Master Curve is printed as a 2-D barcode on the outer label of the kit package and on the QC-Datasheet and has to be scanned with the external barcode scanner prior to the first use of the respective kit lot.

Recalibration is recommended:

- if a new kit lot is used. Each new lot should be verified by running the kit internal re-calibrators and controls before routine use.
- if one or both assay controls are found outside the specified range.
- after 4 weeks of use of the same reagent kit on the unit.

6.4 Calculation of Results

The analyte concentrations are calculated automatically by the DRG:HYBRiD-XL's system software.

7 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

It is also recommended to participate in national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analysing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following: expiration dates and storage conditions of reagents, operational reliability of the analyser. In addition, it is indicated to perform a Recalibration.

In case of further questions please contact your local distributor or DRG directly.

7.1 Internal Controls

For Quality Control it is necessary to use the two internal controls provided with each kit.

Acceptance ranges for both internal controls (*Control 1 & 2*) were established by the manufacturer and are summarized in the QC-Datasheet added to the kit. Note that the expected values and acceptance ranges stated in the QC-Datasheet always refer to the current kit lot.

Internal controls should be run in single determination:

- on a routine basis (e.g. once per 24 h)
- if re-calibration is required (if one or both internal controls are out of range)
- if a new kit lot is used (in order to avoid any negative impact on the kit performance by improper transport or to detect improper storage during transport).

7.2 External Controls

Use controls at both normal and pathological levels.

The control intervals and control ranges for external controls should be adapted to the individual requirements of each laboratory. All results must be within the defined limits.

Each laboratory should establish corrective measures to be taken if values of external controls are not found in the acceptance range.

8 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and pathological values.

In a study conducted with healthy subjects, using the DRG:HYBRID-XL CRP-HS the following data were observed:

| Population | n | Mean (mg/L) | 5 th Percentile (mg/L) | 95 th Percentile (mg/L) |
|------------|----|-------------|-----------------------------------|------------------------------------|
| Adults | 43 | 0.84 | 0.14 | 8.66 |

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated with other clinical observations and diagnostic tests.

9 LIMITATIONS OF USE

Reliable and reproducible results will be obtained, when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10 PERFORMANCE CHARACTERISTICS

10.1 Assay Dynamic Range

The dynamic range of the assay is defined by the limit of detection and the maximum value of the Master Curve.

Values found below the measuring range are indicated as “< 0.054 mg/L”.

Values found above the measuring range are indicated as “> 10 mg/L”. The measuring range of the assay is between 0.054 mg/L – 10 mg/L.

10.2 Specificity of Antibodies (Cross-Reactivity)

Due to its antibodies, the DRG CRP-HS test is a specific immunoassay for human CRP.

10.3 Sensitivity

The sensitivity study was designed according to CLSI guideline EP17-A2.

The Limit of Blank (LoB) is 0.02 mg/L.

The Limit of Detection (LoD) is 0.054 mg/L.

The Limit of Quantification (LoQ) is 0.369 mg/L.

10.4 Precision Performance

The precision study was designed according to CLSI guideline EP5-A2.

10.4.1 Within-Run (Intra Assay)

Samples were analysed 2 times daily in duplicates for 20 days on 3 different instruments. CV (%) is given as mean of 40 duplicate CVs. The within-run variability is shown below:

| Sample | n | Mean (mg/L) | CV (%) |
|--------|----|-------------|--------|
| 1 | 40 | 0.29 | 6.8 |
| 2 | 40 | 1.21 | 2.7 |
| 3 | 40 | 2.09 | 3.1 |
| 4 | 40 | 4.93 | 2.2 |

10.4.2 Total Precision (Inter Assay)

Samples were analysed 2 times daily in duplicates for 20 days (n=80) on 3 different instruments. The total precision is shown below:

| Sample | n | Mean (mg/L) | CV (%) |
|--------|----|-------------|--------|
| 1 | 80 | 0.29 | 11.6 |
| 2 | 80 | 1.21 | 3.7 |
| 3 | 80 | 3.72 | 4.4 |
| 4 | 80 | 4.93 | 4.8 |

10.4.3 Inter-Lot

The between-lots variation was determined by 6 measurements of 4 samples with 3 different kit lots.

| Sample | n | Mean (mg/L) | CV (%) |
|--------|----|-------------|--------|
| 1 | 18 | 0.27 | 4.4 |
| 2 | 18 | 0.24 | 5.0 |
| 3 | 18 | 1.67 | 3.9 |
| 4 | 18 | 5.37 | 1.6 |

10.5 Recovery

Recovery was determined by adding four increasing concentrations of the analyte to four different patient samples containing different amounts of endogenous analyte. Each sample (non-spiked and spiked) was assayed and analyte concentrations of the samples were calculated from the Master Curve. The percentage recoveries were determined by comparing expected and measured values of the samples.

| Sample | 1 | 2 | 3 | 4 | |
|-----------------------|------|-------|-------|-------|-------|
| Concentration (mg/L) | 0.28 | 1.23 | 2.24 | 5.50 | |
| Average Recovery (%) | 97.5 | 101.4 | 99.5 | 102.6 | |
| Range of Recovery (%) | from | 91.5 | 97.9 | 97.8 | 102.0 |
| | to | 102.0 | 103.6 | 101.7 | 103.5 |

10.6 Linearity

Four samples containing different amounts of analyte were serially diluted with *Sample Diluent*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

| Sample | 1 | 2 | 3 | 4 | |
|-----------------------|-------|-------|-------|-------|-------|
| Concentration (mg/L) | 11.05 | 5.19 | 10.80 | 10.55 | |
| Average Recovery (%) | 109.0 | 97.2 | 108.8 | 107.6 | |
| Range of Recovery (%) | from | 106.1 | 95.3 | 106.9 | 88.6 |
| | to | 110.6 | 104.9 | 111.9 | 112.8 |

10.7 Interfering Substances

Haemoglobin (up to 5 mg/mL), Bilirubin (up to 0.3 mg/mL) and Triglyceride (up to 10.0 mg/mL) have no influence on the assay results. Until today no substances (drugs) are known to us, which have an influence to the measurement of CRP in a sample.


10.8 High-Dose-Hook Effect

Hook effect was not observed in this test up to a concentration of 122 mg/L of CRP.

11 METHOD COMPARISON

A comparison of DRG:HYBRID-XL CRP-HS Test HYC-5613 (y) and DRG ELISA EIA-3954 (x) using clinical samples gave the following correlation:

$$\begin{aligned} n &= 35 \\ r &= 0.989 \\ y &= 1.021x - 0.3402 \end{aligned}$$

12 LEGAL ASPECTS

Only for countries where the declaration of European Conformity (CE mark) is applicable.

12.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different kit lots could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

For further information please refer to the User Manual of the DRG:HYBRID-XL, analyser-specific application sheets, product information and package inserts of all necessary components.

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4 REAGENZIEN

4.1 Kitinhalt

4.1.1 Reagent Cartridges (Reagenzien-Cartridges)

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
13 REFERENCES / LITERATURE / LITERATUR / BIBLIOGRAFIA / LITERATURA

1. Agrwal A, Kilpatrick JM, Volanakis JE. Structure and function of human C-reactive protein. In: Mackiewicz A, Kushner I, Baumann H. eds. Acute phase proteins. London: CRC Press, 1993, 79-92.
2. Whicher JT, Rifai N, Biasucci LM. Markers of the acute phase response in cardiovascular disease: an update. Clin Chem Lab Med 2001, 39, 85-88.
3. Thomson D, Whicher JT, Banks RE. Acute phase reactions in predicting disease outcome. In Emery P, ed. Management of early inflammatory arthritis. Balliere's Clinical Rheumatology. London; Balliere Tynhall, 1992, 393-404.
4. Chamers RE, Stross P, Barry RE, Whicher JT. Serum amyloid A protein compared with C-reactive protein, alpha 1-antitrypsin and alpha 1-acid-glycoprotein a monitor of inflammatory bowel disease. Eur J Clin Invest. 1987, 17, 460-467.
5. Clyne B, Olshaker Js. The C-Reactive Protein. The Journal of Emergency Medicine. 1999, 17, 1019-1025.
6. Clinical Laboratory Diagnostics/ Lothar Thomas. 8th ed. Frankfurt: TH-Books Verlagsgesellschaft; 2012. 1278 ff.
7. Henry's Clinical Diagnosis and Management. 22nd ed. Philadelphia: W.B. Saunders Company; 2011. 254-255.
8. Tietz Fundamentals of Clinical Chemistry/Carl A. Burtis, Edward R. Ashwood - 5th ed. Philadelphia: W.B. Saunders Company; 2001. 325 ff.
9. Clinical Chemistry/Carl A. Burtis, Edward R. Ashwood - 5th ed. Philadelphia: W.B. Saunders Company; 2001. 325 ff.

SYMBOLS USED

| Symbol | English | Deutsch | Italiano | Español | Français | Polski |
|--|--|--------------------------------|-------------------------------------|--|--|---|
|  | European Conformity | CE-Konformitätskennzeichnung | Conformità europea | Conformidad europea | Conforme aux normes européennes | Zgodność z normami europejskimi |
|  | Consult instructions for use* | Gebrauchsanweisung beachten* | Consultare le istruzioni per l'uso | Consulte las Instrucciones | Consulter les instructions d'utilisation | Zapoznać się z instrukcją użytkowania |
|  | <i>In vitro</i> diagnostic medical device* | <i>In-vitro</i> -Diagnostikum* | Per uso Diagnostica in vitro | Para uso Diagnóstico in vitro | Usage Diagnostic in vitro | Wyrób medyczny do diagnostyki in vitro |
|  | For research use only | Nur für Forschungszwecke | Solo a scopo di ricerca | Sólo para uso en investigación | Seulement dans le cadre de recherches | Tylko do użytku w badaniach |
|  | Catalogue number* | Artikelnummer* | No. di Cat. | Número de catálogo | Référence | Numer katalogowy |
|  | Batch code* | Chargencode* | Lotto no | Número de lote | No. de lot | Numer LOT |
|  | Contains sufficient for <n> tests * | Ausreichend für <n> Prüfungen | Contenuto sufficiente per "n" saggi | Contenido suficiente para <n> ensayos | Contenu suffisant pour "n" tests | Zawartość przeznaczona na <n> testów |
|  | Temperature limit* | Temperaturbegrenzung* | Temperatura di conservazione | Temperatura de conservación | Température de conservation | Przechowywać w temperaturze |
|  | Use-by date* | Verwendbar bis* | Data di scadenza | Fecha de caducidad | Date limite d'utilisation | Data przydatności |
|  | Manufacturer* | Hersteller* | Fabbricante | Fabricante | Fabricant | Producent |
|  | Caution* | Achtung* | | | | UWAGA |
| <i>Distributed by</i> | Distributor | Vertreiber | Distributore | Distribuidor | Distributeur | Dystrybutor |
| <i>Content</i> | Content | Inhalt | Contenuto | Contenido | Contenu | Zawartość |
| <i>Volume/No.</i> | Volume / No. | Volumen/Anzahl | Volume/Quantità | Volumen/Número | Volume/Numéro | Objętość / Numer |
| <i>Reagent Cartridge</i> | Reagent Cartridge | Reagenzien-Cartridge | Cartucce di reagente | Cartucho de reactivo | Cartouche de réactif | Kartridż z odczynnikami |
| <i>Re-Calibrator</i> | Re-Calibrator | Re-Kalibrator | Re-Calibratore | Re-Calibrador | Re-calibrateurs | Re-kalibrator |
| <i>Control</i> | Control | Kontrolle | Controllo | Control | Contrôle | Kontrola |
| <i>Enzyme Conjugate</i> | Enzyme Conjugate | Enzymkonjugat | Tracciante enzimatico | Conjugado enzimático | Conjugué enzymatique | Koniugat enzymatyczny |
| <i>Enzyme Complex</i> | Enzyme Complex | Enzymkomplex | Complesso enzimatico | Complex enzimático | Complexe enzymatique | Enzymu kompleks |
| <i>Substrate Solution</i> | Substrate Solution | Substratlösung | Soluzione di substrato | Solución de sustrato | Solution substrat | Roztwór Substratu |
| <i>Sample Diluent</i> | Sample Diluent | Probenverdünnungsmedium | Diluente dei campioni | Solución para dilución de la muestra | Diluant d'échantillon | Rozcieńczalnik Próbkki |
| <i>Reagent 1</i> | Reagent 1 | Reagenz 1 | Reagente 1 | Reactivo 1 | Réactif 1 | Reagent 1 |
| <i>Reagent 2</i> | Reagent 2 | Reagenz 2 | Reagente 2 | Reactivo 2 | Réactif 2 | Reagent 2 |
| <i>Wash Buffer</i> | Wash Buffer | Waschpuffer | Tampone di lavaggio | Tampón de lavado | Tampon de lavage | Bufor Płuczący |
| <i>System Solution</i> | System Solution | Systemlösung | | Solución de sistema | Solution système | Roztwór systemowy |
| <i>Needle Cleaning Solution</i> | Needle Cleaning Solution | Nadel-Reinigungslösung | Soluzione lavaggio ago | Solución de lavado de la aguja | Solution de nettoyage des aiguilles | Roztwór Czyszczący Igły |
| <i>Denaturation Buffer</i> | Denaturation Buffer | Denaturierungspuffer | Tampone die denaturazione | Tampón de denaturalización | Tampon de dénaturation | |
| <i>Neutralization Buffer</i> | Neutralization Buffer | Neutralisierungspuffer | Tampone di neutralizzazione | Tampón de neutralización | Tampon de neutralisation | |
| <i>Assay Buffer</i> | Assay Buffer | Assaypuffer | Tampone del test | Tampón de ensayo | Tampon d'essai | Bufor testowy |
| | | | | | | |
| | | | | | | |
| <i>Secondary Sample Holder</i> | Secondary Sample Holder | Sekundärprobenhalter | Sostegno secondario dei campioni | Soporte de muestras para tubos secundarios | Support d'échantillons secondaires | Statyw na próbówki drugorzędowe |
| <i>Sample Tubes (Secondary)</i> | Sample Tubes (Secondary) | Sekundärprobenröhrchen | Tubetti di campioni secondari | Tubos de muestra secundarios | Tubes échantillon secondaires | Probówki drugorzędowe |
| <i>Dilution Cartridge</i> | Dilution Cartridge | | | | Cartouches de dilution | Kartridż do rozcieńczeń |
| <i>Vial Adapter</i> | Vial Adapter | Fläschchen-Adapter | l'adattatore dei tubetti | Adaptador de tubos | l'Adaptateur de flacons | Nakładki na butelki dla kontroli i kalibratorów |