Instructions for Use

Salivary Testosterone ELISA

IVD

REF SLV-3013

Σ 96

DRG

DRG Instruments GmbH, Germany
Frauenbergstraße. 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos / Sommaire

1	INTRODUCTION	2	1	INTRODUCCIÓN	26
2	PRINCIPLE OF THE TEST	2	2	FUNDAMENTO DEL ENSAYO	26
3	WARNINGS AND PRECAUTIONS	3	3	PRECAUCIONES	26
4	REAGENTS.....	4	4	COMPONENTES DEL KIT.....	27
5	SPECIMEN COLLECTION AND PREPARATION.....	5	5	MUESTRAS	28
6	ASSAY PROCEDURE.....	5	6	PROCEDIMIENTO DE ENSAYO	28
7	EXPECTED NORMAL VALUES.....	7	7	VALORES ESPERADOS	30
8	QUALITY CONTROL.....	7	8	CONTROL DE CALIDAD	30
9	PERFORMANCE CHARACTERISTICS.....	8	9	CARACTERÍSTICAS DEL ENSAYO.....	31
10	LIMITATIONS OF USE.....	10	10	LIMITACIONES DE USO	31
11	LEGAL ASPECTS	11	11	ASPECTOS LEGALES	32
1	EINLEITUNG	12	12	REFERENCES / LITERATURE	33
2	TESTPRINZIP	12		SYMBOLS USED.....	34
3	VORSICHTSMAßNAHMEN.....	12			
4	BESTANDTEILE DES KITS	13			
5	PROBENVORBEREITUNG.....	14			
6	TESTDURCHFÜHRUNG	14			
7	ERWARTETE WERTE	16			
8	QUALITÄTSKONTROLLE	16			
9	ASSAY CHARACTERISTIKA.....	17			
10	GRENZEN DES TESTS	17			
11	RECHTLICHE GRUNDLAGEN	18			
1	INTRODUZIONE	19			
2	PRINCIPIO DEL TEST	19			
3	PRECAUZIONI.....	19			
4	COMPONENTI DEL KIT.....	20			
5	CAMPIONI.....	21			
6	ATTUAZIONE DEL TEST.....	21			
7	VALORI NORMALI	23			
8	CONTROLLO QUALITÀ.....	23			
9	CARATTERISTICHE DEL TEST.....	24			
10	LIMITAZIONE DEL TEST.....	24			
11	ASPETTI LEGALI	25			

1 INTRODUCTION

1.1 Intended Use

The **Salivary Testosterone ELISA** is an enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of free active testosterone in saliva.

Measurement of testosterone is used in the diagnosis and treatment of disorders involving the male sex hormones (androgens), including primary and secondary hypogonadism, delayed or precocious puberty, impotence in males and, in females hirsutism (excessive hair) and virilization (masculinization) due to tumors, polycystic ovaries, and adrenogenital syndromes.

1.2 Summary and Explanation

Testosterone is the major androgenic steroid hormone. It is responsible for the development of the male external genitalia and secondary sexual characteristics. In females, its main role is as an estrogen precursor. In both genders, it also exerts anabolic effects and influences behavior.

In men, testosterone is secreted by the testicular Leydig cells and, to a minor extent, by the adrenal cortex. At puberty, testosterone increases dramatically in boys and progressively declines starting between the fourth and sixth decades of life. In premenopausal women, the ovaries are the main source of testosterone with minor contributions by the adrenals and peripheral tissues. After menopause, ovarian testosterone production is significantly diminished. Testosterone production in testes and ovaries is regulated via pituitary-gonadal feedback involving luteinizing hormone (LH) and, to a lesser degree, inhibins and activins. 40-50% of the total testosterone in blood is bound with high affinity to sex hormone-binding globulin (SHBG), while 40-50% is bound with low affinity to albumin, and only 1-2% can be detected as unbound or free testosterone.

During childhood, excessive production of testosterone induces premature puberty in boys and masculinization in girls. In adult women, excess testosterone production results in varying degrees of virilization, including hirsutism, acne, oligo-amenorrhea. Furthermore, it can cause polycystic ovaries, ovarian and adrenal tumors, adrenal hyperplasia and infertility. In men, mild-to-moderate testosterone elevations are usually asymptomatic while high levels of testosterone are associated with hypothalamic-pituitary-unit dysfunction, testicular tumors, congenital adrenal hyperplasia and prostate cancer.

Decreased testosterone in females causes subtle symptoms. These may include some decline in libido and nonspecific mood changes. In males, testosterone levels below the reference range indicate partial or complete hypogonadism, caused either by primary or secondary/tertiary (pituitary/hypothalamic) testicular failure. Low levels of testosterone are encountered in male patients with the following diseases: Klinefelter's syndrome, hypopituitarism, testicular feminization, orchidectomy, cryptorchidism, enzymatic defects and some autoimmune diseases.

2 PRINCIPLE OF THE TEST

The **Salivary Testosterone ELISA** Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding.

The microtiter wells are coated with a monoclonal (mouse) antibody directed towards an antigenic site of the testosterone molecule. Endogenous free testosterone present in a patient sample competes with a testosterone-horseradish peroxidase conjugate for binding to the coated antibody. After incubation the unbound conjugate is washed off.

The amount of bound peroxidase conjugate is inversely proportional to the concentration of testosterone in the sample. After addition of the substrate solution, the intensity of colour developed is inversely proportional to the concentration of testosterone in the patient sample.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with *Stop Solution* containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

4 REAGENTS

4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;
Wells coated with anti-testosterone antibody (monoclonal).
22. **Standard (Standard 0 - 6)**, 7 vials, 1 mL, ready to use;
Concentrations: 0 – 10 – 40 – 80 – 160 – 400 – 1000 pg/mL
Conversion: pg/mL x 3.47 = pmol/L
Standards are calibrated against mass spectrometry.
Contain non-mercury preservative.
23. **Control Low & High**, 2 vials, 1 mL each, ready to use;
For control values and ranges please refer to vial label or QC-Datasheet.
Contain non-mercury preservative.
24. **Enzyme Conjugate**, 1 vial, 14 mL, ready to use,
Testosterone conjugated to horseradish peroxidase;
Contains non-mercury preservative.
25. **Substrate Solution**, 1 vial, 25 mL, ready to use,
Tetramethylbenzidine (TMB).
26. **Stop Solution**, 1 vial, 14 mL, ready to use,
contains 0.5 M H₂SO₄,
Avoid contact with the stop solution. It may cause skin irritations and burns.
27. **Wash Solution**, 1 vial, 30 mL (40X concentrated),
See "Reagent Preparation".

Note: Additional *Standard 0* for sample dilution is available upon request.

4.2 Materials required but not provided

- A microtiter plate calibrated reader (450 ± 10 nm)
- Calibrated variable precision micropipettes (100 µL and 200 µL)
- Absorbent paper
- Distilled or deionized water
- Timer (60 min. range)
- Semi logarithmic graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C - 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C - 8 °C. Microtiter wells must be stored at 2 °C - 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for two months if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature prior to use.

Wash Solution

Add deionized water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL deionized water to a final volume of 1200 mL.

The diluted Wash Solution is stable for 2 weeks at room temperature.

4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet.

4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Eating, drinking, chewing gums or brushing teeth should be avoided for 30 minutes before sampling. Otherwise, it is recommended to rinse mouth thoroughly with cold water 5 minutes prior to sampling.

Do not collect samples when oral diseases, inflammation or lesions exist (blood contamination).

If there is visible blood contamination the patient specimen, it should be discarded, rinse the sampling device with water, wait for 10 minutes and take a new sample.

Note: Samples containing sodium azide should not be used in the assay.

5.1 Specimen Collection

Saliva samples should be collected only using special saliva sampling devices (vial and straw), e.g. SALI-TUBES 100 [REF] SLV-4158, available from DRG.

Do not use Salivette tubes for sampling.

Due to the cyclic secretion pattern of steroid hormones it is important to care for a proper timing of the sampling.

In order to avoid arbitrary results we recommend that 5 samples always be taken within a period of 2 – 3 hours (*multiple sampling*) preferably before a meal.

As food might contain significant amounts of steroid hormones samples preferably should be taken while fasting. If fasting should be a problem the collection period should be timed just before lunch or before dinner.

5.2 Specimen Storage and Preparation

The saliva samples may be stored at 2 °C to 8 °C up to one week, and should be frozen at -20 °C for longer periods; repeated thawing and freezing is no problem.

Each sample has to be frozen, thawed, and centrifuged at least once in order to separate the mucins by centrifugation.

Upon arrival of the samples in the lab the samples have to stay in the deep freeze at least overnight. Next morning the frozen samples are warmed up to room temperature and mixed carefully.

Then the samples have to be centrifuged for 5 to 10 minutes (at 3000 - 2000 × g).

Now the clear colorless supernatant is easy to pipette.

If a set of multiple samples is to be tested, the lab (after at least one freezing, thawing, and centrifugation cycle) has to mix the 5 single samples in a separate sampling device and perform the testing from this mixture.

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Standard 0* and re-assayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

a) Dilution 1:10: 10 µl saliva + 90 µL *Standard 0* (mix thoroughly)

b) Dilution 1:100: 10 µl of dilution a) + 90 µL *Standard 0* (mix thoroughly).

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Test Procedure

Each run must include a standard curve.

The assay should be performed at room temperature (21 °C to 24 °C). Lower temperatures will lead to clearly decreased OD values.

1. Secure the desired number of Microtiter wells in the frame holder.
2. Dispense **75 µL** of each **Standard, Control** and **samples** with new disposable tips into appropriate wells.
3. Dispense **100 µL Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Incubate for **60 minutes** at room temperature.
5. Briskly shake out the contents of the wells.
Rinse the wells **5 times** with diluted *Wash Solution* (400 µL per well). Strike the inverted wells sharply on absorbent paper towel to remove residual droplets.
Important note:
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
6. Add **200 µL** of **Substrate Solution** to each well.
7. Incubate for **30 minutes** at room temperature.
8. Stop the enzymatic reaction by adding **100 µL** of **Stop Solution** to each well.
9. Determine the absorbance (OD) of each well at **450 ± 10 nm** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Using semi-logarithmic graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 1000 pg/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Units (450 nm)
Standard 0 (0 pg/mL)	1.95
Standard 1 (10 pg/mL)	1.75
Standard 2 (40 pg/mL)	1.48
Standard 3 (80 pg/mL)	1.15
Standard 4 (160 pg/mL)	0.80
Standard 5 (400 pg/mL)	0.43
Standard 6 (1000 pg/mL)	0.15

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In order to determine the normal range of Salivary Testosterone, saliva samples from apparently healthy subjects were collected in the morning and analyzed using the Salivary Testosterone ELISA kit.

The following range was calculated from this study.

Table 1:

Population	n	Age (years)	Mean (pg/mL)	Median (pg/mL)	2.5 th - 97.5 th Percentile (pg/mL)	Range (Min. - Max.) (pg/mL)
Males	174	13 - 71	71.5	64.1	10.0 - 169.7	6.1 - 230.9
Females	67	15 - 73	19.8	18.9	3.5 - 44.9	3.1 - 50.2

Table 2:

Population	n	Age (years)	Mean (pg/mL)	Median (pg/mL)	2.5 th - 97.5 th Percentile (pg/mL)	Range (Min. - Max.) (pg/mL)
Males	8	14	30.7	19.6	13.8 - 59.0	13.5 - 59.02
	22	15	47.3	43.8	12.0 - 125.5	8.23 - 131.1
	10	16	69.9	68.9	7.5 - 152.3	6.1 - 167.8
	14	17	100.3	113.4	54.1 - 157.3	53.7 - 170.6
	12	18	125.1	135.0	61.3 - 209.7	53.7 - 230.9
	5	19	95.5	87.2	50.8 - 156.1	46.8 - 162.9
	15	21 - 25	96.4	93.3	56.1 - 168.1	52.9 - 186.5
	4	26 - 30	96.7	81.6	70.1 - 148.9	69.9 - 153.7
	2	31 - 35	74.7	74.7	54.4 - 95.0	53.3 - 96.1
	5	36 - 40	81.8	80.2	48.9 - 118.4	45.8 - 122.2
	12	46 - 50	64.3	56.8	28.6 - 127.4	23.03 - 127.4
	12	51 - 55	61.1	72.4	28.1 - 81.4	28.1 - 81.4
	4	56 - 60	79.7	79.7	68.5 - 90.8	68.5 - 90.8

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

Since testosterone levels show diurnal cycles, we recommend that the samples be obtained the same hour each day. Furthermore, we recommend that each laboratory determine its own range for the population tested.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 0.94 - 1000 pg/mL.

9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross reactivity of the assay:

Steroid	Cross reaction (%)
Testosterone	100
5 α -Dihydrotestosterone	0.80
Androstenedione	0.90
11 β -hydroxysterone	3.30
17 α -methyltestosterone	0.10
19-Nortestosterone	3.30
Epitestosterone	0.10
Estradiol	0.10
Progesterone	< 0.10
Cortisol	< 0.10
Estrone	< 0.10
Danazol	< 0.10

9.3 Sensitivity

The sensitivity study was designed according to CLSI guideline EP17-A2.

The Limit of Blank (LoB) was calculated by:

$$\text{LoB} = \text{mean}_{\text{blank}} + 1.645(\text{SD}_{\text{blank}}) = 0.94 \text{ pg/mL.}$$

The Limit of Detection (LoD); lowest analyte concentration likely to be reliably distinguished from LoB) was calculated by:

$$\text{LoD} = \text{LoB} + 1.645(\text{SD}_{\text{low concentration sample}}) = 2.63 \text{ pg/mL.}$$

The Limit of Quantification (LoQ); lowest actual amount of an analyte that can be reliably detected with CV < 20%) was determined as:

$$\text{LoQ} = 10.07 \text{ pg/mL.}$$

9.4 Reproducibility

9.4.1 Intra Assay

The intra-assay variation was determined by 20 replicate measurements of 3 saliva samples within one run. The within-assay variability is shown below:

Sample	n	Mean (pg/mL)	CV (%)
1	20	109.7	6.2
2	20	228.9	3.6
3	20	556.4	4.2

9.4.2 Inter Assay

The inter-assay (between-run) variation was determined by duplicate measurements of 4 saliva samples over 10 days. The between-assay variability is shown below:

Sample	n	Mean (pg/mL)	CV (%)
1	40	80.2	10.7
2	40	231.5	7.5
3	40	574.5	5.5
4	40	923.2	6.8

9.4.3 Inter-Lot

The inter-lot (between-lots) variation was determined by triplicate measurements of 3 saliva samples in three different kit lots. The between lot variability is shown below:

Sample	n	Mean (pg/mL)	CV (%)
1	18	146.9	4.6
2	18	267.9	14.4
3	18	586.7	11.1

9.5 Recovery

Samples have been spiked by adding testosterone solutions with known concentrations in a 1:1 ratio to 5 different saliva samples containing different amounts of endogenous analyte. Each sample (native and spiked) was assayed and analyte concentrations of the samples were calculated from the standard curve.

The % recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100 (expected value = (endogenous testosterone + added testosterone) / 2; because of a 1:2 dilution of serum with spike material).

		Sample 1	Sample 2	Sample 3	Sample 4	Sample 5
Concentration [pg/mL]		84.0	154.5	280.3	564.7	968.6
Average Recovery		99.7	103.3	100.9	104.5	102.0
Range of Recovery [%]	from	89.8	98.8	93.0	88.8	92.4
	to	110.4	109.1	112.2	110.2	106.7

9.6 Linearity


Five saliva samples containing different amounts of analyte were serially diluted with zero standard and assayed with the DRG ELISA. The percentage recovery was calculated by comparing the expected and measured values for testosterone.

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	
Concentration [pg/mL]	142.7	231.3	362.0	611.6	867.0	
Average Recovery	104.5	97.7	103.8	100.2	102.1	
Range of Recovery [%]	from	101.6	94.3	95.4	94.6	98.7
	to	109.2	101.3	111.6	109.0	110.8

9.7 Comparison Studies

A study was performed that evaluated 17 saliva samples collected from adult men and women. The samples were analyzed with the DRG ELISA and a LC-MS/MS mass spectrometer to determine the concentration of testosterone in the saliva samples.

A correlation of $r = 0.9904$ and regression formula of $y = 1.015x - 2.8203$ were obtained.

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

The patient should not eat, drink, chew gum or brush teeth for 30 minutes before sampling. Otherwise rinse mouth thoroughly with cold water 5 min prior to sample collection. Do not collect samples when oral diseases, inflammation or lesions exist (blood contamination).

10.1 Interfering Substances

Blood contamination of more than 0.16 % in saliva samples will affect results, and usually can be seen by eye. Concentrations of Sodium Azide ≥ 0.02 % interferes in this assay and may lead to false results.

10.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of testosterone in a sample.

10.3 High-Dose-Hook Effect

No hook effect was observed in this test.

11 LEGAL ASPECTS

Only for countries where the declaration of European Conformity (CE mark) is applicable.

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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
12 REFERENCES / LITERATURE

1. Mulhall JP, et al. (2004); Effect of testosterone supplementation on sexual function in hypogonadal men with erectile dysfunction, *Urology* 63, 348-352
2. Grinspoon S et al. (2000); Effects of hypogonadism and testosterone administration on depression indices in HIV-infected men, *J Clin. Endocrinol Metab.* 85, 60-65
3. Leder BZ, et al. (2004): Effect of aromatase inhibition in elderly men with low or borderline-low serum testosterone levels, *J Clin Endocrinol Metab.* 89, 1174-1180
4. Bunch T. et al. (2002): Pituitary radiographic abnormalities and clinical correlates of hypogonadism in elderly males presenting with erectile dysfunction, *Aging Male.* 5, 38-46
5. Ruutiainen K. et al. (1987): Salivary testosterone in hirsutism: correlations with serum testosterone and the degree of hair growth, *J Clin Endocrinol Metab.* 64, 1015-1020
6. Cserepes E. et al. (2002): Ovarian steroid cell tumor and a contralateral Ovarian thecoma in a postmenopausal woman with severe hyperandrogenism, *Gynecol Endocrinol.* 16, 213-216
7. Baldwin L. J. et al. (2004): Virilisation in a case of transitional cell carcinoma of the ovary, *J Clin Pathol.* 57, 1331-1332
8. Blumberg J. et al. (1991): Testicular tumors in congenital adrenal hyperplasia: steroid measurements from adrenal and spermatic veins, *J Clin Endocrinol Metab.* 73, 1129-1133.
9. Cardoso EM et al. (2011): Salivary testosterone for the diagnosis of androgen deficiency in end-stage renal disease. *Nephrol Dial Transplant.* 26; 677-83
10. Arregger AL et al. 2007: Salivary testosterone: a reliable approach to the diagnosis of male hypogonadism. *Clin Endocrinol (Oxf)* 67; 656-62
11. Ostatníková D et al. 2002: Salivary testosterone levels in preadolescent children. *BMC Pediatr.* 3; 2:5
12. Büttler RM et al. 2016: Reference values for salivary testosterone in adolescent boys and girls determined using Isotope-Dilution Liquid-Chromatography Tandem Mass Spectrometry (ID-LC-MS/MS). *Clin Chim Acta.* 456; 15-8
13. Karrer-Voegeli S. et al. 2009: Androgen dependence of hirsutism, acne, and alopecia in women: retrospective analysis of 228 patients investigated for hyperandrogenism. *Medicine (Baltimore).* 88; 32-45

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitätskennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Diagnostica in vitro	Diagnóstico in vitro	Diagnostic in vitro
	Catalogue number *	Artikelnummer *	No. di Cat.	No de catálogo	Référence
	Batch code *	Chargencode *	Lotto no	Número de lote	No. de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Caution *	Achtung *			
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement
<i>Volume/No.</i>	Volume / No.	Volumen / Anzahl	Volume / Quantità	Volumen / Número	Volume / Quantité
<i>Microtiterwells</i>	Microtiterwells	Mikrotiterwells	Micropozzetti	Placas multipocillo	Plaques de micro-titration
<i>Antiserum</i>	Antiserum	Antiserum	Antisiero	Antisero	Antisérum
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Tracciante enzimatico	Conjugado enzimático	Conjugué enzymatique
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complesso enzimatico	Complejo enzimático	Complexe enzymatique
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat
<i>Stop Solution</i>	Stop Solution	Stopplösung	Soluzione d'arresto	Solución de parada	Solution d'arrêt
<i>Zero Standard</i>	Zero Standard	Nullstandard	Standard zero	Estándar cero	Zero Standard
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampone del test	Tampón de ensayo	Tampon d'essai
<i>Wash Solution</i>	Wash Solution	Waschlösung	Soluzione di lavaggio	Solución de lavado	Solution de lavage
<i>1N NaOH</i>	1N NaOH	1N NaOH	1N NaOH (idrossido di sodio 1N)	1N NaOH	1N NaOH
<i>1 N HCl</i>	1 N HCl	1 N HCl		1 N HCl	1N HCl
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungs-medium	Diluyente dei campioni	Solución para dilución de la muestra	Solution pour dilution de l'échantillon
<i>Conjugate Diluent</i>	Conjugate Diluent	Konjugatverdünnungs-medium	Diluyente del tracciante	Solución para dilución del conjugado	Solution pour dilution du conjugué