# DRG ELISA

# Fasciola

# DRG

# DRG

# DRG ELISA

# Fasciola IgG (human)

## Superspecific ELISA for detection of Fasciola IgG antibodies

### Advantages of the DRG Fasciola IgG ELISA

- Low sample volume: 10  $\mu$ L Serum
- Detection of Fasciola hepatica and Fasciola gigantica IgG antibodies
- Easy to perform
- Overall incubation time: 1h 45 min
- High Sensitivity: 100%
- High Specificity: 100%
- No cross reaction with other parasites, e.g. Trichinella or Giardia
- This ELISA could be performed manually or with an open automated system

### Essential steps in DRG Fasciola IgG ELISA

1. During incubation Fasciola-specific IgG antibodies from the serum are bound to the immobilized antigen
2. After a washing step to remove unbound material, horseradish peroxidase conjugated anti-human IgG antibodies is added
3. Unbound conjugate is washed away
4. Substrate (TMB) is added for chromogenesis
5. Stop solution is added to terminate color development
6. Absorbance at 450 nm is read using an ELISA microtiter plate reader

### Fasciola IgG ELISA

#### EIA-4503

- Sample volume: 10  $\mu$ L Serum  
Inc. time: 60' 37° C, 30'/15' r.t.  
Sensitivity: 100%  
Specificity: 100%

## **Fasciola Diagnosis: The DRG Advantage**

Fasciolosis is caused by trematodes belonging to the genus *Fasciola* (*F. hepatica* and *F. gigantica*).

In the past, infection was limited to specific and typical geographical areas, but now widespread throughout the world. With human cases being increasingly reported from Europe, the Americas and Oceania (where only *F. hepatica* is transmitted) and from Africa and Asia (where the two species overlap).

As a consequence, human fasciolosis should be considered as a zoonosis of major global and regional importance.

Globally, the estimated number of human infections ranges from 2.4 million to 17 million.

*Fasciola hepatica* causes liver rot in sheep and cattle. Snails are the first intermediate host and encystation then occurs on aquatic vegetation. Humans usually acquire infection by eating contaminated freshwater plants, but can occasionally be infected by drinking unboiled contaminated water.

## **Diagnosis:**

In generally made via microscopy by identifying characteristic eggs in fecal samples or bile specimen. (However egg production does not begin until approx. 3 months after infection.)

Serological testing becomes positive during the early phase of migration through the liver and is therefore useful in diagnosing early symptoms before the appearance of eggs in the feces.

## **Clinical aspects:**

Adult worms usually reside in the bile duct, where they can live for many years and produce eggs that pass out with bile into the feces.

The early phase of migration of parasites through the liver can cause liver parenchymal destruction and be associated with fever, pain and hepatomegaly (6-12 weeks following ingestion).

## **DRG Fasciola IgG ELISA**

meets the claim for the early detection of Fasciolosis with the highest sensitivity (100%) and specificity (100%), no cross reaction was verifiable detected to other parasites, like: *Opisthorchis felineus*, *Toxocara canis*, *Schistosoma mansoni*, *Echinococcus granulosus/multilocularis*, *Trichinella spiralis*, *Strongyloides stercoralis*, *Ascaris lumbricoides*, *Giardia lamblia*. DRG Fasciola IgG ELISA is also applicable for the follow-up of treatment.

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
## DRG Diagnostics

DRG Instruments GmbH, founded in 1973 by Dr. Geacintov, subsidiary of DRG Intl. Inc., USA, is a diagnostics manufacturer of ELISAS. The DRG Group operates through a network of DRG subsidiaries in Germany, Poland, Russia, China and the Czech Republic and through distributors worldwide.

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## ELISAS that perform

DRG develops and manufactures diagnostic ELISA test kits for use in clinical and research laboratories. The experience of our production and management team guarantees to provide high quality products, competitive prices and excellent customer service.

DRG works to EN ISO 9001:2008 and EN ISO 13485:2003 + AC:2007 and ISO 13485:2003 under CMDCAS, certified by TÜV Rheinland LGA Products GmbH, an indication of our commitment to customer service, quality control and improved health care.

DRG Instruments GmbH, Germany  
Frauenbergstraße 18  
D-35039 Marburg  
Tel. +49 (0) 64 21/170 00,  
Fax +49 (0) 64 21/17 00 50  
Internet: www.drg-diagnostics.de  
E-mail: drg@drg-diagnostics.de

Distributed by

DRG International Inc. USA  
1167 U.S. Highway 22 East  
Mountainside, N.J. 07092 USA  
Phone: +1 (908) 233-2079  
Fax +1 (908) 233-0758  
Internet: www.drg-international.com  
E-mail: corp@drg-international.com