## Instructions for Use

# Thyroglobulin ELISA

IVD

CE

REF EIA-3377

Σ 96

DRG Instruments GmbH, Germany  
Frauenbergstraße. 18, D-35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

DRG International, Inc., USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.**  
**Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.**  
**Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.**  
**Por favor, se usa solo la versión válida de la metodico técnico incluido aquí en el kit.**  
**Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**  
**Por favor, use apenas a versão válida das instruções de utilização fornecidas com o kit.**

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## 1 INTENDED PURPOSE

Thyroglobulin ELISA is a test system for the quantitative measurement of thyroglobulin in human serum or plasma. This product is intended for professional in vitro diagnostic use only.

Elevated thyroglobulin levels are found in patients with pathological changes to the thyroid gland, e.g. hyperthyroidism, thyroiditis, thyroid carcinoma. The Thyroglobulin ELISA is primarily used for postoperative monitoring of patients with differentiated thyroid carcinoma.

## 2 PRINCIPLE OF THE TEST

Highly specific anti-human-thyroglobulin antibodies are bound to microwells.

The reaction is based on indirect enzyme immuno assay (ELISA) method with these steps:

thyroglobulin present in a patient sample binds to the antibody coated forming an antigen-antibody-complex. Washing of the microwells removes unbound unspecific serum and plasma components. During incubation with enzyme-conjugate immunologically a conjugate/antibody/antigen complex is formed. Washing of the microwells removes unbound conjugate. An anti-thyroglobulin enzyme substrate in the presence of bound conjugate hydrolyzes to form a blue colour. The addition of an acid stops the reaction forming a yellow end-product. The intensity of this yellow colour is measured photometrically at 450 nm. The amount of colour is directly proportional to the concentration of antibodies present in the original sample.

This assay includes a recovery test.

## 3 WARNINGS AND PRECAUTIONS

- All reagents of this kit are intended for professional in vitro diagnostic use only.
- Components containing human serum were tested and found negative for HBsAg, HCV, HIV1 and HIV2 by FDA approved methods. No test can guarantee the absence of HBsAg, HCV, HIV1 or HIV2, and so all human serum based reagents in this kit must be handled as though capable of transmitting infection.
- Bovine serum albumin (BSA) used in components has been tested for BSE and found negative.
- Avoid contact with the substrate TMB (3,3',5,5'-Tetramethyl-benzidine).
- Stop solution contains acid, classification is non-hazardous. Avoid contact with skin.
- Control, sample buffer and wash buffer contain sodium azide 0.09% as preservative. This concentration is classified as non-hazardous.
- Enzyme conjugate contains ProClin 300 0.05% as preservative. This concentration is classified as non-hazardous.

During handling of all reagents, controls and serum samples observe the existing regulations for laboratory safety regulations and good laboratory practice:

- First aid measures:  
In case of skin contact, immediately wash thoroughly with water and soap. Remove contaminated clothing and shoes and wash before reuse. If system fluid comes into contact with skin, wash thoroughly with water. After contact with the eyes carefully rinse the opened eye with running water for at least 10 minutes. Get medical attention if necessary.
- Personal precautions, protective equipment and emergency procedures:  
Observe laboratory safety regulations. Avoid contact with skin and eyes. Do not swallow. Do not pipette by mouth. Do not eat, drink, smoke or apply makeup in areas where specimens or kit reagents are handled. When spilled, absorb with an inert material and put the spilled material in an appropriate waste disposal.
- Exposure controls / personal protection:  
Wear protective gloves of nitril rubber or natural latex. Wear protective glasses. Used according to intended use no dangerous reactions known.
- Conditions to avoid:  
Since substrate solution is light-sensitive. Store in the dark.
- For disposal of laboratory waste the national or regional legislation has to be observed.

Observe the guidelines for performing quality control in medical laboratories by assaying control sera.

#### 4 CONTENTS OF THE KIT

Symbols		Sufficient for 96 determinations
<b>MICROPLATE</b>	1	One divisible <b>microplate</b> consisting of 12 modules of 8 wells each. Ready to use.
<b>CALIBRATOR A</b>	1x 1.5 mL	<b>Calibrator A, 0 ng/mL</b> , containing serum/buffer matrix (PBS, BSA, detergent, NaN <sub>3</sub> 0.09%), yellow. Ready to use.
<b>CALIBRATOR B</b>	1x 1.5 mL	<b>Calibrator B, 3 ng/mL</b> , containing thyroglobulin in a serum/buffer matrix (PBS, BSA, detergent, NaN <sub>3</sub> 0.09%), yellow. Ready to use.
<b>CALIBRATOR C</b>	1x 1.5 mL	<b>Calibrator C, 10 ng/mL</b> , containing thyroglobulin in a serum/buffer matrix (PBS, BSA, detergent, NaN <sub>3</sub> 0.09%), yellow. Ready to use.
<b>CALIBRATOR D</b>	1x 1.5 mL	<b>Calibrator D, 30 ng/mL</b> , containing thyroglobulin in a serum/buffer matrix (PBS, BSA, detergent, NaN <sub>3</sub> 0.09%), yellow. Ready to use.
<b>CALIBRATOR E</b>	1x 1.5 mL	<b>Calibrator E, 100 ng/mL</b> , containing thyroglobulin in a serum/buffer matrix (PBS, BSA, NaN <sub>3</sub> 0.09%), yellow. Ready to use.
<b>CALIBRATOR F</b>	1x 1.5 mL	<b>Calibrator F, 300 ng/mL</b> , containing thyroglobulin in a serum/buffer matrix (PBS, BSA, detergent, NaN <sub>3</sub> 0.09%), yellow. Ready to use.
<b>CONTROL +</b>	1x 1.5 mL	<b>Control positive</b> , containing thyroglobulin in a serum/buffer matrix (PBS, BSA, detergent, NaN <sub>3</sub> 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
<b>CONTROL -</b>	1x 1.5 mL	<b>Control negative</b> , containing thyroglobulin in a serum/buffer matrix (PBS, BSA, detergent, NaN <sub>3</sub> 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis.
<b>RECOVERY</b>	1x 3 mL	<b>Recovery</b> , 50 ng/mL, containing thyroglobulin in a serum/buffer matrix (PBS, BSA, detergent, NaN <sub>3</sub> 0.09%), yellow. Ready to use.
<b>DILUENT</b>	20 mL	<b>Sample Buffer STP</b> , containing PBS, BSA, detergent, preservative NaN <sub>3</sub> 0.09%, yellow, Ready to use.
<b>CONJUGATE</b>	15 mL	<b>Enzyme Conjugate</b> ; containing anti-human thyroglobulin antibodies, HRP labelled; PBS, BSA, detergent, preservative Proclin 0.05%, light red. Ready to use.
<b>TMB</b>	15 mL	<b>TMB Substrate</b> ; containing 3,3', 5,5'-Tetramethylbenzidin, colourless. Ready to use.
<b>STOP</b>	15 mL	<b>Stop Solution</b> ; contains acid. Ready to use.
<b>WASH</b>	20 mL	<b>Wash Buffer</b> , containing Tris, detergent, preservative NaN <sub>3</sub> 0.09%; 50X conc.
	1	Instruction for Use
	1	Certificate of Analysis

#### 5 MATERIALS REQUIRED

- Microplate reader capable of endpoint measurements at 450 nm; optional: reference filter at 620 nm
- Data reduction software
- Multi-channel dispenser or repeatable pipette for 100 µL
- Vortex mixer
- Pipettes for 50 µL, 100 µL and 1000 µL
- Laboratory timing device
- Distilled or deionised water
- Measuring cylinder for 1000 mL and 100 mL
- Plastic container for storage of the wash solution

This ELISA is suitable for use on open automated ELISA processors. Each assay has to be validated on the respective automated system. Detailed information is provided upon request.

## 6 SPECIMEN COLLECTION, STORAGE AND HANDLING

- Collect whole blood specimens using acceptable medical techniques to avoid hemolysis.
- Allow blood to clot and separate the serum or plasma by centrifugation.
- Test serum should be clear and non-hemolyzed. Contamination by hemolysis or lipemia should be avoided, but does not interfere with this assay.
- Specimens may be refrigerated at 2 °C - 8 °C for up to five days or stored at -20 °C up to six months.
- Avoid repetitive freezing and thawing of serum or plasma samples. This may result in variable loss of antibody activity.
- Testing of heat-inactivated sera is not recommended.

## 7 STORAGE AND STABILITY

- Store test kit at 2 °C - 8 °C in the dark.
- Do not expose reagents to heat, sun, or strong light during storage and usage.
- Store microplate sealed and desiccated in the clip bag provided.
- Shelf life of the unopened test kit is 18 months from day of production. Unopened reagents are stable until expiration of the kit. See labels for individual batch.
- Diluted Wash Buffer and Sample Buffer are stable for at least 30 days when stored at 2 °C - 8 °C.  
We recommend consumption on the same day.

## 8 PROCEDURAL NOTES

- Do not use kit components beyond their expiration dates.
- Do not interchange kit components from different lots and products.
- All materials must be at room temperature (20 °C - 28 °C) prior to use.
- Prepare all reagents and samples. Once started, perform the test without interruption.
- Double determinations may be done. By this means pipetting errors may become obvious.
- Perform the assay steps only in the order indicated.
- Always use fresh sample dilutions.
- Pipette all reagents and samples into the bottom of the wells.
- To avoid carryover or contamination, change the pipette tip between samples and different kit controls.
- Wash microwells thoroughly and remove the last droplets of wash buffer.
- All incubation steps must be accurately timed.
- Do not re-use microplate wells.

## 9 PREPARATION OF REAGENTS

### ***Wash Buffer (WASH)***

Dilute the contents of one vial of the buffered wash solution concentrate (50x) with distilled or deionised water to a final volume of 1000 mL prior to use.

### ***Sample Buffer STP (DILUENT)***

Sample buffer STP is ready to use.

## 10 PREPARATION OF SAMPLES

Use undiluted sample.

**Note:** Calibrators / Controls are ready to use and need not be diluted.

## 11 TEST PROCEDURE

Prepare enough microplate modules for all calibrators / controls and patient samples.

1. Pipette **50 µL** of calibrators, controls and patient samples into the wells.  
Calibrators and controls: add **50 µL** sample buffer  
Patient samples: add **50 µL** samples buffer (unspiked) /  
add **50 µL** RECOVERY (spiked)
2. Incubate for **60 minutes** at room temperature (20 °C - 28 °C).
3. Discard the contents of the microwells and **wash 3 times with 300 µL** of wash solution.
4. Dispense **100 µL** of enzyme conjugate into each well.
5. Incubate for **60 minutes** at room temperature.
6. Discard the contents of the microwells and **wash 3 times with 300 µL** of wash solution.
7. Dispense **100 µL** of TMB substrate solution into each well.
8. Incubate for **15 minutes** at room temperature
9. Add **100 µL** of stop solution to each well of the modules
10. Incubate for **5 minutes** at room temperature.
11. **Read** the optical density at 450 nm (reference 600-690 nm) and calculate the results. The developed colour is stable for at least 30 minutes. Read during this time.

Example for a pipetting scheme:

	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>
<b>A</b>	A	P1	P1+R									
<b>B</b>	B	P2	P2+R									
<b>C</b>	C	P3	P3+R									
<b>D</b>	D	P4	P4+R									
<b>E</b>	E	P5	P5+R									
<b>F</b>	F	P6	P6+R									
<b>G</b>	C+	P7	P7+R									
<b>H</b>	C-	P8	P8+R									

P1, ... patient sample (unspiked), P1+R, ... patient sample + **RECOVERY**, A-F calibrators, C+, C- controls

### 11.1 RECOVERY Test

The presence of autoantibodies against thyroglobulin (anti-TG) can interfere with the determination of human thyroglobulin (hTG) in patient samples: anti-TG can attach to epitopes of hTG molecules and thus cause false negative results in the determination of hTG. Therefore, it is necessary to prove the presence of anti-TG autoantibodies in patient samples. This can be done either by direct quantitative measurement with an anti-TG test (e.g. EIA-3561) or indirectly by recovery experiments in combination with the quantitative thyroglobulin determination.

In this Thyroglobulin assay a recovery test is included:

A patient sample is determined twice, unspiked and spiked with exogenous hTG called **RECOVERY** which contains exactly 50 ng/mL hTG. The recovery test provides evidence as to the presence of anti-TG autoantibodies.

Neither the correct anti-TG concentration nor the exact thyroglobulin concentration in presence of anti-TG can be calculated with this measurement.

## 12 VALIDATION

Test results are valid if the optical densities at 450 nm for calibrators / controls and the results for controls comply with the reference ranges indicated on the Certificate of Analysis enclosed in each test kit.

If these quality control criteria are not met the assay run is invalid and should be repeated.

## 13 CALCULATION OF RESULTS

For quantitative results plot the optical density of each calibrator versus the calibrator concentration to create a calibration curve. The concentration of patient samples may then be estimated from the calibration curve by interpolation.

Using data reduction software a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice.

The percentage recovery is calculated:

$$\% \text{ recovery} = (\text{ng/mL hTG spiked}/\text{ng/mL hTG non-spiked} + 50 \text{ ng/mL}) \times 100$$

Recovery should be expected in the range of 80-120 %.

If percentage of thyroglobulin recovery is below or above this range, thyroglobulin values for the respective patient sample should be excluded for further assessment.

## 14 PERFORMANCE CHARACTERISTICS

### 14.1 Calibration

The assay system is calibrated against the international Certified Reference Material CRM 457 from BCR, Brussels for human Thyroglobulin.

### 14.2 Measuring range

The calculation range of this ELISA assay is 0 - 300 ng/mL

### 14.3 Expected values

In a normal range study with samples from healthy blood donors the following ranges have been established with this ELISA assay:

Cut-off: 2 - 50 ng/mL

### 14.4 Interpretation of results

According to literature cut-off values for serum thyroglobulin of around 60 ng/mL, with a median of 5 to 10 ng/mL. In newborn babies as well as in pregnant woman of the 3rd trimester higher thyroglobulin concentrations may be detected. For comprehensive interpretation of thyroglobulin concentrations knowledge of alimentary iodine supply is indispensable. In regions with endemic goiter hTG values tend to be higher. In patients with total thyroidectomy no detectable thyroglobulin should be present after an X-ray therapy. Every increase of thyroglobulin to detectable serum concentrations is indicative for recidivation on thyroglobulin producing metastasis.

### 14.5 Linearity

Patient samples containing high levels of thyroglobulin were serially diluted in sample buffer to demonstrate the dynamic range of the assay and the upper / lower end of linearity. Activity for each dilution was calculated from the calibration curve using a 4-Parameter-Fit with lin-log coordinates.

Sample	Dilution	Observed [ng/mL]	Expected [ng/mL]	O/E [%]
1	1:1	268.0	268.0	100
	1:2	141.0	134.0	105
	1:4	65.0	67.0	97
	1:8	31.0	34.0	91
2	1:1	207.0	207.0	100
	1:2	101.0	104.0	97
	1:4	50.0	52.0	96
	1:8	24.0	26.0	92

### 14.6 Limit of detection

Functional sensitivity was determined to be: 1 ng/mL

## 14.7 Reproducibility

### Intra-assay precision:

Coefficient of variation (CV) was calculated for each of three samples from the results of 24 determinations in a single run. Results for precision-within-assay are shown in the table below.

### Inter-assay precision:

Coefficient of variation (CV) was calculated for each of three samples from the results of 6 determinations in 5 different runs. Results for run-to-run precision are shown in the table below.

Intra-Assay		
Sample	Mean [ng/mL]	CV [%]
1	33.0	1.9
2	93.0	2.4
3	227.0	3.2

Inter-Assay		
Sample	Mean [ng/mL]	CV [%]
1	31.0	1.7
2	88.0	1.7
3	212.0	1.1

## 14.8 Interfering substances

No interference has been observed with haemolytic (up to 1000 mg/dL) or lipemic (up to 3 g/dL triglycerides) sera or plasma, or bilirubin (up to 40 mg/dL) containing sera or plasma.

Nor have any interfering effects been observed with the use of anticoagulants (Citrate, EDTA, Heparin).

However for practical reasons it is recommended that grossly hemolyzed or lipemic samples should be avoided.

## 15 LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire clinical picture of the patient. Also every decision for therapy should be taken individually.

The above pathological and normal reference ranges for antibodies in patient samples should be regarded as recommendations only. Each laboratory should establish its own ranges according to ISO 15189 or other applicable laboratory guidelines.

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Echantillon	Dilution	Observé [ng/mL]	Escompté [ng/mL]	O/E [%]
1	1:1	268,0	268,0	100
	1:2	141,0	134,0	105
	1:4	65,0	67,0	97
	1:8	31,0	34,0	91
2	1:1	207,0	207,0	100
	1:2	101,0	104,0	97
	1:4	50,0	52,0	96
	1:8	24,0	26,0	92

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1. Marriq, S. et al. Polypeptide chains of 19-S thyroglobulin from several mammalian species and of porcine 27-S iodoprotein. Eur. J. Biochem. 1977; 79: 143- 149.
2. Gärtner, D.F. et al. Evidence for autonomous thyroglobulin release from euthyroid and hyperthyroid nodular goiters - Thyroglobulin, a possible helpful parameter in diagnosis of non-malignant thyroid disorders. Klin. Wochenschr. 1983; 61: 737 - 741.
3. Gebel, F. et al. The site of leakage of intrafollicular thyroglobulin into the blood stream in simple human goiter. J. Clin. Endocrinol. Metab. 1983; 57: 915 -919.
4. Uller, R.P. and van Herle, A.J. Effect of therapy on serum thyroglobulin levels in patients with Graves` disease. J. Clin. Endocrinol. Metab. 1978; 46: 747 -755.
5. Gardner, et al. Serum thyroglobulin in normal subjects and patients with hyperthyroidism due to Graves` disease: effects of T3, iodine, 131J, and antithyroid drugs. Clin. Endocr. (Oxf.) 1979; 11: 585 - 594.
6. Kawamura, S. et al. Serum thyroglobulin changes in patients with Graves` disease treated with long term antithyroid drug therapy. J. Clin. Endocrinol. Metab. 1983; 56: 507 - 512.
7. Czernichow, P. et al. Plasma thyroglobulin measurements help determine the type of thyroid defect in congenital hypothyroism. J. Clin. Endocrinol. Metab. 1983; 56: 242 - 245.
8. Mariotti, S. et al. Low serum thyroglobulin as a clue to the diagnosis of thyrotoxicosis factitia. New Engl. J. Med. 1982; 307: 410 - 412.
9. Reiners, C. Klinische Wertigkeit der Thyreoglobulinbestimmung im Serum. Akt. Endokr. Stoffw. 1984; 5: 76 - 82.
10. Kastrup, J. et al. An enzyme linked immunosorbent assay for measurement of human serum thyroglobulin. Evaluation of the influence of thyroglobulin autoantibodies. Scand. J. Lab. Invest. 1985; 45: 471 - 476.
11. Feldt-Rasmussen, U. et al. Serum thyroglobulin (Tg) in presence of thyroglobulin autoantibodies (TgAb). Clinical and methodological relevance of the interaction between Tg and TgAb in vitro and in vivo. J. Endocrinol. Invest. 1985; 8: 571 - 576.
12. van Herle, A.J. and Uller, R.P Elevated serum thyroglobulin: a marker of metastases in differentiated thyroid carcinomas. J. Clin. Invest. 1975; 56: 272 -277.

## SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français	Português
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes	Conformidade Europeia
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation	Consultar as instruções de uso
	In vitro diagnostic medical device *	In-vitro-Diagnostikum *	Diagnistica in vitro	Diagnóstico in vitro	Diagnostic in vitro	Dispositivo médico para diagnóstico in vitro
	Catalogue number *	Katalognummer *	No. di Cat.	No de catálogo	Référence	Número de catálogo
	Batch code *	Chargen-bezeichnung *	Lotto no	Número de lote	No. de lot	Código do lote
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests	Suficiente para <n> determinações
	Temperature limit *	Temperaturgrenzwerte *	Temperatura di conservazione	Temperatura de conservacion	Temperature de conservation	Limites de temperatura
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation	Prazo de validade
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant	Fabricante
	Distributor *	Vertriebspartner *	Distributore	Distribuidor	Distributeur	Distribuidor
	Date of manufacture *	Herstellungsdatum *	Data di produzione	Fecha de fabricación	Date de production	Data de fabricação
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques	Riscos biológicos
	Caution *	Achtung *	Attenzione	Precaución	Attention	Cuidado
	Unique device Identifier *	eindeutige Produktidentifizierung *	Identificativo unico del dispositivo*	Identificación exclusiva del dispositivo *	Identifiant de dispositif unique*	Identificador único do dispositivo *
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches	
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement	Conteúdo
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité	Volume / Quantidade