## Instructions for Use

# ASCA IgA ELISA

**REF** EIA-4279

 **96**

**DRG Instruments GmbH**, Germany  
Frauenbergstraße 18, 35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

**Distributed by:**

**DRG International, Inc.**, USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.  
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.  
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.  
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.**

### Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos

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## 1 INTENDED PURPOSE

The ASCA IgA ELISA is a quantitative immunoassay for the determination of IgA antibodies against *Saccharomyces cerevisiae* in human serum.

The ASCA IgA ELISA is intended as an aid in the diagnosis of Crohn's disease in conjunction with other clinical and laboratory findings.

The immunoassay is designed for manual professional in vitro diagnostic use.

## 2 DIAGNOSTIC RELEVANCE

Non-specific inflammatory bowel diseases including Crohn's disease (Enteritis regionalis) and ulcerative colitis (UC) are characterized by unknown etiology as well as chronic-remitting inflammatory processes of the intestine. Whereas the inflammation of ulcerative colitis is restricted to the mucosa and submucosa of colon and rectum, Crohn's disease (CD) shows a wide spread inflammation of the gastro-intestinal tract with granuloma formation.

The risk developing one of these diseases is strongly influenced by immunologic, genetic, infectious and environmental factors.

The differential diagnosis of inflammatory bowel diseases to chronic diarrhea, recurrent abdominal dolor, infectious colitis, anorexia as well as the differentiation of CD to ulcerative colitis is still a high challenge.

The determination of IgA and IgG antibodies to *Saccharomyces cerevisiae* (baker's yeast) has been described as one important serological marker for the differential diagnosis of Crohn's disease recently. Up to 70 % of patients with CD show antibody levels to *Saccharomyces cerevisiae*. Although the cause for their occurrence has been unclear, antibodies to *Saccharomyces cerevisiae* (ASCA) are strongly associated with inflammatory processes of the intestine.

In combination with the detection of autoantibodies to atypical anti-neutrophil cytoplasmic antigens (aANCA) which are mainly found in patients with ulcerative colitis, ASCA are a valid parameter for the differentiation of Crohn's disease and ulcerative colitis.

## 3 TEST PRINCIPLE

The ELISA (Enzyme Linked Immunosorbent Assay) is an immunoassay for the determination of specific antibodies. The strips of the microtiter plate are coated with test-specific antigens. If antibodies are present in the patient's sample, they bind to the antigens. A secondary antibody conjugated with the enzyme peroxidase detects the generated immune complex. A colorless substrate is converted into the colored product. The signal intensity of the reaction product is proportional to the antibody activity in the sample. After stopping the signal intensity of the reaction product is measured photometrically.

**4 TEST COMPONENTS**

Component	Description
<b>Microtiter plate</b> <b>A</b> Ag   96 1 piece	12 breakable microtiter strips (ready-to-use), 8 wells per strip, each well coated with mannan from <i>Saccharomyces cerevisiae</i>
<b>Calibrator</b> <b>0 – 4</b> CAL 5 x 1 mL, white cap	Colored dilutions of human serum (ready-to-use; contains ProClin 950) <i>The antibody activities are indicated on the quality control certificate.</i>
<b>Negative control</b> <b>N</b> CONTROL   - 1 x 1 mL, green cap	Colored dilution of human serum (ready-to-use; contains ProClin 950) <i>The antibody activity is indicated on the quality control certificate.</i>
<b>Positive control</b> <b>P</b> CONTROL   + 1 x 1 mL, red cap	Colored dilution of human serum (ready-to-use; contains ProClin 950) <i>The antibody activity is indicated on the quality control certificate.</i>
<b>Sample diluent</b> <b>C</b> DIL 1 x 100 mL, black cap	Colored solution (ready-to-use; contains ProClin 950)
<b>Wash buffer</b> <b>B</b> WASH   10x 1 x 100 mL, white cap	Concentrated solution (10x; contains ProClin 950)
<b>Conjugate IgA</b> <b>D</b> CONJ 1 x 15 mL, violet cap	Colored solution of polyclonal anti-human IgA antibody conjugated to horseradish peroxidase (ready-to-use; contains ProClin 950)
<b>Substrate</b> <b>E</b> SOLN   TMB 1 x 15 mL, blue cap	3,3',5,5'-Tetramethylbenzidin (ready-to-use)
<b>Stop solution</b> <b>F</b> H2SO4   0.25M 1 x 15 mL, yellow cap	0.25 M Sulfuric acid (ready-to-use)
<b>Adhesive Foil</b> 2 pieces	-
<b>QC Certificate</b> 1 piece	-
<b>Instructions for Use</b> 1 piece	-

**5 MATERIALS REQUIRED IN ADDITION**

- Common laboratory equipment
- Precision pipettes (5 – 1000 µL), multi-channel pipettes (100 – 1000 µL) and disposable pipette tips
- Graduated cylinders (100 – 1000 mL)
- Sample tubes for the preparation of dilutions
- Vortex mixer or other rotators
- Microtiter plate washer or wash comb
- Microtiter plate reader with optical filters for 450 nm and 620 nm or 690 nm
- Adsorbent paper or paper towel
- Distilled or de-ionized water

**6 STORAGE AND STABILITY**

Upon receipt, all test components must be stored at 2 °C to 8 °C, preferably in the original kit box. If stored properly in their original containers, all components are stable until their expiry date. All components are stable for at least 2 months after opening when stored properly at 2 °C to 8 °C.

## 7 GENERAL INFORMATION

This product is for in vitro diagnostic use only. The instructions for use must be carefully read before use. They are valid only for the present product with the given composition and must be strictly followed to ensure reliable test results. Deviations can lead to erroneous test results. Components must not be exchanged by test reagents of different lots or of other manufacturers.

Contamination of reagents must be avoided by use of aseptic techniques when removing aliquots from the vials. After use, reagent vials must be tightly closed with their corresponding caps.

Cross-contamination of samples or reagents can lead to inconsistent test results and must be avoided by use of consistent pipetting techniques.

Exposure of reagents to strong light must be avoided throughout the entire test procedure and storage.

Insufficient washing will result in poor precision and elevated measurement signals. After each washing step any residual fluid has to be removed completely.

## 8 PREPARATION

### 8.1 Preparation of Reagents

All components including the microtiter plate must be brought to room temperature (RT: 18 °C to 25 °C) before use for at least 30 min. All liquid components must be mixed gently to ensure homogeneity.

#### 8.1.1 Microtiter Plate

The microtiter plate is sealed in an aluminium bag. Unused test strips should always be stored refrigerated and protected from moisture with the desiccant in the properly sealed aluminum bag. Carefully resealed, the test strips can be used for 8 weeks after opening.

#### 8.1.2 Calibrators

The calibrators are ready-to-use and must not be diluted any further. Calibrators must be used in each test run.

#### 8.1.3 Controls

The positive and the negative controls are ready-to-use and must not be diluted any further. Controls must be used in each test run. Laboratories can also validate their own control samples and use them alternatively.

#### 8.1.4 Sample Diluent

The sample diluent is ready-to-use.

#### 8.1.5 Wash Buffer

The wash buffer is concentrated and must be diluted **1:10** with distilled water before use (e. g. 100 mL + 900 mL). A sufficient amount of washing solution must be prepared. The diluted washing solution can be stored at 2 °C to 8 °C up to 30 days.

#### 8.1.6 Conjugate

The conjugate is ready-to-use and stable up to 8 weeks after opening when stored at 2 °C to 8 °C.

#### 8.1.7 Substrate

The substrate is ready-to-use. Exposure of the substrate solution to strong light should be avoided.

#### 8.1.8 Stop Solution

The stop solution is ready-to-use.

### 8.2 Preparation of Samples

#### 8.2.1 Sample Material

The use of freshly collected **serum** from blood taken by venipuncture is recommended.

The use of icteric, lipemic, hemolytic or bacterially contaminated samples should be avoided. Insoluble substances must be removed from the sample by centrifugation. Samples must not be thermally inactivated.

#### 8.2.2 Sample Dilution

The samples must be **diluted 1:101** (e. g. 10 µL + 1000 µL) with sample diluent and mixed thoroughly. Building of foam should be avoided.

#### 8.2.3 Sample Storage

Samples may be kept at 2 °C to 8 °C up to three days.

Long-term storage requires -20 °C. Repeated freezing and thawing should be avoided. For multiple use, samples should be aliquoted and kept at -20 °C.

## 9 TEST PERFORMANCE

### 9.1 Pipetting Scheme

The following pipetting scheme is recommended:

	1	2	3	4
A	CAL 0	Sample 2		
B	CAL 1	Sample 3		
C	CAL 2	Sample 4		
D	CAL 3	Sample 5		
E	CAL 4	...		
F	N	...		
G	P	...		
H	Sample 1	...		

### 9.2 Procedure

The indicated incubation times and temperatures must be adhered to, and significant time shifts during pipetting samples and reagents must be avoided. The microtiter plate should be shortly shaken after addition of reagents.

Step	Description
1. Addition of calibrators, controls and diluted samples	Add 100 µL ready-to-use calibrators, controls and diluted samples per well
2. Incubation	Cover the plate and incubate for 60 min. at RT
3. Wash cycle	Aspirate the solution and wash 3 times with 300 µL washing solution with at least 5 seconds soaking time each; dry by tapping the microtiter plate on a paper towel to remove any residual droplets
4. Addition of conjugate	Add 100 µL ready-to-use conjugate to each well
5. Incubation	Cover the plate and incubate for 30 min. at RT
6. Wash cycle	Aspirate the solution and wash 3 times with 300 µL washing solution with at least 5 seconds soaking time each; dry by tapping the microtiter plate on a paper towel to remove any residual droplets
7. Addition of substrate	Add 100 µL ready-to-use substrate to each well
8. Incubation	Cover the plate and incubate for 15 min. in the dark at RT
9. Addition of Stop Solution	Add 100 µL ready-to-use stop solution to each well
10. Analysis	Read optical density (OD) at 450 nm versus 620 or 690 nm within 30 min. after stopping the reaction

### 9.3 Automation

Automated processing of the immunoassays must be performed analogous to manual use and validated by the user.

## 10 TEST EVALUATION

### 10.1 Metrological Traceability

The immunoassay is calibrated using an internal reference sample. Quantitative results are expressed in U/mL.

### 10.2 Standard Curve

For generation of a standard curve, the optical signals (optical density, OD) of the calibrators are plotted against their antibody activities and correlated by a 4-parameter logistic (4 PL) fit. Antibody activities of unknown samples can be derived directly from their optical signals by use of the generated standard curve.

### 10.3 Criteria of Validity

Test runs are only valid if the following criteria of validity are fulfilled:

- OD CAL 0 < CAL 1 < CAL 2 < CAL 3 < CAL 4
- OD CAL 4 > 1.2
- The negative control must be evaluated negative.
- The positive control must be evaluated positive and present an antibody activity within the validity range indicated on the quality control certificate.

If these criteria are not met, the test is not valid and must be repeated.

### 10.4 Troubleshooting

In case of an invalid test run, the expiry dates and storage conditions, incubation times and temperatures, and precise calibration of all instruments used should be verified. If no reason for an invalid test run could be identified, please contact the supplier or manufacturer of the product.

### 10.5 Reference Ranges

The reference ranges are indicated below:

	Interpretation
Antibody activity < 20 U/mL	negative
Antibody activity ≥ 20 U/mL	positive

As a result of different seroprevalences in individual regions, each laboratory should verify the reference ranges by own analysis and adapt, if necessary.

### 10.6 Interpretation of Test Results

A positive test result indicates the presence of specific antibodies.

A negative result indicates the absence of specific antibodies, but does not exclude the possibility of an autoimmune reaction.

In case of a borderline test result, a reliable evaluation is not possible.

### 10.7 Limitations of the Method

The interpretation of test results must always be considered in combination with the clinical picture of the patient. The diagnosis should not be based on the results of a sole diagnostic method. All clinical and laboratory findings should be evaluated to state a diagnosis. For confirmation, further investigations should be carried out.

## 11 PERFORMANCE CHARACTERISTICS

### 11.1 Analytical Performance Characteristics

#### 11.1.1 Precision

The precision of test results was assessed by the determination of the intra- and inter-assay variation by the analysis of multiple samples with different antibody activities.

	Intra-assay Precision		Inter-assay Precision	
	U/mL	CV (%)	U/mL	CV (%)
Sample 1	199	6.6	215	8.2
Sample 2	65	4.4	65	6.8
Sample 3	18	6.7	19	8.6

### 11.1.2 Measurement Range

Reliable accuracy, trueness, precision, linearity and recovery of test results have been observed within the measurement range of the assay from the LoQ to the upper calibrator in comprehensive studies.

Samples with test results above the upper calibrator should be reported as > max.

Samples with test results below the LoQ should be reported as < min.

If test results above the upper calibrator are observed, the samples may be tested at a higher dilution. The resulting antibody activity must be multiplied with the additional dilution factor.

## 11.2 Diagnostic Performance Characteristics

### 11.2.1 Diagnostic Sensitivity and Specificity

Sensitivity and specificity were assessed by the analysis of 116 serum samples from patients with Crohn`s disease and 75 serum samples from patients with ulcerative colitis:

	Diagnostic Performance
Sensitivity	42.0 %
Specificity	92.0 %

## 12 WARNINGS AND PRECAUTIONS

The product is designed exclusively for *in vitro* diagnostic use by qualified, authorized and trained personnel. All test components and human samples should be handled with care as potentially hazardous. Good laboratory practices (GLP) and all relevant regulations should be adhered to.

In case the product is damaged or product information including labelling is wrong or incorrect, please contact the manufacturer or supplier.

This product contains preparations of human and / or animal origin. Any material derived from human body fluids or organs used for the preparation of components were tested and found negative for HBsAg (Hepatitis B-Virus-surface Antigen) and anti-HIV as well as anti-HCV antibodies. However, all components and all patient samples should be handled as potentially hazardous in accordance with national laws and appropriate guidelines on biological safety.

As the product contains potentially hazardous materials, the following precautions should be followed: Do not smoke, eat or drink while handling kit material or samples. Avoid direct contact to kit material or samples by wearing protective gloves, laboratory coat and safety glasses. Never pipette material by mouth. Wipe up spills promptly and wash the affected surface thoroughly with a decontaminant. Wash hands thoroughly after use.

Some of the reagents contain ProClin (< 1.0 %) as a preservative, may cause skin sensitization (H317) and must not be swallowed or allowed to come into contact with skin or mucosa (P280, P333+P313).

The information in the safety data sheet on possible hazards, first aid measures, measures in the event of the unintentional release of large quantities, handling and storage, personal protective equipment, information on disposal as well as information on toxicology must be observed.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the member state in which the user and/or the patient is established.

## 13 DISPOSAL

For decontamination and disposal the recommendations of the CDC as well as the relevant local and national environmental guidelines and regulations should be adhered to.

Samples, potentially contaminated materials and infectious waste must be decontaminated, e.g. by autoclaving for 20 minutes at 121 °C.

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## 14 REFERENZEN

- Conrad K, Schmechta H, Klafki A, Lobeck G, Uhlig HH, Gerdi S, Henker J: Serological differentiation of inflammatory bowel diseases. Eur J Gastro & Hepatol. 2002, 14, 129 – 35.
- Vermeire S: Serological Diagnosis in IBD. IBDM 2002, 3, 82 – 9.
- Degenhardt F, Dirmeier A, Lopez R, Lang S, Kunst C, Roggenbuck D, Reinhold D, Szymczak S, Rogler G, Klebl F, Franke A, Rieder F. Serologic Anti-GP2 Antibodies Are Associated with Genetic Polymorphisms, Fibrostenosis, and Need for Surgical Resection in Crohn's Disease. Inflamm. Bowel Dis. 2016, 22, 2648 – 57.

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## 14 RIFERIMENTI

- Conrad K, Schmechta H, Klafki A, Lobeck G, Uhlig HH, Gerdi S, Henker J: Serological differentiation of inflammatory bowel diseases. *Eur J Gastro & Hepatol.* 2002, 14, 129 – 35.
- Vermeire S: Serological Diagnosis in IBD. *IBDM* 2002, 3, 82 – 9.
- Degenhardt F, Dirmeier A, Lopez R, Lang S, Kunst C, Roggenbuck D, Reinhold D, Szymczak S, Rogler G, Klebl F, Franke A, Rieder F. Serologic Anti-GP2 Antibodies Are Associated with Genetic Polymorphisms, Fibrostenosis, and Need for Surgical Resection in Crohn's Disease. *Inflamm. Bowel Dis.* 2016, 22, 2648 – 57.

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## 14 REFERENCIAS

1. Conrad K, Schmechta H, Klafki A, Lobeck G, Uhlig HH, Gerdi S, Henker J: Serological differentiation of inflammatory bowel diseases. *Eur J Gastro & Hepatol.* 2002, 14, 129 – 35.
2. Vermeire S: Serological Diagnosis in IBD. *IBDM* 2002, 3, 82 – 9.
3. Degenhardt F, Dirmeier A, Lopez R, Lang S, Kunst C, Roggenbuck D, Reinhold D, Szymczak S, Rogler G, Klebl F, Franke A, Rieder F. Serologic Anti-GP2 Antibodies Are Associated with Genetic Polymorphisms, Fibrostenosis, and Need for Surgical Resection in Crohn's Disease. *Inflamm. Bowel Dis.* 2016, 22, 2648 – 57.

**SYMBOLS USED**

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
<b>IVD</b>	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico In vitro	Dispositif médical de diagnostic in vitro
<b>REF</b>	Catalogue number *	Artikelnummer *	Numero di Catalogo	Número de catálogo	Référence de catalogue
<b>LOT</b>	Batch code *	Fertigungslosnummer, Charge *	Codice del lotto	Codigo de lote	Numéro de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis *	Utilizzare prima del	Establa hasta	Utiliser jusque
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Biological risks*	Biologische Risiken*	Rischi biologici	Riesgos biológicos	Risques biologiques
	Caution *	Achtung *	Attenzione	Precaución	Attention
<b>RUO</b>	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité