## Instructions for Use

# 25-OH Vitamin D (total) ELISA

IVD

REF EIA-5396

$\Sigma$  96

**DRG Instruments GmbH**, Germany  
Frauenbergstraße. 18, D-35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

Distributed by:

**DRG International, Inc.**, USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.  
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.  
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.  
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.**

**Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos**

1	INTRODUCTION.....	2
2	PRINCIPLE OF THE TEST .....	2
3	WARNINGS AND PRECAUTIONS .....	3
4	REAGENTS .....	4
5	SPECIMEN COLLECTION AND PREPARATION.....	6
6	ASSAY PROCEDURE .....	6
7	EXPECTED NORMAL VALUES.....	8
8	QUALITY CONTROL .....	9
9	PERFORMANCE CHARACTERISTICS.....	9
10	LIMITATIONS OF USE .....	12
11	LEGAL ASPECTS.....	12
12	REFERENCES / LITERATURE.....	12
1	EINLEITUNG.....	13
2	TESTPRINZIP.....	13
3	VORSICHTSMAßNAHMEN .....	14
4	BESTANDTEILE DES KITS .....	15
5	PROBENVORBEREITUNG .....	17
6	TESTDURCHFÜHRUNG .....	17
7	ERWARTETE WERTE.....	19
8	QUALITÄTS-KONTROLLE .....	20
9	ASSAY CHARACTERISTIKA.....	20
10	GRENZEN DES TESTS.....	20
11	RECHTLICHE GRUNDLAGEN .....	21
12	REFERENZEN / LITERATUR .....	21
1	INTRODUZIONE .....	22
2	PRINCIPIO DEL TEST .....	22
3	PRECAUZIONI.....	22
4	COMPONENTI DEL KIT .....	23
5	CAMPIONI .....	25
6	ATTUAZIONE DEL TEST .....	25
7	VALORI NORMALI.....	27
8	CONTROLLO QUALITÀ.....	28
9	CARATTERISTICHE DEL TEST.....	28
10	LIMITAZIONE DEL TEST.....	28
11	ASPETTI LEGALI.....	29
12	BIBLIOGRAFIA .....	29
	SYMBOLS USED.....	30

## 1 INTRODUCTION

### 1.1 Intended Use

The **DRG 25-OH Vitamin D (total) ELISA** is an enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of total 25-OH Vitamin D (Vitamin D<sub>2</sub> and Vitamin D<sub>3</sub>) **in serum**.

### 1.2 Summary and Explanation

Vitamin D is a steroid hormone involved in the intestinal absorption of calcium and the regulation of calcium homeostasis. The two major forms of Vitamin D, named Vitamin D<sub>3</sub> (cholecalciferol) and Vitamin D<sub>2</sub> (ergocalciferol), have isomeric structures, but D<sub>2</sub> is supposed to be less active than D<sub>3</sub><sup>1</sup>.

Physiological Vitamin D<sub>3</sub> levels result not only from dietary uptake but can also be produced from a cholesterol precursor, 7-dehydrocholesterol, in the skin during sun exposure. D<sub>2</sub> is obtained from plant sources and only represents less than 5% of the total Vitamin D in the body<sup>2</sup>. In the liver, the Vitamin D is hydroxylated to 25-hydroxyvitamin D (25-OH D), the major circulating metabolite of Vitamin D. Vitamin D and 25-OH D enter the circulation bound to Vitamin D binding protein (VDBP). Upon request, a small portion of 25-OH D is further hydroxylated in the kidney to form the biologically active hormone 1,25 dihydroxyvitamin D (1,25-(OH)<sub>2</sub> D)<sup>3</sup>. This process is tightly regulated by the concentration of 1,25-(OH)<sub>2</sub> D, parathyroid hormone, hypophosphatemia and ionized calcium levels. Concentrations of 1,25-(OH)<sub>2</sub> are about 1000-fold lower than that of 25-OH D<sup>4</sup>. Although 1,25-(OH)<sub>2</sub> D portrays the biological active form of Vitamin D, it is widely accepted that the measurement of circulating 25-OH D provides better information with respect to patients Vitamin D status and allows its use in diagnosis of hypovitaminosis<sup>5</sup>.

The concentration of 25-OH D decreases during winter time (reduced sun exposure), with dark skin colour and with age<sup>6</sup>. Determination of 25-OH D in serum or plasma will support the diagnosis and therapy control of postmenopausal osteoporosis, rickets in children, osteomalacia, renal osteodystrophy, neonatal hypocalcemia and hyperparathyroidism. In addition, the effects of prevailing subclinical Vitamin D deficiency in different European countries is critically discussed<sup>6</sup>. Vitamin D intoxication mostly occurs during a large intake of pharmaceutical preparations of Vitamin D and may lead to hypercalcemia and nephrocalcinosis in susceptible infants.

## 2 PRINCIPLE OF THE TEST

The **DRG 25-OH Vitamin D total ELISA Kit** is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding.

In the first step, samples have to be pretreated in separate vials with denaturation buffer to extract the analyte, since most circulating 25-OH Vitamin D is bound to VDBP *in vivo*. After neutralization, biotinylated 25-OH Vitamin D (enzyme conjugate) and peroxidase-labeled streptavidin- (enzyme complex) are added.

After careful mixing, the solution is transferred to the wells of the microtiter plate. Endogenous 25-OH Vitamin D of a patient sample competes with a 25-OH Vitamin-D<sub>3</sub>-biotin conjugate for binding to the VDBG that is immobilized on the plate. Binding of 25-OH Vitamin D-biotin is detected by peroxidase-labeled streptavidin. Incubation is followed by a washing step to remove unbound components. The color reaction is started by addition of enzyme substrate and stopped after a defined time. The colour intensity is inversely proportional to the concentration of 25-OH Vitamin D in the sample.

**3 WARNINGS AND PRECAUTIONS**

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with *Stop Solution* containing 0.5 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

## 4 REAGENTS

### 4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;  
Wells coated with Vitamin D binding protein (VDBG).
2. **Standard (Standard 0-5)**, 6 vials, 1 mL, ready to use;  
Concentrations: 0 – 4 – 10 – 25 – 60 – 130 ng/mL  
Conversion: 1 ng/mL = 2.5 nmol/L.  
*The standards are calibrated against the NIST Standard Reference Material (SRM) 2972;*  
Contain non-mercury preservative.
3. **Control Low & High**, 2 vials, 1 mL each, ready to use;  
For control values and ranges please refer to vial label or QC-Datasheet.  
Contain non-mercury preservative.
4. **Denaturation Buffer**, 1 vial, 10 mL, ready to use
5. **Neutralization Buffer**, 1 vial, 25 mL, ready to use,  
Contains non-mercury preservative.
6. **Enzyme Conjugate**, 1 vial, 7 mL, ready to use,  
Vitamin D<sub>3</sub> conjugated to biotin;  
Contains non-mercury preservative.
7. **Enzyme Complex**, 1 vial, 7 mL, ready to use,  
Streptavidin-Peroxidase Conjugate  
Contains non-mercury preservative.
8. **Substrate Solution**, 1 vial, 25 mL, ready to use,  
Tetramethylbenzidine (TMB).
9. **Stop Solution**, 1 vial, 14 mL, ready to use,  
contains 0.5 M H<sub>2</sub>SO<sub>4</sub>.  
Avoid contact with the stop solution. It may cause skin irritations and burns.
10. **Wash Solution**, 1 vial, 30 mL (40X concentrated),  
see „Preparation of Reagents“.
11. **Cover Foil**, 1 sheet of adhesive cover foil

**Note:** Additional *Standard 0* for sample dilution is available upon request.

### 4.2 Materials required but not provided

- Vials for the Vitamin D release step  
(e.g. rack with 96 vials, 0.65 mL each, 5 pieces; [REF](#) 781565-5; (BRAND GmbH))
- Single vials (e.g. single vials 0.65 mL each, 1000 pieces ([REF](#): EIA-5396-VIALS))
- Incubator 37 °C (98.6 °F)
- A microtiter plate calibrated reader (450 ± 10 nm)
- Calibrated variable precision micropipettes.
- Absorbent paper.
- Distilled or deionized water
- Timer
- Semi logarithmic or linear graph paper or software for data reduction

### 4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for two months if stored as described above.

#### 4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature prior to use.

##### **Wash Solution**

Add deionized water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL deionized water to a final volume of 1200 mL.

*The diluted Wash Solution is stable for 2 weeks at room temperature.*

##### **Working Conjugate Solution**

Prepare an adequate volume of Working Conjugate Solution by mixing *Enzyme Conjugate* with *Enzyme Complex* in a 1:1 ratio **at least 30 minutes before use**.

Shorter pre-incubation of this conjugate-complex-mixture will lead to increased Vitamin D concentrations.

*Stability of the prepared Working Conjugate Solution: 1 week at 2 °C to 8 °C in a sealed container.*

##### Example:

If the whole plate is used, mix 6 mL of *Enzyme Conjugate* with 6 mL of *Enzyme Complex* to a total volume of 12 mL.

If the whole plate is not used at once prepare the required quantity of *Working Conjugate Solution* by mixing 0.5 mL of *Enzyme Conjugate* with 0.5 mL of *Enzyme Complex* per strip (see table below):

No. of strips	<i>Enzyme Conjugate</i> (mL)	<i>Enzyme Complex</i> (mL)
1	0.5	0.5
2	1.0	1.0
3	1.5	1.5
4	2.0	2.0
5	2.5	2.5
6	3.0	3.0
7	3.5	3.5
8	4.0	4.0
9	4.5	4.5
10	5.0	5.0
11	5.5	5.5
12	6.0	6.0

#### 4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet.

#### 4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

## 5 SPECIMEN COLLECTION AND PREPARATION

**Serum** can be used in this assay.

Do not use haemolytic, icteric or lipaemic specimens.

Please note: Samples containing sodium azide should not be used in the assay.

### 5.1 Specimen Collection

#### **Serum:**

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

### 5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 3 days at 2 °C to 8 °C prior to assaying.

Specimens held for a longer time (up to two months) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

### 5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Standard 0* and reassayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

#### Example:

- a) dilution 1:10:           10 µL sample + 90 µL *Standard 0* (mix thoroughly)
- b) dilution 1:100:        10 µL dilution a) 1:10 + 90 µL *Standard 0* (mix thoroughly).

## 6 ASSAY PROCEDURE

### 6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

## 6.2 Test Procedure

Each run must include a standard curve.

All standards, samples, and controls should be run in duplicate. All standards, samples, and controls should be run concurrently so that all conditions of testing are the same.

### 6.2.1 Release Procedure and Pretreatment

1. Prepare an adequate volume of **Working Conjugate Solution** (see chapter 4.4)
2. Secure the desired number of appropriate vials or uncoated plates for the Vitamin D release step (not included in the kit).
3. Dispense **25 µL** of each **Standard, Control** and **sample** with new disposable tips into the vials.
4. Dispense **50 µL Denaturation Buffer** into each vial.
5. Seal vials and incubate for **30 minutes** at 37 °C.
6. Add **200 µL of Neutralization Buffer** to each vial.
7. Add **100 µL of Working Conjugate Solution** to each vial.
8. Thoroughly mix for 10 seconds. It is important to have a complete mixing of the solution in this step. Use 150 µL of this mixed solution for the ELISA.

### 6.2.2 ELISA Procedure

1. Secure the desired number of Microtiter wells in the frame holder.
2. Transfer **150 µL** of the mixed solution of each **Standard, Control** and **sample** with new disposable tips into the appropriate wells.  
(E.g. if in step 6.2.1 a rack with 96 vials/wells is used, an 8-channel pipette can be used for the transfer.)
3. Seal wells carefully and incubate for **60 minutes** at 37 °C.
4. Briskly shake out the contents of the wells.  
Rinse the wells **4 times** with diluted *Wash Solution* (300 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.  
**Important note:**  
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
5. Add **150 µL of Substrate Solution** to each well.
6. Incubate for **15 minutes** at room temperature.
7. Stop the enzymatic reaction by adding **100 µL of Stop Solution** to each well.
8. Determine the absorbance (OD) of each well at **450 ± 10 nm** with a microtiter plate reader.  
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

## 6.3 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Using semi-logarithmic or linear graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 130 ng/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

### 6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Units (450 nm)
Standard 0 (0 ng/mL)	2.13
Standard 1 (4 ng/mL)	1.91
Standard 2 (10 ng/mL)	1.63
Standard 3 (25 ng/mL)	1.12
Standard 4 (60 ng/mL)	0.57
Standard 5 (130 ng/mL)	0.24

## 7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently normal healthy adults, using the DRG 25-OH Vitamin D ELISA the following values were observed:

Population	Valid N	Age (years)	Mean Age	Mean Conc. (ng/mL)	5. Percentile (ng/mL)	95. Percentile (ng/mL)
Total	120	21 - 76	44	23.7	8.69	51.0
Males	64	24 - 76	47	25.8	9.78	51.9
Females	56	21 - 75	40	21.2	7.70	41.3
Caucasian	61	21 - 75	44	28.2	12.6	52.5
Hispanic	24	21 - 76	40	23.2	12.1	40.8
Afro-American	35	24 - 65	45	16.1	5.95	36.5
Northern	40	21 - 75	43	20.3	11.8	32.3
Central	40	24 - 62	43	16.1	6.36	29.6
Southern	40	24 - 76	45	34.6	18.9	56.0
Summer	60	24 - 76	47	29.7	12.7	52.7
Winter	60	21 - 66	40	17.6	6.67	29.5

The samples were collected from subjects with different skin tones, at 3 different geographical locations (North, South and Central United States), during summer and winter time.

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

A review of the literature suggests the following ranges for the classification of 25-OH Vitamin D status<sup>3</sup>:

Vitamin D status	25-OH Vitamin D (ng/mL)	25-OH Vitamin D (nmol/L)
Deficiency	< 10	< 25
Insufficiency	10 – 29	25 – 72.5
Sufficiency	30 – 100	75 – 250
Toxicity	> 100	> 250

## 8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs like the international Vitamin D Quality Assessment Scheme (DEQAS) in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

## 9 PERFORMANCE CHARACTERISTICS

### 9.1 Assay Dynamic Range

The range of the assay is between 3.5 – 130 ng/mL.

### 9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross reactivity of the assay:

25-OH Vitamin D <sub>3</sub> :	102.4 %
25-OH Vitamin D <sub>2</sub> :	69.5 %
1,25 (OH) <sub>2</sub> Vitamin D <sub>3</sub> :	< 0.1 %
1,25 (OH) <sub>2</sub> Vitamin D <sub>2</sub> :	< 0.1 %
3-Epi-25-OH-Vitamin D <sub>3</sub> :	66.3 %
Vitamin D <sub>3</sub> :	3.8 %
Vitamin D <sub>2</sub> :	3.2 %

### 9.3 Sensitivity

The analytical sensitivity of the DRG ELISA was calculated according to CLSI EP17-A by subtracting 1.645 standard deviations from the mean of 60 replicate analyses of the Zero Standard (S0) and was found to be < 2.5 ng/mL.

The functional sensitivity of the DRG ELISA was evaluated according to CLSI EP17-A, is defined as the Vitamin D concentration that approximates a CV of 20% and was found to be < 3.5 ng/mL.

### 9.4 Reproducibility

Reproducibility of the DRG ELISA was evaluated according to CLSI EP5-A. Six samples, containing different concentrations of analyte, were assayed in duplicate in two assays per day over 20 operating days to determine the within- and between-assay variability.

#### 9.4.1 Intra Assay

The within assay variability is shown below:

Sample	n	Mean (ng/mL)	CV (%)
1	80	6.0	4.4
2	80	10.9	3.3
3	80	22.5	1.7
4	80	29.2	1.6
5	80	57.3	1.6
6	80	88.9	1.4

### 9.4.2 Inter Assay

The between assay variability is shown below:

Sample	n	Mean (ng/mL)	CV (%)
1	80	6.0	8.3
2	80	10.9	5.2
3	80	22.5	4.1
4	80	29.2	3.6
5	80	57.3	3.2
6	80	88.9	2.1

### 9.4.3 Inter-Lot

The inter-assay (between-lots) variation was determined by repeated measurements of 6 samples in 3 different kit lots.

Sample	n	Mean (ng/mL) Lot 1	Mean (ng/mL) Lot 2	Mean (ng/mL) Lot 3	CV (%)
1	18	7,3	7,6	6,4	8,7
2	18	16,8	15,7	14,9	6,2
3	18	33,0	30,7	29,4	5,9
4	18	52,7	50,1	47,1	5,6
5	18	80,9	74,3	74,1	5,1
6	18	108,8	105,7	102,3	3,1

### 9.5 Recovery

Samples have been spiked by adding 4 solutions with known concentrations in a 1:1 ratio.

The % recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100 (expected value = (endogenous value + added value) / 2; because of a 1:2 dilution of serum with spike material).


	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	
<b>Concentration [ng/mL]</b>	9.0	17.4	36.4	60.5	82.6	
<b>Average Recovery</b>	94.7	93.9	104.1	105.2	92.5	
<b>Range of Recovery [%]</b>	from	91.4	92.3	92.9	96.9	88.2
	to	99.7	96.7	112.4	111.3	94.1

### 9.6 Linearity

According to CLSI EP6-A samples have been diluted with Zero Standard (S0). The results were analyzed as a linear regression of the Expected against Observed values. The resulting regression equation is:

Observed = Expected 1.01 + 0.30; R<sup>2</sup> = 0.997

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	
<b>Concentration [ng/mL]</b>	69.4	88.7	104.8	107.4	128.3	
<b>Average Recovery</b>	90.8	98.4	107.5	103.2	102.7	
<b>Range of Recovery [%]</b>	from	85.8	94.7	105.7	97.9	97.5
	to	97.8	103.3	109.2	108.5	107.8


### 9.7 Comparison Studies

For method comparison according to CLSI EP09-A2-IR, 145 samples spanning the assay range were tested by the DRG 25-OH Vitamin D (total) ELISA and by the DiaSorin LIAISON® 25-OH Vitamin D TOTAL Assay.

The resulting regression equation was:

DRG = DiaSorin  $0.99 + 1.72$ ;  $R = 0.967$ .

## 10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

### 10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL), Triglyceride (up to 7.5 mg/mL), Cholesterol (up to 2.8 mg/mL), Albumin (up to 75 mg/mL) and Rheumatoid Factor (up to 100 U/mL) have no influence on the assay results.

### 10.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of 25-OH Vitamin D in a sample.

### 10.3 High-Dose-Hook Effect

No hook effect was observed in this test.

## 11 LEGAL ASPECTS

**Only for countries where the declaration of European Conformity (CE mark) is applicable.**

### 11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

### 11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

## 12 REFERENCES / LITERATURE

1. Armas LAG., Hollis M., Heaney RP. Vitamin D<sub>2</sub> is much less effective than vitamin D<sub>3</sub> in humans. *J. Clin. Endocrinol. Metab.* 2004; 89(11) 5387-91.
2. Houghton LA., Vieth R. The case against ergocalciferol (vitamin D<sub>2</sub>) as a vitamin supplement. *Am. J. Nutr.* 2006; 84, 694-97.
3. Holick M. Vitamin D: the underappreciated D-lightful hormone that is important for skeletal and cellular health. *Curr. Opin. Endocrinol. Diabetes* 2002; 9(1) 87-98.
4. Pilz S. et al. Vitamin D: clinical implications beyond musculoskeletal diseases. *J. Lab. Med.* 2011; 35(4) 211-16.
5. Visser M. et al. Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *Am. J. Clin. Nutr.* 2006; 84(3) 616-22.
6. Souberbielle JC. Et al. Vitamin D and musculoskeletal health, cardiovascular disease, autoimmunity and cancer: Recommendations for clinical practice. *Autoimmun Rev.* 2010; 9 709-15.

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