Instructions for Use

Leptospira Ab vet ELISA

VET

REF EIA-6072

Σ 96

DRG Instruments GmbH, Germany
Frauenbergstraße. 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la versión válida de la metodico técnico incluido aquí en el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.

Table of Contents / Inhaltsverzeichnis / Tabla de Contenidos

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1 INTRODUCTION

Leptospirosis is probably the most widespread zoonosis in the world. It is caused by infection with spirochete bacteria of the genus *Leptospira* and affects humans as well as a broad spectrum of animal hosts (e.g. dogs, pigs, cattle). The incidence is significantly higher in warm climate countries than in temperate regions. The disease is seasonal, with peak incidence occurring in summer or fall in temperate regions, where temperature is the limiting factor in survival of leptospires, and during rainy seasons in warm climate regions, where rapid desiccation would otherwise prevent survival. Natural reservoirs for the pathogenic *Leptospira interrogans* include rodents as well as a large variety of domesticated mammals. Leptospires occupy the lumen of nephritic tubules in their natural host and are shed into the urine.

Transmission can occur when humans or animals are directly or indirectly exposed to the urine of infected animals or a urine-polluted environment. Leptospires gain entry into the blood stream via cuts, skin abrasions or mucous membranes through contact with moist soil, vegetation, contaminated water and also by handling infected animal tissues or ingestion of food and water. The incubation period is usually 5-14 days, with a range of 2-30 days. The spectrum of clinical symptoms is extremely wide. The vast majority of leptospiral infections are either subclinical or result in very mild illness and recover without any complications. Clinical manifestations of leptospirosis range from mild influenza-like symptoms to severe life-threatening disease forms, characterized by jaundice, renal failure, bleeding and severe pulmonary hemorrhage. Acute kidney injury (AKI) is the most commonly recognized disease in dogs, accounting for more than 90% of reported cases of leptospirosis. Hepatic disease occurs concurrently in 10% - 20% of dogs with AKI. *Leptospira pomona* and *Leptospira hardjo* cause a febrile syndrome, loss of appetite, and, in dairy cows, a large drop in - even the suppression of - milk secretion for a period ranging from two to ten days. *L. pomona* causes important reproductive problems in female breeding pigs spreading slowly through the herd. *L. tarassovi* causes a similar syndrome as *L. pomona* but tends to be milder and to spread more slowly. The clinical presentation of leptospirosis is biphasic, with the acute or septicemic phase lasting about a week, followed by the immune phase, characterized by antibody production and excretion of leptospires in the urine. Most of the complications of leptospirosis are associated with localization of leptospires within the tissues during the immune phase.

The presence of pathogen resp. infection may be identified by

- Pathogen detection: dark-field microscopy of a culture from blood, urine, cerebrospinal fluid or tissue PCR
- Serology: microscopic agglutination test (MAT), ELISA

2 INTENDED USE

The Leptospira Ab vet ELISA is intended for the qualitative determination of antibodies against Leptospira in **veterinary serum**.

3 PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplates reader.

4 MATERIALS

4.1 Reagents supplied

- **Microtiterplate:**
12 break-apart 8-well snap-off strips coated with Leptospira antigens; in resealable aluminium foil.
- **Sample Dilution Buffer:**
1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015 % (v/v) CMIT/MIT (3:1).
- **Stop Solution:**
1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):**
1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:**
1 bottle containing 20 mL of peroxidase Protein A/G; coloured yellow, ready to use; white cap; ≤ 0.02 % (v/v) MIT.
- **TMB Substrate Solution:**
1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap
- **Positive Control**:**
1 vial containing 2 mL; coloured yellow; ready to use; red cap; ≤ 0.02 % (v/v) MIT.
- **Cut-off Control**:**
1 vial containing 3 mL; coloured yellow; ready to use; green cap; ≤ 0.02 % (v/v) MIT.
- **Negative Control**:**
1 vial containing 2 mL; coloured yellow; ready to use; blue cap; ≤ 0.0015 % (v/v) CMIT/MIT (3:1).

For hazard and precautionary statements see 12.1

For potential hazardous substances please check the safety data sheet.

4.2 Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3 Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplate wells
- Pipettes to deliver volumes between 10 µL and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes
- Timer

5 STABILITY AND STORAGE

Store the kit at 2 °C - 8 °C.

The opened reagents are stable up to the expiry date stated on the label when stored at 2 °C - 8 °C.

6 REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20 °C - 25 °C) and mix them before starting the test run!

6.1 Microtiterplate

The break-apart snap-off strips are coated with Leptospira antigens.

Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2 °C - 8 °C.

6.2 Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water.

The diluted buffer is stable for 5 days at room temperature (20 °C - 25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

6.3 TMB Substrate Solution

The reagent is ready to use and has to be stored at 2 °C - 8 °C, away from the light.

The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7 SPECIMEN COLLECTION AND PREPARATION

Use **canine serum** samples with this assay.

If the assay is performed within 5 days after sample collection, the samples should be kept at 2 °C - 8 °C; otherwise they should be aliquoted and stored deep-frozen (-70 °C to -20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

Heat inactivation of samples is not recommended.

7.1 Sample Dilution

Before assaying all samples should be diluted **1 + 100** with Sample Dilution Buffer.

Dispense 10 µL sample and 1 mL Sample Dilution Buffer into tubes to obtain a 1 + 100 dilution and thoroughly mix with a Vortex.

8 ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 °C ± 1 °C.

1. Dispense 100 µL controls and diluted samples into their respective wells.
Leave well A1 for the substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 C ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the substrate blank well (A1).
6. **Incubate for 30 min at room temperature (20 °C - 25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells
9. **Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.1 Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

9 RESULTS

9.1 Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate blank** Absorbance value < 0.100
- **Negative Control** Absorbance value < 0.200 and < Cut-off
- **Cut-off Control** Absorbance value 0.150 – 1.300
- **Positive control** Absorbance value > Cut-off.

If these criteria are not met, the test is not valid and must be repeated.

9.2 Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43
Cut-off = 0.43

9.2.1 Results in Units [DU]

Sample (mean) absorbance value x 10 = [DRG Units = DU]
Cut-off

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ DU}$

9.3 Interpretation of Results

Normal value ranges for this ELISA should be established by each laboratory based on its own sample populations in the geographical areas serviced.

The following values should be considered as a guideline:

Cut-off	10 DU	-
Positive	> 11 DU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 DU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 DU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.		

10 SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

The performance data have been established with canine samples. Due to the nature of the Protein A/G conjugate this ELISA should react with other mammalian species also. More detailed information is available on request.

10.1 Precision

Intra-assay	n	Mean (E)	CV (%)
#1	24	0.176	11.78
#2	24	1.237	7.50
#3	24	0.324	8.47
Inter-assay	n	Mean (DU)	CV (%)
#1	12	3.63	3.68
#2	12	10.89	4.93
#3	12	30.98	3.85

10.2 Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte.

Diagnostic Specificity canine: 100.0 % (95% confidence interval: 59.04 % - 100.0 %)

10.3 Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte.

Diagnostic Sensitivity canine: 88.89 % (95% confidence interval: 51.75 % - 99.72 %)

10.4 Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

10.5 Cross-Reactivity

Cross-reactions cannot be excluded.

11 LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12 PRECAUTIONS AND WARNINGS

- Only for veterinary in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

12.1 Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

	<p>H317 May cause an allergic skin reaction. P261 Avoid breathing spray P280 Wear protective gloves/protective clothing. P302+P352 IF ON SKIN: Wash with plenty of soap and water. P333+P313 If skin irritation or rash occurs: Get medical advice/attention. P362+P364 Take off contaminated and Wash it before reuse.</p>
---	--

Further information can be found in the safety data sheet.

12.2 Disposal Considerations

1 RESIDUES OF CHEMICALS AND PREPARATIONS ARE GENERALLY CONSIDERED AS HAZARDOUS WASTE. THE DISPOSAL OF THIS KIND OF WASTE IS REGULATED THROUGH NATIONAL AND REGIONAL LAWS AND REGULATIONS. CONTACT YOUR LOCAL AUTHORITIES OR WASTE MANAGEMENT COMPANIES WHICH WILL GIVE ADVICE ON HOW TO DISPOSE HAZARDOUS WASTE.

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13 BIBLIOGRAPHY / LITERATUR / BIBLIOGRAFÍA

1. Adler, B.; Murphy, A.M.; Locarnini, S.A.; Faine, S. Detection of specific anti-leptospiral immunoglobulins M and G in human serum by solid-phase enzyme-linked immunosorbent assay. *J. Clin. Microbiol.* May 1980 vol. 11 no. 5 452-457.
2. Ahmad, S.N.; Shah, S.; Ahmad, F.M. Laboratory diagnosis of leptospirosis. *J Postgrad Med* 2005;51:195-200.
3. Cruz, L.S.; Vargas, R.; Lopes, A.A. Leptospirosis: a worldwide resurgent zoonosis and important cause of acute renal failure and death in developing nations. *Ethn Dis.* 2009 Spring;19(1 Suppl 1):S1-37-41.
4. Doungchawee, G.; Sirawaraporn, W.; Icksang-Ko, A.; Kongtim, S.; Naigowit, P.; Thongboonkerd, V. Use of immunoblotting as an alternative method for serogrouping *Leptospira*. *J Med Microbiol.* 2007 May;56(Pt 5):587-592.
5. Jansen, A.; Stark, K. (2006): Leptospirose (Kap. VIII - 1.25). In: F. Hofmann (Hrsg), *Handbuch der Infektionskrankheiten: Epidemiologie, Diagnostik, Therapie, Prophylaxe, Gesetzliche Regelungen*, 2. Aufl., 17. Erg.Lfg. 11/06. Landsberg/Lech: EcoMed, S. 1-13.
6. Johnson, R. C. (2001) *Leptospira*. In: *Medical Microbiology*. edited by Baron, S. The University of Texas Medical Branch, Levett, P.N. *Leptospirosis*. *Clin Microbiol Rev.* 2001 Apr;14(2):296-326.
7. Sanders, E.J.; Rigau-Pérez, J.G.; Smits, H.L.; Deseda, C.C.; Vorndam, V.A.; Aye, T.; Spiegel, R.A.; Weyant, R.S.; Bragg, S.L. Increase of leptospirosis in dengue-negative patients after a hurricane in Puerto Rico in 1996 [correction of 1966]. *Am J Trop Med Hyg.* 1999 Sep;61(3):399-404.
8. Terpstra, W.J.; Ligthart, G.S.; Schoone, G.J. ELISA for the detection of specific IgM and IgM in human leptospirosis. *J Gen Microbiol.* 1985 Feb;131(2):377-385.
9. Utzinger, J.; Becker, S.L.; Knopp, S.; Blum, J.; Neumayr, A.L.; Keiser, J.; Hatz, C.F. Neglected tropical diseases: diagnosis, clinical management, treatment and control. *Swiss Med Wkly.* 2012 Nov 22;142:w13727.
10. Vijayachari P, Sugunan AP, Shriram AN. Leptospirosis: an emerging global public health problem. *J Biosci.* 2008 Nov;33(4):557-69.
11. World Health Organization. Human leptospirosis: guidance for diagnosis, surveillance and control. WHO (2003)
12. Major A, Schweighauser A, Francy T. Increasing incidence of canine leptospirosis in Switzerland. I. *Int J Environ Res Public Health.* 2014 Jul 16;11(7):7242-60. doi: 10.3390/ijerph110707242.
13. Alton GD, Berke O, Reid-Smith R, et al. Increase in seroprevalence of canine leptospirosis and its risk factors, Ontario 1998–2006. *Can J Vet Res.* 2009;73:167–175.
14. Sykes JE, Hartmann K, Lunn KF, Moore GE, Stoddard RA, Goldstein RE. 2010 ACVIM small animal consensus statement on leptospirosis: diagnosis, epidemiology, treatment, and prevention. *J Vet Intern Med.* 2011;25:1–13.

Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas

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14 SCHEME OF THE ASSAY**Test Preparation**

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all samples and standards/controls on the plate layout supplied in the kit.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate Blank (A1)	Negative Control	Cut-off Control	Positive Control	Sample (diluted 1+100)
Negative Control	-	100 µL	-	-	-
Cut-off Control	-	-	100 µL	-	-
Positive Control	-	-	-	100 µL	-
Sample (diluted 1+100)	-	-	-	-	100 µL
Cover wells with foil supplied in the kit Incubate for 1 h at 37±1 °C					
Wash each well three times with 300 µL of Washing Buffer					
Conjugate	-	100 µL	100 µL	100 µL	100 µL
Incubate for 30 min at room temperature (20 °C - 25 °C) Do not expose to direct sunlight Wash each well three times with 300 µL of Washing Buffer					
TMB Substrate Solution	100 µL	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark					
Stop Solution	100 µL	100 µL	100 µL	100 µL	100 µL
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnóstico in vitro	Producto sanitario para diagnóstico In vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	Numero di Catalogo	Número de catálogo	Référence de catalogue
	Batch code *	Chargencode *	Codice del lotto	Codigo de lote	Numéro de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservacion	Température de conservation
	Use-by date *	Verwendbar bis *	Utilizzare prima del	Establa hasta	Utiliser jusque
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Caution *	Achtung *			
	For veterinary in vitro diagnostic use only				
Distributed by	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
Content	Content	Inhalt	Contenuto	Contenido	Contenu
Volume/No.	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité
MTP	Microtiterplate	Mikrotiterplatte	Piastre di Microtitolazione	Placa de Microtitulación	Plaque de Microtitrage
CONJL	Conjugate	Konjugat	Coniugato	Conjugado	Conjugué
CONTROL -	Negative Control	Negativkontrolle	Controllo Negativo	Control Negativo	Contrôle Négatif
CONTROL +	Positive Control	Positivkontrolle	Controllo Positivo	Control Positivo	Contrôle Positif
CUT OFF	Cut off Control	Cut off-Kontrolle	Controllo Cut-off	Control Cut-off	Contrôle Cut-off
CAL	Standard or Calibrator	Standard oder Kalibrator	Standard o Calibratore	Estándar o Calibrador	Standard o Etalon
DIL	Sample Dilution Buffer	Probenverdünnungspuffer	Tampone di Diluizione del Campione	Tampón de Dilución de Muestras	Tampon de Dilution d'Échantillon
SOLN STOP	Stop solution	Stoplösung	Soluzione Bloccante	Solución de Parada	Solution d'Arrêt
SUB TMB	TMB Substrate solution	TMB-Substratlösung	Soluzione Substrato TMB	Solución de Sustrato de TMB	Solution de Substrat TMB
WASH BUF 20x	Washing Buffer 20x concentrated	Waschpuffer 20x konzentriert	Tampone di Lavaggio concentrazione x20	Tampón de Lavado concentrado x20	Tampon de Lavage concentré 20 x