Instructions for Use

Bordetella pertussis Toxin (PT) IgA ELISA

IVD

CE

REF EIA-6132

96

Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la versión válida de la metodico técnico incluido aquí en el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos

1	INTRODUCTION.....	3
2	INTENDED USE.....	3
3	PRINCIPLE OF THE ASSAY	3
4	MATERIALS.....	4
5	STABILITY AND STORAGE	4
6	REAGENT PREPARATION	4
7	SAMPLE COLLECTION AND PREPARATION.....	5
8	ASSAY PROCEDURE	5
9	RESULTS.....	6
10	SPECIFIC PERFORMANCE CHARACTERISTICS	7
11	LIMITATIONS OF THE PROCEDURE	8
12	PRECAUTIONS AND WARNINGS	8

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1 INTRODUCTION

Bordetella pertussis is a respiratory pathogen that causes pertussis, commonly known as whooping cough, a localized infection of the ciliated epithelium of the bronchial tree. Pertussis is characterized by a prolonged paroxysmal cough often accompanied by an inspiratory whoop.

The disease affects mainly children, but adults have also been increasingly reported to be affected. The pathogen produces toxins which cause local damage to the cilia of epithelial cells, which leads to prolonged illness and pertussis. Disease presentation varies with age and history of previous exposure or vaccination. Severe disease is infrequent in healthy, vaccinated persons. Infants, particularly those who have not received the primary vaccination series against pertussis, are at risk for complications and mortality.

In addition to *B. pertussis*, three other *Bordetella* species can cause disease in humans: *B. parapertussis*, *B. holmesii*, and *B. bronchiseptica*. *B. parapertussis* causes a pertussis-like illness that is generally milder than pertussis because the bacteria do not produce pertussis toxin. Co-infection of *B. pertussis* and *B. parapertussis* is not unusual.

B. pertussis is of worldwide prevalence. Globally, 20 - 40 million cases of pertussis occur each year, 90 % of which are in developing countries, and there are up to 400 000 fatalities each year, mostly in young infants.

Transmission of *B. pertussis* occurs primarily via close direct contact with an infected person or inhalation of airborne droplets. Symptoms develop following inhalation of the airborne pathogen. The organism is highly contagious, with up to 90 % of household contacts developing the disease. Infected persons are most contagious in the catarrhal and the paroxysmal stages.

The incubation period is usually seven to 10 days, with a range of 4 - 21 days.

Typical pertussis symptoms occur in three different stages: catarrhal, paroxysmal, and convalescent.

The catarrhal stage lasts for about 1 - 2 weeks, and is characterized by non-specific symptoms such as rhinorrhea, sneezing, low-grade fever and cough. The second stage is the paroxysmal stage, lasting for about 4 - 6 weeks, and is characterized by various pathognomonic symptoms of pertussis such as episodes of paroxysmal cough with a characteristic whooping sound. The final stage is the convalescent stage. During this stage, the respiratory symptoms gradually decrease although cough can continue for months.

Many factors can alter the usual course of pertussis, causing an atypical presentation. Previously vaccinated adolescents and adults may have less severe paroxysmal symptoms.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Bordetella pertussis</i>	Pertussis whooping cough	<p><u>1. Stadium catarrhale:</u> symptoms of a cold with slight fever (1-2 weeks)</p> <p><u>2. Stadium convulsivum:</u> severe, spasmodic coughing; after deep inspiration follows a coughing staccato (2-6 weeks)</p> <p><u>3. Stadium decrementi:</u> Ease of disease with symptoms of a bronchitis (up to 6 weeks)</p>	Highly contagious droplet infection

Infection or presence of pathogen may be identified by

- Microscopy
- PCR
- Serology: e. g. by ELISA

2 INTENDED USE

The *Bordetella pertussis* Toxin (PT) IgA ELISA is intended for the quantitative determination of IgA class antibodies against *B. pertussis* toxin (PT) antigens in human serum or plasma (citrate, heparin).

3 PRINCIPLE OF THE ASSAY

The quantitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

4 MATERIALS

4.1 Reagents supplied

- **Microtiterplate:**
12 break apart 8-well snap-off strips coated with B. pertussis toxin (PT) antigens; in resealable aluminium foil.
- **IgA Sample Dilution Buffer:**
1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2 ; coloured yellow; ready to use; white cap; $\leq 0.0015\% \text{ (v/v)}$ CMIT/ MIT (3:1). **Stop Solution:**
1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):**
1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2 , for washing the wells; white cap.
- **Conjugate:**
1 bottle containing 20 mL of peroxidase labelled antibody to human IgA in phosphate buffer (10 mM); coloured violet; ready to use; black cap.
- **TMB Substrate Solution:**
1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap
- **Standards:** 4 vials, each containing 2 mL standard; coloured yellow; ready to use; yellow cap. $\leq 0.02\% \text{ (v/v)}$ MIT.
 Standard A: 0 IU/mL
 Standard B: 10 IU/mL
 Standard C: 25 IU/mL
 Standard D: 50 IU/mL

The standards are calibrated in accordance with the "1st WHO International Standard Pertussis antiserum (human)", Code 06/140, of the National Institute for Biological Standards and Control (NIBSC), Potters Bar, UK.

• Control Low:

1 vial containing 2 mL control; coloured yellow; ready to use; blue cap.
Concentration in IU/mL is indicated on the label; $\leq 0.02\% \text{ (v/v)}$ MIT.

• Control High:

1 vial containing 2 mL control; coloured yellow; ready to use; red cap.
Concentration in IU/mL is indicated on the label; $\leq 0.02\% \text{ (v/v)}$ MIT.

For hazard and precautionary statements see 12.1

For potential hazardous substances please check the safety data sheet.

4.2 Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3 Materials and Equipment needed

- ELISA Microtiterplate plate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 µL and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

5 STABILITY AND STORAGE

Store the kit at 2 °C - 8 °C.

The opened reagents are stable up to the expiry date stated on the label when stored at 2 °C - 8 °C.

6 REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20 °C - 25 °C) and mix them before starting the test run!

6.1 Microtiterplate

The break-apart snap-off strips are coated with *B. pertussis* toxin (PT) antigen.

Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2 °C - 8 °C.

6.2 Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water.

The diluted buffer is stable for 5 days at room temperature (20 °C - 25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

6.3 TMB Substrate Solution

The reagent is ready to use and has to be stored at 2 °C - 8 °C, away from the light.

The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7 SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2 °C - 8 °C; otherwise, they should be aliquoted and stored deep-frozen (-70 °C to -20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

Heat inactivation of samples is not recommended.

7.1 Sample Dilution

Before assaying, all samples should be diluted **1+100** with IgA Sample Dilution Buffer.

Dispense 10 µL sample and 1 mL IgA Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8 ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 °C ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 °C ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20 °C - 25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.1 Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

9 RESULTS

9.1 Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value < **0.100**
- **Standard A:** Absorbance value < **0.200**
- **Standard B** Absorbance value > **Standard A**
- **Standard C** Absorbance value > **Standard B**
- **Standard D:** Absorbance value > **1.000**
- **Control Low:** Result in IU/mL within range indicated on the label
- **Control High:** Result in IU/mL within range indicated on the label

Standard A < Standard B < Standard C < Standard D

If these criteria are not met, the test is not valid and must be repeated.


9.2 Calculation of Results

In order to obtain **quantitative results in IU/mL** plot the (mean) absorbance values of the 4 Standards A - D on (linear/Linear) graph paper in a system of coordinates against their corresponding concentrations (0, 10, 25 and 50 IU/mL) and draw a standard curve (absorbance values on the y-axis, concentrations on the x-axis).

Read results from this calibration curve employing the (mean) absorbance values of each patient sample and control.

For the calculation of the standard-curve mathematical Point to Point function should be used.

9.3 Typical Standard Curve

9.4 Interpretation of Results

According to recent literature and recommendations from reference laboratories the following interpretation of results is recommended:

IgG		IgA
< 40 IU/mL	not indicative of recent contact	

$\geq 40 - < 100$ IU/mL	testing of a second sample after 7-10 days or further analysis of the IgA antibody response is recommended	< 12 IU/mL	not indicative of recent contact
		≥ 12 IU/mL	indicative of recent contact
≥ 100 IU/mL	indicative of recent contact		
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.			

9.4.1 Antibody Isotypes and State of Infection

Serology	Significance
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection
IgA	Produced in mucosal linings throughout the body (⇒ protective barrier) Usually produced early in the course of the infection

10 SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact DRG.

10.1 Precision

Intra assay	n	Mean (E)	CV (%)
#1	23	0.547	8.01
#2	24	1.954	2.03
#3	24	1.048	2.66

Inter assay	n	Mean (IU/mL)	CV (%)
#1	12	5.19	7.28
#2	12	16.29	4.37
#3	12	6.23	8.17

10.2 Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 90.91% (95% confidence interval: 75.67% - 98.08%).

10.3 Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100.0% (95% confidence interval: 69.15% - 100.0%).

10.4 Analytical Sensitivity

The analytical sensitivity (according to CLSI EP17-A) is defined as the apparent concentration of the analyte that can be distinguished from the zero calibrator. It is 1.56 IU/mL.

10.5 Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

10.6 Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters (including antibodies to several respiratory pathogens) did not reveal evidence of false-positive results due to cross-reactions.

10.7 Measurement range

The measurement range is 1.56 IU/mL - 50 IU/mL.

11 LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12 PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following standards of good laboratory practice(GLP).
- For further internal quality control each laboratory should additionally use known samples.

12.1 Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning	H317	May cause an allergic skin reaction.
	P261	Avoid breathing spray
	P280	Wear protective gloves/ protective clothing.
	P302+P352	IF ON SKIN: Wash with plenty of soap and water.
	P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
	P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet

12.2 Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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13 BIBLIOGRAPHY / LITERATUR / BIBLIOGRAFIA / BIBLIOGRAFÍA

1. Baughman, Andrew L.; Bisgard, Kristine M.; Edwards, Kathryn M.; Guris, Dalya; Decker, Michael D.; Holland, Kathy et al. (2004): Establishment of diagnostic cutoff points for levels of serum antibodies to pertussis toxin, filamentous hemagglutinin, and fimbriae in adolescents and adults in the United States.
In *Clinical and diagnostic laboratory immunology* 11 (6), pp. 1045–1053. DOI: 10.1128/CDLI.11.6.1045-1053.2004.
2. Guiso, N.; Berbers, Guy A.; Fry, N. K.; He, Q.; Riffelmann, M.; Wirsing von König, Carl Heinz (2011): What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories.
In *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology* 30 (3), pp. 307–312. DOI: 10.1007/s10096-010-1104-y.
3. Le, Thuan; Cherry, James D.; Chang, Swei-Ju; Knoll, Maria Deloria; Lee, Martin L.; Barenkamp, Steve et al. (2004): Immune responses and antibody decay after immunization of adolescents and adults with an acellular pertussis vaccine: the APERT Study. In *The Journal of Infectious Diseases* 190 (3), pp. 535–544. DOI: 10.1086/422035.
4. Melker, H. E. de; Versteegh, F. G.; Conyn-Van Spaendonck, M. A.E.; Elvers, L. H.; Berbers, Guy A.; van der Zee, A.; Schellekens, J. F. (2000): Specificity and sensitivity of high levels of immunoglobulin G antibodies against pertussis toxin in a single serum sample for diagnosis of infection with Bordetella pertussis.
In *Journal of Clinical Microbiology* 38 (2), pp. 800–806.
5. Pebody, R. G.; GAY, N. J.; Giannamico, A.; Baron, S.; Schellekens, J.; Tischer, A. et al. (2005): The seroepidemiology of Bordetella pertussis infection in Western Europe.
In *Epidemiology and infection* 133 (1), pp. 159–171. DOI: 10.1017/S0950268804003012.
6. Poynten, I. M.; Hanlon, M.; Irwig, L.; Gilbert, G. L. (2002): Serological diagnosis of pertussis. Evaluation of IgA against whole cell and specific Bordetella pertussis antigens as markers of recent infection.
In *Epidemiol. Infect.* 128 (02). DOI: 10.1017/S0950268801006598.
7. Riffelmann, M.; Thiel, K.; Schmetz, J.; Wirsing von König, Carl Heinz (2010): Performance of commercial enzyme-linked immunosorbent assays for detection of antibodies to Bordetella pertussis.
In *Journal of Clinical Microbiology* 48 (12), pp. 4459–4463. DOI: 10.1128/JCM.01371-10.
8. Riffelmann, Marion; Littmann, Martina; Hülße, Christel; Hellenbrand, Wiebke; Wirsing von König, Carl Heinz (2008): Pertussis: not only a disease of childhood.
In *Deutsches Arzteblatt international* 105 (37), pp. 623–628. DOI: 10.3238/arztebl.2008.0623.
9. Robert Koch Institut (RKI) (2013): Pertussis (Keuchhusten). In *RKI-Ratgeber für Ärzte*.
10. Ward, Joel I.; Cherry, James D.; Chang, Swei-Ju; Partridge, Susan; Keitel, Wendy; Treanor, John et al. (2006): Bordetella Pertussis Infections in Vaccinated and Unvaccinated Adolescents and Adults, as Assessed in a National Prospective Randomized Acellular Pertussis Vaccine Trial (APERT).
In *Clinical Infectious Diseases* 43, pp. 151–157.
11. Wong, Emily J.; Hewlett, Erik L. (2006): Pertussis. In Richard L. Guerrant, David H. Walker, Peter F. Weller (Eds.): Tropical infectious diseases. Principles, pathogens & practice. 2nd ed. Philadelphia: Churchill Livingstone, 369-373.
12. Xing, Dorothy; Markey, Kevin; Newland, Penny; Rigsby, Peter; Hockley, Jason; He, Qiushui (2011): EUVAC.NET collaborative study: evaluation and standardisation of serology for diagnosis of pertussis.
In *Journal of immunological methods* 372 (1-2), pp. 137–145. DOI: 10.1016/j.jim.2011.07.005.
13. Xing, Dorothy; Wirsing von König, Carl Heinz; Newland, Penny; Riffelmann, Marion; Meade, Bruce D.; Corbel, Michael J.; Gaines-Das, Rose (2009): Characterization of reference materials for human antiserum to pertussis antigens by an international collaborative study.
In *Clinical and vaccine immunology : CVI* 16 (3), pp. 303–311. DOI: 10.1128/CVI.00372-08.

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14 SCHEME OF THE ASSAY**Test Preparation**

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all samples and standards/controls on the plate layout supplied in the kit.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate Blank (A1)	Standards A - D	Control Low	Control High	Sample (diluted 1+100)
Standards A - D	-	100 µL	-	-	-
Control Low	-	-	100 µL	-	-
Control High	-	-	-	100 µL	-
Sample (diluted 1+100)	-	-	-	-	100 µL
Cover wells with foil supplied in the kit Incubate for 1 h at 37 °C Wash each well three times with 300 µL of Washing Buffer					
Conjugate	-	100 µL	100 µL	100 µL	100 µL
Incubate for 30 min at room temperature (20 °C - 25 °C) Do not expose to direct sunlight Wash each well three times with 300 µL of Washing Buffer					
TMB Substrate solution	100 µL	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark					
Stop Solution	100 µL	100 µL	100 µL	100 µL	100 µL
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français	Português
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes	Conformidade Europeia
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation	Consultar as instruções de uso
	In vitro diagnostic medical device *	In-vitro-Diagnostikum *	Diagnistica in vitro	Diagnóstico in vitro	Diagnostic in vitro	Dispositivo médico para diagnóstico in vitro
	Catalogue number *	Katalognummer *	No. di Cat.	No de catálogo	Référence	Número de catálogo
	Batch code *	Chargen-bezeichnung *	Lotto no	Número de lote	No. de lot	Código do lote
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests	Suficiente para <n> determinações
	Temperature limit *	Temperaturgrenzwerte *	Temperatura di conservazione	Temperatura de conservacion	Temperature de conservation	Limites de temperatura
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation	Prazo de validade
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant	Fabricante
	Distributor *	Vertriebspartner *	Distributore	Distribuidor	Distributeur	Distribuidor
	Date of manufacture *	Herstellungsdatum *	Data di produzione	Fecha de fabricación	Date de production	Data de fabricação
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques	Riscos biológicos
	Caution *	Achtung *	Attenzione	Precaución	Attention	Cuidado
	Unique device Identifier *	eindeutige Produktidentifizierung *	Identificativo unico del dispositivo*	Identificación exclusiva del dispositivo *	Identifiant de dispositif unique*	Identificador único do dispositivo *
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches	Apenas para uso em investigação
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement	Conteúdo
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité	Volume / Quantidade
<i>MTP</i>	Microplate	Mikrotiterplatte	Micrioplastra	Microplaca	Microplaqué	Microplaca
<i>CONJ</i>	Conjugate	Konjugat	Coniugato	Conjugado	Conjugué	Conjugado
<i>CONTROL -</i>	Negative Control	Negativkontrolle	Controllo negativo	Control positivo	Contrôle négatif	Controlo Negativo
<i>CONTROL +</i>	Positive Control	Positivkontrolle	Controllo positivo	Control positivo	Contrôle positif	Controlo positivo
<i>CUT OFF</i>	Cut off control	Cut off-Kontrolle	Controllo cut-off	Control cut-off	Contrôle cut-off	Controlo de corte
<i>CAL</i>	Standard or Calibrator	Standard oder Kalibrator	Standard o Calibratore	Estándar o Calibrador	Standard o Etalon	Padrão ou Calibrador
<i>DIL A</i>	IgA Sample Diluent	IgA-Proben-verdünnungspuffer	Tampon diluente IgA	Diluyente para IgA de la muestra	Diluant pour échantillon IgA	Diluente de amostra IgA
<i>DIL G</i>	IgG Sample Diluent	IgG-Proben-verdünnungspuffer	Tampon diluente IgG	Diluyente para IgG de la muestra	Diluant pour échantillon IgG	Diluente de Amostra IgG
<i>DIL M</i>	IgM Sample Diluent	IgM-Proben-verdünnungspuffer	Tampon diluente IgM	Diluyente para IgM de la muestra	Diluant pour échantillon IgM	Diluente de amostra IgM
<i>SOLN STOP</i>	Stop solution	Stoplösung	Soluzione bloccante	Solución de parada	Solution d'arrêt	Solução de paragem
<i>SUB TMB</i>	TMB Substrate solution	TMB-Substratlösung	Soluzione substrato TMB	Solución de sustrato de TMB	Solution de substrat TMB	Solução de substrato TMB
<i>WASH BUF 20x</i>	Washing Buffer 20x concentrated	Waschpuffer 20x konzentriert	Tampon di lavaggio concentrazione x20	Tampón de lavado concentrado x20	Tampon de lavage concentré 20 x	Tampão de lavagem 20x concentrado