# DRG:HYBRID•XL<sup>®</sup>

Instructions for Use

## Calprotectin

IVD

CE

REF

HYE-5767

40 tests


**DRG Instruments GmbH**, Germany  
Frauenbergstraße. 18, D-35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

**DRG International, Inc.**, USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**ASSAY PROTOCOL BARCODE (APB)**

(Version 2.61 Software or later / ab Softwareversion 2.61 oder höher / Versione Software 2.61 o superiore /  
Versión de Software 2.61 o posterior / Version de logiciel 2.61 et supérieure /  
Wersja oprogramowania 2.61 lub późniejsza)

**HYE-5767 – v02.70.04**

The barcode must be used to install the assay protocol into the DRG:HYBRiD-XL software via the SCAN NEW LOT page.  
Der Barcode muss in dem Menü „NEUES LOT SCANNEN“ eingelesen werden, um das Protokoll in der DRG:HYBRiD-XL-Software zu installieren.

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Please refer to section **3: Routine Procedures: "Installing a new assay product"** of the User Manual v2.60 or later.

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## 1 INTRODUCTION

### 1.1 Intended Use

The **DRG:HYBRiD-XL Calprotectin** is an enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of Calprotectin in stool.

Only for use with the DRG:HYBRiD-XL Analyzer.

### 1.2 Summary and Explanation

Calprotectin is a calcium-binding protein secreted predominantly by neutrophils and monocytes.

Fecal calprotectin is a marker for acute inflammation in general and for neoplastic and inflammatory gastrointestinal diseases (1-3). It is often difficult to distinguish between irritable bowel syndrome and chronic inflammatory bowel disease. This can lead to extensive and unnecessary colonoscopic examinations. In contrast, the fecal calprotectin test allows clear differentiation between patients with inflammatory bowel disease (acute Morbus Crohn 's disease and ulcerative colitis) and irritable bowel syndrome (4-9). Fecal calprotectin levels correlate significantly with histologic and endoscopic assessment of disease activity in Morbus Crohn 's disease and ulcerative colitis as well as with the fecal excretion of indium-111-labelled neutrophilic granulocytes that has been suggested as the "gold standard" of disease activity in inflammatory bowel disease (10). Elevated levels of calprotectin are a much better predictor of relapse of inflammatory bowel disease than standard inflammatory markers (11). Comparing this marker with standard fecal occult blood screening in colorectal cancer demonstrates clearly the diagnostic advantages of the fecal calprotectin test. The parameter is of a high diagnostic value: if the calprotectin level in stool is low, the probability is high that no organic intestinal disease exists.

## 2 PRINCIPLE OF THE TEST

The DRG:HYBRiD-XL Calprotectin Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **sandwich principle**.

The antibody coated wells (ACW) of the reagent cartridges are coated with a monoclonal (mouse) antibody directed towards a unique antigenic site of the Calprotectin molecule.

An aliquot of sample containing endogenous Calprotectin is incubated in the coated well with enzyme conjugate, which is an anti-Calprotectin antibody conjugated with horseradish peroxidase.

After incubation unbound substances are washed off.

Having added the substrate solution, the intensity of color developed is proportional to the concentration of Calprotectin in the patient sample. The concentrations of unknown samples are determined from the master curve.

## 3 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only.  
For professional use only.
- This kit can only be used in combination with the DRG:HYBRiD-XL Analyzer
- Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is clear and understood.
- Do not remove, exchange, discard or damage any of the barcode labels provided with each kit and its components. All barcodes build an integral system for the kit lot.**
- Respect the general safety measures for use of laboratory reagents.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear appropriate disposable gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may cause false results.
- Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Unused reagent cartridges must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant provided.
- Optimal test results are only obtained when using calibrated pipettes.

- Do not mix or use components from kits with different lot numbers. It is advised not to interchange reagent cartridges of different kits even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the wells in the reagent cartridges may differ slightly.
- Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
- TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant amount of water and skin with soap and plenty of water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.
- For information on hazardous substances included in the kit please refer to Safety Data Sheets.  
For professional users the Safety Data Sheet for this product is available upon request directly from DRG.

## 4 REAGENTS

### 4.1 Reagents provided

#### 4.1.1 Reagent Cartridges

40 pieces containing the following:

- Antibody Coated Well (ACW)**  
coated with anti-Calprotectin antibody (monoclonal).
- Enzyme Conjugate**, 240 µL,  
Anti-Calprotectin antibody conjugated with horseradish peroxidase;  
Contains non-mercury preservative.
- Substrate Solution**, 260 µL  
Tetramethylbenzidine (TMB).

#### 4.1.2 Re-Calibrator 1 & 2

2 vials, 1 mL each, lyophilized, see "Reagent Preparation ";  
For re-calibration of the quantitative DRG:HYBRiD-XL Calprotectin test.

Concentrations are lot-specific.

Contain non-mercury preservative.

#### 4.1.3 Control 1 & 2

2 vials, 1 mL each lyophilized, see "Reagent Preparation ";  
For control values and ranges please refer to the vial label or to the CoA.

Contain non-mercury preservative.

### 4.2 Materials required but not provided

- General needed laboratory equipment
- Ultra-pure water  
DRG recommends to use Clinical Laboratory Reagent Water (CLRW) according to CLSI guideline 3C-A4 with the following specifications:  
Resistivity at 25 °C [MΩ·cm]: > 10  
Conductivity at 25 °C [µS·cm<sup>-1</sup>]: < 0.1  
Total Organic Carbon/p.p.b.[µg/L] : < 50  
Colloids [µg/mL]: <0.05
- REF** HYB-5252 DRG:HYBRiD-XL Analyzer
- REF** HYI-5392: *System Solution 5L*, 5000 mL;  
(Instrument Feed Water according to CLSI guideline 3C-A4 with the following specification can also be used:  
Resistivity at 25 °C [MΩ·cm]: > 1  
Conductivity at 25 °C [µS·cm<sup>-1</sup>]: < 1  
Total Organic Carbon/p.p.b [µg/L] : < 200  
Colloids [µg/mL]: <0.1)
- REF** HYI-6234: *Wash Buffer*, 10x concentrate, 100 mL
- REF** HYI-5395: *Needle Cleaning Solution*, 30 mL. Cleaning solution for the pipetting needle (daily and weekly maintenance, see also user manual)
- REF** HYI-5387: *Cuvettes*, (2 x 360 pieces)

For use of the *Secondary Sample Holder* for secondary tubes the following tubes are required:

- REF** HYI-5391: *Sample Tubes (Secondary)*, 2500 pcs.

### 4.3 Storage Conditions

All kit components should be stored at 2 °C to 8 °C to ensure product performance until the defined expiry date.

When stored at 2 °C to 8 °C, **unopened kits** will retain reactivity until expiration date. Do not use reagents beyond this date.

- Cartridges (stored at 2 °C to 8 °C) in the supplied and unopened zip/foil bags will retain reactivity until expiration date.
- Unopened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity until expiration date.

**Opened reagents and the reagent cartridges** must be stored at 2 °C to 8 °C.

Once the plastic bag has been opened, care should be taken to tightly close it again along with the supplied desiccant bag.

**Immediately after end of each** run the Re-Calibrator and Control vials have to be removed from the instrument, tightly capped and stored at 2 °C to 8 °C.

- Unused cartridges in opened zip/foil bags (stored at 2 °C to 8 °C) will retain reactivity until expiration date, if stored as described above.
  - Pierced or open cartridges must be disposed of immediately.
- Opened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity for 8 weeks.

#### 4.3.1 On-board Stability

For Re-Calibrators and Controls the on-board stability has been evaluated under controlled laboratory conditions at room temperature (20 °C to 25 °C).

Due to the differences in laboratory environmental conditions and reagent volumes, the on-board stability may deviate from the declared value.

On-board stability:	8 h
---------------------	-----

### 4.4 Reagent Preparation

Bring all reagents, such as controls and re-calibrators, to room temperature (20 °C to 25 °C) prior to use. Reagent Cartridges can be used directly after storage in the refrigerator.

#### Re-Calibrator 1 & 2

Reconstitute the lyophilized content of each vial with 1 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

Stability after reconstitution:	at 2 °C to 8 °C	1 week
	at -20 °C (in aliquots)	8 weeks

#### Control 1 & 2

Reconstitute the lyophilized content of each vial with 1 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

Stability after reconstitution:	at 2 °C to 8 °C	1 week
	at -20 °C (in aliquots)	8 weeks

#### Wash Buffer (10x) (HYI-6234; not included in the kit)

For Wash Buffer (1x) dilute 100 mL of Wash Buffer (10x) with 900 mL ultra-pure water to a final volume of 1000 mL.

stability after dilution:	20 °C - 26 °C	4 weeks
---------------------------	---------------	---------

### 4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet.

### 4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

## 5 SPECIMEN COLLECTION AND PREPARATION

Stool samples can be used in this assay.

**A minimum of 80 µL of sample is needed for one determination.** This includes 20 µL sample and 60 µL dead volume.

- Do not use standards or external controls stabilized with azide.

### 5.1 Specimen Collection

We **strongly** recommend the use of the **Fecal Extraction System 2.6 (FES; REF: HYI-5989)**.

**Only the Fecal Extraction System 2.6 (FES; REF: HYI-5989)** can be used for stool sample preparation in the DRG:HYBRID-XL Calprotectin test.

No other stool sample extraction system has been tested for the DRG Calprotectin test (HYE-5767).

For sample preparation please follow the Instructions for Use included in the package of HYI-5989.

After finishing the sample preparation in the FES tubes the stool samples are ready for use in the DRG:HYBRID-XL Calprotectin test.

**Place the opened extraction tube into the adapter for the primary arc.**

### 5.2 Specimen Storage

#### Raw Stool

Calprotectin in stool is described to be stable for at least 3 days at room temperature (8,17,18).

Nevertheless, we recommend storing the samples for no more than 48 h at 2 °C - 8 °C.

Long term storage up to 12 months is recommended at -20 °C.

Allow frozen samples to thaw slowly, preferably at 2 °C - 8 °C, and warm the samples to room temperature before analysis. Avoid repeated freezing and thawing of the sample.

Freezing can cause neutrophil granulocytes in the stool sample to burst and release calprotectin. Therefore frozen samples can be expected to contain slightly elevated concentrations of calprotectin compared to fresh samples.

#### Stool Extract

Stool extract is stable for seven days at 2 °C - 8 °C or up to 12 months at -20 °C.

Avoid more than one freeze-thaw cycles.

### 5.3 Specimen Dilution

#### 5.3.1 Manual Sample Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the **pretreated stool sample (1:200 dilution in FES, see chapter 5.1)** can be diluted with *Sample Diluent\** and measured again as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

a) dilution 1:10: 10 µL sample + 90 µL *Sample Diluent* (mix thoroughly)

\* *Sample Diluent for manual dilution is not included in this kit, but can be ordered on request (REF HYE-5767-DIL, 20 mL).*

## 6 ASSAY PROCEDURE

### 6.1 General Remarks

- All reagents, such as controls and re-calibrators and specimens, must be allowed to come to room temperature (20 °C to 25 °C) before use.

All reagents and samples must be mixed without foaming.

Reagent Cartridges can be used directly after storage in the refrigerator.

- Samples, controls and re-calibrators should be measured within 2 hours in order to avoid possible evaporation effects.
- The *Secondary Sample Holder* (HYI-5437) for secondary tubes has the capacity for a maximum of 20 samples including controls and re-calibrators. They all have to be pipetted into the secondary tubes, and the respective barcodes of control/re-calibrator vials and, if available, the sample barcodes have to be read with the external barcode scanner.

## 6.2 Test Procedure

- The total assay time for DRG:HYBRiD-XL Calprotectin is 60 minutes.
- To ensure proper operation of the test, the instructions in the user manual for the DRG:HYBRiD-XL should strictly be followed.
- All test specific information required for the correct operation is included in the respective barcodes of the reagents.  
**Take care not to damage these bar codes!**
- It is recommended to tap the bottom of the Cartridge Segments containing the reagent cartridges once on the bench before placing them on the rotor. This is to avoid foam and adhering of the liquid on the sealing of the reagent cartridge.
- Place reagent cartridges on the rotor of the unit. The heating to 37 °C incubation temperature is performed automatically in the unit.

## 6.3 Calibration

Each DRG:HYBRiD-XL reagent contains a barcode with the specific information for recalibration of the reagent lot. The Master Curve is printed as a 2-D barcode on the outer label of the kit package and on the QC-Datasheet and has to be scanned with the external barcode scanner prior to the first use of the respective kit lot.

### Recalibration is recommended:

- if a new kit lot is used. Each new lot should be verified by running the kit internal re-calibrators and controls before routine use.
- if one or both assay controls are found outside the specified range.
- after 1 week of use of the same reagent kit on the unit.

## 6.4 Calculation of Results

The analyte concentrations are calculated automatically by the DRG:HYBRiD-XL's system software.

## 7 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results.

It is also recommended to participate in national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analysing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following: expiration dates and storage conditions of reagents, operational reliability of the analyser. In addition, it is indicated to perform a Recalibration.

In case of further questions please contact your local distributor or DRG directly.

### 7.1 Internal Controls

For Quality Control it is necessary to use the two internal controls provided with each kit.

Acceptance ranges for both internal controls (*Control 1 & 2*) were established by the manufacturer and are summarized in the QC-Datasheet added to the kit. Note that the expected values and acceptance ranges stated in the QC-Datasheet always refer to the current kit lot.

Internal controls should be run in single determination:

- on a routine basis (e.g. once per 24 h)
- if re-calibration is required (if one or both internal controls are out of range)
- if a new kit lot is used (in order to avoid any negative impact on the kit performance by improper transport or to detect improper storage during transport).

### 7.2 External Controls

Use controls at both normal and pathological levels.

The control intervals and control ranges for external controls should be adapted to the individual requirements of each laboratory. All results must be within the defined limits.

Each laboratory should establish corrective measures to be taken if values of external controls are not found in the acceptance range.

## 8 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and pathological values.

Study 1: In a clinical study, faecal Calprotectin concentration was determined in 302 apparently healthy children (12).

Age	n	Median faecal Calprotectin (µg/g)
0 - 3 months	14	345 (195 - 621)
3 - 6 months	13	278 (85 - 988)
6 - 12 months	27	183 (109 - 418)
1 - 4 years	89	75 (53 - 119)
4 - 12 years	159	28 (25 - 35)

Study 2: In this clinical study, faecal Calprotectin concentration was determined in 117 apparently healthy children (13).

Age	n	Median faecal Calprotectin (µg/g)
4 - 6 years	27	28.2
7 - 10 years	30	13.5
11 - 14 years	27	9.9
15 - 17 years	33	14.6

These values correspond well to other reported clinical values (14-15). A suggested cut-off level for adults (50 µg/g) can be used for children aged 4 - 17 years.

In a study conducted with apparently healthy adults using the DRG:HYBRiD-XL HYE-5767 the following data were observed:

Population	healthy adults
n	50
Range (µg/g)	2.6 - 103
Mean (µg/g)	17.78
2.5 <sup>th</sup> - 97.5 <sup>th</sup> Percentile (µg/g)	4.19 - 71.87
Median (µg/g)	25.27

### 8.1 Interpretation of results

- The median value in healthy adults is about 18 µg/g (5).
- Samples with a calprotectin concentration < 50 µg/g are regarded as negative.
- Samples with a calprotectin concentration between 50 µg/g and 200 µg/g are regarded as borderline positive. We recommend repeating the measurement at a later time point in order to confirm the result.
- Samples with a calprotectin concentration > 200 µg/g are regarded as positive.

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated with other clinical observations and diagnostic tests.

## 9 LIMITATIONS OF USE

Reliable and reproducible results will be obtained, when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

## 10 PERFORMANCE CHARACTERISTICS

### 10.1 Assay Dynamic Range

The dynamic range of the assay is defined by the limit of detection and the maximum value of the Master Curve.

Values found below the measuring range are indicated as "< 4.502 µg/g".

Values found above the measuring range are indicated as "> 1500 µg/g".

The measuring range of the assay is between 4.50 µg/g - 1500 µg/g.

**10.2 Specificity of Antibodies (Cross-Reactivity)**

The following substances were tested for cross-reactivity of the assay:

Lysozyme	0%
PMN-Elastase	0%
Myeloperoxidase	0%
Lactoferrin	0%

**10.3 Sensitivity**

The sensitivity study was designed according to CLSI guideline EP17-A2.

The Limit of Blank (LoB) is 1.260 µg/g.

The Limit of Detection (LoD) is 4.502 µg/g.

The Limit of Quantification (LoQ) is 11.048 µg/g.

**10.4 Precision Performance**

The precision study was designed on the basis of CLSI guideline EP5-A2.

**10.4.1 Within-Run (Intra Assay)**

The Intra-Assay precision was determined with 4 patient samples covering the measuring range in 5 independent runs on 5 days with 1 device in 5 replicates per run. CV was calculated as mean CV of 5 runs.

Sample	n	Mean (µg/g)	CV (%)
1	5	26.53	4.4
2	5	87.73	4.1
3	5	156.30	5.6
4	5	359.39	7.5

**10.4.2 Total Precision (Inter Assay)**

The Inter-Assay precision was determined with 4 patient samples covering the measuring range in 5 independent runs on 5 days with 5 determinations per run. CV was calculated from 25 determinations. The total precision is shown below:

Sample	n	Mean (µg/g)	CV (%)
1	25	26.53	6.5
2	25	87.73	5.3
3	25	156.30	6.7
4	25	359.39	7.6

**10.4.3 Inter-Lot**

The between-lot variation was determined by 6 measurements of 4 samples with 3 different kit lots.

Sample	n	Mean (µg/g)	CV (%)
1	18	31.46	5.6
2	18	155.15	3.0
3	18	327.76	4.3
4	18	910.25	10.6

**10.5 Recovery**

Recovery was determined by adding four increasing concentrations of the analyte to four different patient samples containing different amounts of endogenous analyte. Each sample (non-spiked and spiked) was assayed, and analyte concentrations of the samples were calculated from the Master Curve. The percentage recoveries were determined by comparing expected and measured values of the samples.

Sample	1	2	3	4	
Concentration (µg/g)	29.77	145.64	386.70	1146.49	
Average Recovery (%)	99.1	110.3	106.9	97.4	
Range of Recovery (%)	from	90.0	105.0	101.9	92.8
	to	109.1	113.1	110.8	101.7

**10.6 Linearity**

Four samples containing different amounts of analyte were serially diluted with *Sample Diluent*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

Sample	1	2	3	4	
Concentration (µg/g)	326.47	1199.79	1444.54	26.63	
Average Recovery (%)	100.9	98.1	96.6	89.2	
Range of Recovery (%)	from	95.8	90.8	86.5	85.5
	to	112.8	112.2	107.1	92.2

**10.7 Interfering Substances**

In general, hemolytic, icteric or lipemic samples should be avoided. Unconjugated bilirubin (up to 0.3 mg/mL), conjugated bilirubin (up to 0.1 mg/mL), hemoglobin (up to 4 mg/mL) and triglyceride (up to 7.5 mg/mL) have no influence on the assay results.


**10.8 High-Dose-Hook Effect**

Hook effect was not observed in this test up to a concentration of 6,000 µg/g of Calprotectin.

**11 METHOD COMPARISON**

A comparison of DRG:HYBRiD-XL Calprotectin (HYE-5767) (y) and the reference method (x) using clinical samples gave the following correlation:

n = 76  
 r = 0.91  
 y = 1.058x + 12.68

## 12 LEGAL ASPECTS

**Only for countries where the declaration of European Conformity (CE mark) is applicable.**

### 12.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

### 12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different kit lots could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

**For further information please refer to the User Manual of the DRG:HYBRiD-XL, analyser-specific application sheets, product information and package inserts of all necessary components.**

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1. Fagerhol MK. Calprotectin, a faecal marker of organic gastrointestinal abnormality. *Lancet*. 2000;356(9244):1783-1784. doi:10.1016/S0140-6736(00)03224-4.
2. van Rheenen PF et al. Faecal calprotectin for screening of patients with suspected inflammatory bowel disease: diagnostic meta-analysis. *BMJ*. 2010;341:c3369. doi:10.1136/bmj.c3369.
3. Poullis A et al. Review article: faecal markers in the assessment of activity in inflammatory bowel disease. *Aliment Pharmacol Ther*. 2002;16(4):675-681. doi:10.1046/j.1365-2036.2002.01196.x.
4. Konikoff MR, Denson LA. Role of fecal calprotectin as a biomarker of intestinal inflammation in inflammatory bowel disease. *Inflamm Bowel Dis*. 2006;12(6):524-534. doi:10.1097/00054725-200606000-00013.
5. D'Haens G et al. Fecal calprotectin is a surrogate marker for endoscopic lesions in inflammatory bowel disease. *Inflamm Bowel Dis*. 2012;18(12):2218-2224. doi:10.1002/ibd.22917.
6. Önal IK et al. The value of fecal calprotectin as a marker of intestinal inflammation in patients with ulcerative colitis. *Turk J Gastroenterol*. 2012;23(5):509-514. doi:10.4318/tjg.2012.0421.
7. Kalantari H et al. Fecal calprotectin is a useful marker to diagnose ulcerative colitis from irritable bowel syndrome. *Adv Biomed Res*. 2015; 4:85. doi:10.4103/2277-9175.156647.
8. von Arnim U et al. Fecal calprotectin: a marker for clinical differentiation of microscopic colitis and irritable bowel syndrome. *Clin Exp Gastroenterol*. 2016; 9:97-103. doi:10.2147/CEG.S97701.
9. Wright EK et al. Comparison of Fecal Inflammatory Markers in Crohn's Disease. *Inflamm Bowel Dis*. 2016;22(5):1086-1094. doi:10.1097/MIB.0000000000000671.
10. Røseth AG et al. Correlation between faecal excretion of indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. *Scand J Gastroenterol*. 1999;34(1):50-54. doi:10.1080/00365529950172835.
11. Tibble JA et al. Surrogate markers of intestinal inflammation are predictive of relapse in patients with inflammatory bowel disease. *Gastroenterology*. 2000;119(1):15-22. doi:10.1053/gast.2000.8523.
12. Hestvik, E. et al. Faecal calprotectin concentrations in apparently healthy children aged 0-12 years in urban Kampala, Uganda: a community-based survey. *BMC Pediatr*. 2011; 11:9. doi:10.1186/1471-2431-11-9.
13. Fagerberg UL et al. Fecal calprotectin levels in healthy children studied with an improved assay. *J Pediatr Gastroenterol Nutr*. 2003;37(4):468-472. doi:10.1097/00005176-200310000-00013.
14. Li F et al. Fecal calprotectin concentrations in healthy children aged 1-18 months. *PLoS One*. 2015;10(3):e0119574. doi: 10.1371/journal.pone.0119574.
15. Krzesiek E. Fecal Calprotectin as an Activity Marker of Inflammatory Bowel Disease in Children. *Adv Clin Exp Med*. 2015;24(5):815-822. doi:10.17219/acem/26003.
16. Bohra A et al. The synergy of dual faecal immunochemical and faecal calprotectin testing for accurate assessment of endoscopic and histological activity in inflammatory bowel disease. *Therapeutic advances in gastroenterology* vol. 17. 2024, doi:10.1177/17562848241237895.
17. Hamer HM et al. Impact of Preanalytical Factors on Calprotectin Concentration in Stool: A Multiassay Comparison. *J Appl Lab Med*. 2022;7(6):1401-1411. doi:10.1093/jalm/jfac057.
18. Delia A et al. Fecal Calprotectin: A Comparison of Two Commercial Enzymoimmunoassays and Study of Fecal Extract Stability at Room Temperature. *J Clin Med Res*. 2018;10(5):396-404. doi:10.14740/jocmr3275e.
19. Ghazal et al. Hormone Immunoassay Interference: A 2021 Update. *Ann Lab Med*. 2022;42(1):3-23. doi:10.3343/alm.2022.42.1.3.

## SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français	Polski
	European Conformity	CE-Konformitätskennzeichnung	Conformità europea	Conformidad europea	Conforme aux normes européennes	Zgodność z normami europejskimi
	Consult instructions for use*	Gebrauchsanweisung beachten*	Consultare le istruzioni per l'uso	Consulte las Instrucciones	Consulter les instructions d'utilisation	Zapoznać się z instrukcją użytkowania
	<i>In vitro</i> diagnostic medical device*	<i>In-vitro</i> -Diagnostikum*	Per uso Diagnostica in vitro	Para uso Diagnóstico in vitro	Usage Diagnostic in vitro	Wyrób medyczny do diagnostyki in vitro
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches	Tylko do użytku w badaniach
	Catalogue number*	Artikelnummer*	No. di Cat.	Número de catálogo	Référence	Numer katalogowy
	Batch code*	Chargencode*	Lotto no	Número de lote	No. de lot	Numer LOT
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests	Zawartość przeznaczona na <n> testów
	Temperature limit*	Temperaturbegrenzung*	Temperatura di conservazione	Temperatura de conservación	Température de conservation	Przechowywać w temperaturze
	Use-by date*	Verwendbar bis*	Data di scadenza	Fecha de caducidad	Date limite d'utilisation	Data przydatności
	Manufacturer*	Hersteller*	Fabbricante	Fabricante	Fabricant	Producent
	Caution*	Achtung*				UWAGA
	Distributor	Vertreiber	Distributore	Distribuidor	Distributeur	Dystrybutor
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu	Zawartość
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Numéro	Objętość / Numer
<i>Reagent Cartridge</i>	Reagent Cartridge	Reagenzien-Cartridge	Cartucce di reagente	Cartucho de reactivo	Cartouche de réactif	Kartridż z odczynnikami
<i>Re-Calibrator</i>	Re-Calibrator	Re-Kalibrator	Re-Calibratore	Re-Calibrador	Re-calibrateurs	Re-kalibrator
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle	Kontrola
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Tracciante enzimatico	Conjugado enzimático	Conjugué enzymatique	Koniugat enzymatyczny
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complesso enzimatico	Complex enzimático		Enzymu kompleks
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat	Roztwór Substratu
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungsmedium	Diluente dei campioni	Solución para dilución de la muestra	Diluant d'échantillon	Rozcieńczalnik Próbkki
<i>Reagent 1</i>	Reagent 1	Reagenz 1	Reagente 1	Reactivo 1	Réactif 1	Reagent 1
<i>Reagent 2</i>	Reagent 2	Reagenz 2	Reagente 2	Reactivo 2	Réactif 2	Reagent 2
<i>Wash Buffer</i>	Wash Buffer	Waschpuffer	Tampone di lavaggio	Tampón de lavado	Tampon de lavage	Bufor Płuczący
<i>System Solution</i>	System Solution	Systemlösung		Solución de sistema	Solution système	Roztwór systemowy
<i>Needle Cleaning Solution</i>	Needle Cleaning Solution	Nadel-Reinigungslösung	Soluzione lavaggio ago	Solución de lavado de la aguja	Solution de nettoyage des aiguilles	Roztwór Czyszczący Igły
<i>Denaturation Buffer</i>	Denaturation Buffer	Denaturierungspuffer	Tampone die denaturazione	Tampón de denaturalización	Tampon de dénaturation	
<i>Neutralization Buffer</i>	Neutralization Buffer	Neutralisierungspuffer	Tampone di neutralizzazione	Tampón de neutralización	Tampon de neutralisation	
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampone del test	Tampón de ensayo	Tampon d'essai	
<i>Secondary Sample Holder</i>	Secondary Sample Holder	Sekundärprobenhalter	Sostegno secondario dei campioni	Soporte de muestras para tubos secundarios	Support d'échantillons secondaires	Statyw na probówki drugorzędowe
<i>Sample Tubes (Secondary)</i>	Sample Tubes (Secondary)	Sekundärprobenröhrchen	Tubetti di campioni secondari	Tubos de muestra secundarios	Tubes échantillon secondaires	Probówki drugorzędowe
<i>Dilution Cartridge</i>	Dilution Cartridge				Cartouches de dilution	Kartridż do rozcieńczeń
<i>Vial Adapter</i>	Vial Adapter	Fläschchen-Adapter	l'adattatore dei tubetti	Adaptador de tubos	l'Adaptateur de flacons	Nakładki na butelki dla kontroli i kalibratorów