Instructions for Use

hGH IRMA (CT)

IVD

CE

REF RIA-3650

Σ 96

DRG Instruments GmbH, Germany
Frauenbergstraße 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos

1	INTENDED USE.....	3
2	CLINICAL BACKGROUND	3
3	PRINCIPLES OF THE METHOD	3
4	REAGENTS PROVIDED	3
5	SUPPLIES NOT PROVIDED	4
6	REAGENT PREPARATION	4
7	STORAGE AND EXPIRATION DATING OF REAGENTS	4
8	SPECIMEN COLLECTION AND PREPARATION.....	4
9	PROCEDURE	5
10	CALCULATION OF RESULTS.....	5
11	TYPICAL DATA.....	5
12	PERFORMANCE AND LIMITATIONS	6
13	LIMITATIONS.....	7
14	INTERNAL QUALITY CONTROL.....	7
15	REFERENCE INTERVALS	7
16	PRECAUTIONS AND WARNINGS	8

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1 INTENDED USE

Immunoradiometric assay kit for the invitro quantitative measurement of human growth hormone (hGH) in serum.

2 CLINICAL BACKGROUND

2.1 Biological activities of hGH

hGH is a polypeptide hormone (molecular weight 21,500 Da) produced by the acidophil cells of the anterior pituitary under the control of two main substances from the median eminence : Growth-hormone Releasing Factor (GRF) and an inhibitory agent, somatostatin. Dopaminergic, adrenergic and serotonergic neuroendocrine pathways also play an important role in the control of hGH secretion. Excitatory stimuli of hGH secretion include hypoglycemia, exercise, fasting, meals with a high protein content, deep sleep, stress, glucagon, L Dopa, amino acids, etc. Inhibitory stimuli include glucose, cortisol, hGH and free fatty acids. Because of its short plasma half-life (± 25 minutes) and of the frequent excitatory or inhibitory stimuli, hGH displays frequent and large variations of concentration in serum.

One of the main physiological functions of hGH is to act on the liver and other tissues to produce Somatomedins, which in turn induce growth by direct action on target tissues. In contrast to hGH, the concentration of somatomedin in serum is kept stable by virtue of being largely bound to circulating plasma proteins.

2.2 Clinical application of hGH-IRMA

Growth retardation

hGH hyposecretion is one of the various causes of small stature in children. Serum hGH measurement with a highly sensitive assay, especially following a provocative stimulus (absence of response), is an important way to establish this diagnosis because this group of patients can be treated by administration of hGH.

Hypopituitarism

Serum hGH measurement is also an index of pituitary function when hypopituitarism (either idiopathic or due to tumour and surgery) is suspected.

Gigantism and acromegaly

Serum hGH measurement, especially following a provocative inhibitory test (absence of response), is an important way to establish the diagnosis of hGH hypersecretion due to acidophilic pituitary tumour. This results in gigantism in children and acromegaly in adults. Both of these disorders may be treated by surgery or radiation.

3 PRINCIPLES OF THE METHOD

The hGH-IRMA is an immunoradiometric assay based on coated-tube separation. Mab 1 - the capture antibody is attached to the lower and inner surface of the plastic tube. Standards or samples added to the tubes will at first show low affinity for Mab 1. Addition of Mab 2, the signal antibody labeled with ^{125}I , will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration.

4 REAGENTS PROVIDED

Reagents		96 Test Kit	Reconstitution
Tubes coated with anti hGH (monoclonal antibodies)	TUBES	2 x 48	Ready for use
Tracer: ^{125}I odine labeled anti-hGH (monoclonal antibodies) in phosphate buffer with bovine serum albumin, azide (< 0.1%) and an inert red dye.	Ab ^{125}I	1 vial 5.5 mL 350 kBq	Ready for use
Zero calibrator (0 $\mu\text{IU}/\text{mL}$) in human plasma with thymol	CAL 0	1 vial lyophilized	Add 2 mL distilled water
Calibrators hGH N = 1 to 6 (see exact values on vial labels) in human plasma with thymol	CAL N	6 vials lyophilized	Add 0.5 mL distilled water
Wash Solution (NaCl, Tween 20)	WASH SOLN CONC	1 vial 40 mL	Dilute 20 x with distilled water (use a magnetic stirrer).
Controls - N = 1 or 2 in human plasma with thymol	CONTROL N	2 vials lyophilized	Add 0.5 mL distilled water

Note : Use the zero calibrator for sample dilutions.

1 μIU of the calibrator preparation is equivalent to 1 μIU NIBSC 2nd IS 98/574.

5 SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 50 µL, 500 µL and 2 mL (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. 5 mL automatic syringe (Cornwall type) for washing
6. Aspiration system (optional)
7. Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).

6 REAGENT PREPARATION

Calibrators:

Reconstitute the zero calibrator with 2.0 mL distilled water and other calibrators with 0.5 mL distilled water.

Controls:

Reconstitute the controls with 0.5 mL distilled water.

Working Wash solution:

Prepare an adequate volume of Working Wash solution by adding 19 volumes of distilled water to 1 volume of Wash Solution (20x). Use a magnetic stirrer to homogenize.

Discard unused Working Wash solution at the end of the day.

7 STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 °C to 8 °C.
- After reconstitution, calibrators and controls are stable for one week at 2 °C to 8 °C. For longer storage periods, aliquots should be made and kept at -20 °C for 3 months. Avoid successive freezing and thawing.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well closed vial at 2 °C to 8 °C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

8 SPECIMEN COLLECTION AND PREPARATION

Serum samples must be kept at 2 °C to 8 °C.

If the test is not run within 24 hours, storage at -20 °C is recommended.

Avoid successive freezing and thawing.

9 PROCEDURE

9.1 Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.

High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.

Prepare a calibration curve for each run; do not use data from previous runs.

9.2 Procedure

1. Label coated tubes in duplicate for each calibrator, sample, and control.
For the determination of total counts, label 2 normal tubes
2. Briefly vortex calibrators, controls and samples and dispense 50 µL of each into respective tubes.
3. Dispense 50 µL of ¹²⁵Iodine labeled anti-hGH into each tube, including the uncoated tubes for total counts.
4. Shake the tube rack gently by hand.
5. Incubate for 2 hours at room temperature.
6. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
7. Wash tubes with 3 mL Working Wash solution (except total counts).
Avoid foaming during the addition of the Working Wash solution.
8. Aspirate (or decant) the content of each tube (except total counts).
9. Wash tubes again with 3 mL Wash solution (except total counts) and aspirate (or decant).
10. After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
11. Count tubes in a gamma counter for 60 seconds.

10 CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. On semilogarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of hGH (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.

11 TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

hGH-IRMA	cpm	B/T (%)
Total count	126065	
Standard 0 µIU/mL	241	0.19
1.0 µIU/mL	1254	0.99
2.4 µIU/mL	2743	2.18
8.0 µIU/mL	7817	6.20
25.5 µIU/mL	19515	15.48
59.1 µIU/mL	34310	27.22
119.9 µIU/mL	49965	39.63

12 PERFORMANCE AND LIMITATIONS

12.1 Detection limit

Twenty zero calibrators were assayed along with a set of other calibrators.

The detection limit, defined as the apparent concentration two standard deviations above the average counts at zero binding, was 0.04 µIU/mL.

12.2 Specificity

Cross-reactive hormones were added to a low and to a high hGH value calibrator. The apparent hGH response was measured.

added Hormone	hGH CAL 1 µIU/mL	hGH CAL 5 µIU/mL
-	0.60	65
hCG 300000 mIU/mL	0.63	63
hPL 50000 ng/mL	0.69	68
PRL 2000 ng/mL	1.35	63

Note : this table shows the cross-reactivity for the hGH-IRMA.

12.3 Precision

INTRA-ASSAY				INTER-ASSAY			
Serum	N	$\langle X \rangle \pm SD$ (µIU/mL)	CV (%)	Serum	N	$\langle X \rangle \pm SD$ (µIU/mL)	CV (%)
A	6	6.7 ± 0.3	4.4%	A	20	3.8 ± 0.3	8.1%
B	6	19.0 ± 0.2	1.0%	B	20	17.9 ± 1.2	6.7%

SD : Standard Deviation; CV: Coefficient of variation

12.4 Accuracy

DILUTION TEST

Sample	Dilution	Theoretical Concent. (µIU/mL)	Measured Concent. (µIU/mL)
Serum 1	1/1	-	35.6
	1/2	17.8	17.4
	1/4	8.9	8.8
	1/8	4.5	5.0
	1/16	2.2	2.4
	1/32	1.1	1.2

Sample was diluted with the zero calibrator.

Conversion factor: 1 µIU hGH-IRMA Calibrator = 0.33 ng

RECOVERY TEST

Sample	added hGH (μ IU/mL)	Recovered hGH (μ IU/mL)	Recovered (%)
Serum 1	5	5.1	102
	10	9.5	95
	25	24.7	99
	50	53.7	107
Serum 2	60	63	105
	40	46.8	117
	20	22.1	111
	10	10.2	102
	5	4.8	96

12.5 Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to coated tubes.

TIME DELAY

Serum (μ IU/mL)	0'	15'	30'
Serum 1	4.4	4.5	4.8
Serum 2	18.8	18.1	19.3
Serum 3	27.6	26.2	27.0
Serum 4	36.5	36.8	36.6

12.6 Hook effect

A sample spiked with hGH up to 4000 μ IU/mL gives a signal above the highest calibrator concentration.

13 LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays.
Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies.
If results are not consistent with other clinical observations, additional information should be required before diagnosis.

14 INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practices.

15 REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

In normal subjects, growth hormone (hGH) secretion is pulsatile. During day-time, hGH concentrations range from <0.2 - 10 μ IU/mL. During sleep, hGH concentrations increase consistently (\pm 30 μ IU/mL).

hGH secretion is greatly stimulated by exercise and stress (venous puncture, hypoglycemia, ...), but is decreased by hyperglycemia.

In normal subjects (n=34), two hours after oral glucose load (75 g in adults), hGH levels were lower than 10 μ IU/mL and the hGH response to stimulation tests (insulin, arginine, glucagon administration) exceeded 20 μ IU/mL.

hGH levels were elevated (even after glucose load) in acromegaly (> 10 μ IU/mL).

In hGH deficiency, the response to stimulation test is absent or blunted (short stature by hGH deficiency; hypopituitarism from various origins).

16 PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ^{125}I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, or serum specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

All radioactive handling should be executed in a designated area, away from regular passage. A log book for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiosafety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

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17 BIBLIOGRAPHY

1. CAMANNI F., MASSARA F., BELFORTE L., ROSATELLO A., MOLENATTI G.M. (1977)
Effect of dopamine on plasma growth hormone and prolactin levels in normal and acromegalic subjects.
J. Clin. Endo. and Metab., 44:465.
2. DAUGHADAY W.H. (1981)
The adenohypophysis
In : "Textbook of Endocrinology", R.H. Williams ed., W.B. Saunders, Philadelphia, p;73.
3. DROBNY E.C., AMBURN K., BAUMANN G. (1983)
Circadian variation of basal plasma growth hormone in man.
J. Clin. Endo. and Metab., 57:524-528.
4. HASHIDA S., NAKAWAGA K., ISHIKAWA E., OHTAKI S. (1985)
Basal level of human growth hormone (hGH) in normal serum.
Clin. Chim. Acta, 151:185-186.
5. LEWIS V.J. et al. (1980)
Human growth hormone : A complex of proteins.
Rec. Progr. Horm. Res., 36:488.
6. LI C.H. (1974)
Human growth hormone : perspective on its chemistry and physiology.
In : "Advances in human growth hormone research", S. Raiti ed. U.S. Department of Health, Education and Welfare
Publ. DMEW publ. no. (NJH) 74-612, p. 321.
7. REUTENS AT. (1995)
**Evaluation and application of a highly sensitive assay for serum growth hormone (GH) in the study of adult
GH deficiency.**
J. Clin. Endocrinol. Metab., 80(2): 480-485
8. CROWLEY S. et al. (1995)
Growth and the growth hormone axis in prepubertal children with asthma.
J. Pediatr. 126(2): 297-303.
9. KELLY JT. et al. (1993)
**Effect of different oral oestrogen formulations on insulin-like growth factor I, growth hormone and growth
hormone binding protein in post-menopausal women.**
Clin. Endocrinol. 39(5) : 561-567

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	In vitro diagnostic medical device	In-vitro-Diagnostikum	Dispositivo medico-diagnóstico in vitro	Producto sanitario para diagnóstico In vitro	Dispositif médical de diagnostic in vitro
	Catalogue number	Katalognummer	Numero di Catalogo	Número de catálogo	Référence de catalogue
	Batch code	Chargenbezeichnung	Codice del lotto	Código de lote	Numéro de lot
	Contains sufficient for <n> tests	Ausreichend für <n> Prüfungen	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit	Temperaturgrenzwerte	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date	Verwendbar bis	Utilizzare prima del	Establa hasta	Utiliser jusque
	Manufacturer	Hersteller	Fabbricante	Fabricante	Fabricant
	Distributor *	Vertriebspartner *	Distributore	Distribuidor	Distributeur
	Date of manufacture	Herstellungsdatum	Data di produzione	Fecha de fabricación	Date de production
	Biological risks	Biologische Risiken	Rischi biologici	Riesgos biológicos	Risques biologiques
	Caution	Achtung	Attenzione	Precaución	Attention
	Unique device Identifier	eindeutige Produktidentifizierung			
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité