## Instructions for Use

# 11-Desoxycortisol (CT) RIA

IVD

CE

REF RIA-4777

96

DRG Instruments GmbH, Germany  
Frauenbergstraße 18, 35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

DRG International, Inc., USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.**

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## 1 INTENDED USE

The 11-Desoxycortisol (Compound S) is an intermediate steroid in the glucocorticoids biosynthesis. Precursor of the cortisol, it comes from 17-hydroxyprogesterone after action of the 21-hydroxylase. This parameter is interesting for the diagnosis and follow up of treatment in case of surrenal enzymatic deficiency in 11 beta hydroxylase which is responsible for congenital surrenal hyperplasia in children and hyperandrogenics in adult women

Under hypothalamic-pituitary control via the Adrenocorticotropic hormone (ACTH), the secretion of 11-Desoxycortisol follows a nycthemeral variation: it reaches a maximum in the morning (around 8 a.m.) and a minimal during the night (between 0 to 4 a.m.)


The measure of the 11-Desoxycortisol can be performed in serum using immunological competition methods (RIA).

The metyrapone inhibits the 11 beta hydroxylase and the conversion of the 11-Desoxycortisol in cortisol.

The metyrapone (Metopirone<sup>®</sup>) test is an indicator of the ACTH reserve.

## 2 PRINCIPLE OF THE METHOD

The 11-Desoxycortisol CT RIA obeys the law of mass action according to the following equation:

Since the concentrations of  $^{125}\text{I}$  - 11 DOC and coated antibodies are constant, the advancing state of the equation depends on the concentration of 11 DOC. The amount of  $^{125}\text{I}$  - 11 DOC bound to the coated tube is inversely proportional to the concentration of 11 DOC in the sample.

Following the incubation, the tube is aspirated to remove excess unbound  $^{125}\text{I}$  labelled 11-DOC.

Patient sample concentrations are read from a calibrator curve.

## 3 MATERIAL PROVIDED AND STORAGE

**Stored at 2 °C - 8 °C, the material can be used up to the expiration date printed on each label.**

### 3.1. Coated tubes

2 x 48 **Polystyrene tubes** (12 x 75 mm) coated with anti-11 DOC polyclonal antibodies.

Systematically, allow the coated tubes to reach room temperature (18 °C - 25 °C) before use.

Store at 2 °C - 8 °C.

### 3.2. (Ag $^{125}\text{I}$ ) $^{125}\text{I}$ - 11 DOC E analog tracer (yellow, 52 mL)

1 bottle, containing < 5 µCi (105 kBq) of  $^{125}\text{I}$ -labelled 11 DOC analog in protein-based buffer containing < 0.1% NaN<sub>3</sub> as a preservative. Store at 2 °C - 8 °C.

### 3.3. (CAL N) 11 DOC Calibrators

5 vials, 0.5 mL each, containing 11 DOC calibrators in human serum with < 0.1% sodium azide as a preservative.  
**(See exact values on vial labels.)**

Store at 2 °C - 8 °C. Can be used up to the expiration date printed on each label.

### 3.4. (CAL 0) 11-DOC Zero Calibrator

One vial, 1 mL, containing 0 ng/mL 11-DOC zero calibrator in human serum with < 0.1% sodium azide as a preservative.

Store at 2 °C - 8 °C. Can be used up to the expiration date printed on the label.

### 3.5. (CONTROL N) Control 1 – 2

2 vials, lyophilised, 0.5 mL each, level I and level II containing low and high 11 DOC levels in human plasma with < 0.1% sodium azide as a preservative. **(See exact values on vial labels.)**

Store at 2 °C - 8 °C. Can be used up to the expiration date printed on each label.

*Before use, reconstitute the contents of the controls with 0.5 mL of distilled water.*

### 3.6. (WASH SOLN CONC) Concentrated Wash Buffer

1 bottle concentrated buffered solution containing sodium azide (NaN<sub>3</sub> < 0.1 %).

*Poor the solution in 700 mL of distilled water.*

The diluted washing solution is stable for 2 weeks at 2 °C - 8 °C if covered with adhesive film to avoid contamination.

**Note: Conversion factor:** 1 ng/mL = 2.886 nmol/L  
1 nmol/L = 0.3465 ng/mL

#### 4 MATERIAL REQUIRED BUT NOT PROVIDED

- bench surfaces protected by absorbent paper to reduce the effects of radioactive spillage.
- waste disposal containers appropriately labelled and designed as suitable for solid or liquid radioactive materials.
- manual or automated precision micropipettes for dispensing samples or reagents without cross-contamination.
- absorbent paper.
- vacuum pump connected through a trap for aspiration.
- water bath
- a gamma scintillation counter.

#### 5 METHODOLOGY

##### 5.1 Collection and handling of blood samples

The blood sample may be collected into a dry tube.

After separation from the red blood cells, serum samples may be assayed immediately, within 24 hours if stored at 2 °C - 8 °C, or later, after period up to several months if stored at -20 °C.

Repeating freezing and thawing must be avoided.

For samples from children of less than 2 years old, the 11-Desoxycortisol must be assayed immediately after an extraction step.

##### 5.2 Assay Procedure

Reagents stored at 2 °C - 8 °C must be brought at room temperature (18 °C - 25 °C) prior to use.

Do not mix reagents of different lots.

Label the tubes for T (« Total Counts » do not use coated tubes) calibrators, samples and control sera.

Perform the assay in duplicate. Calibrators, controls and samples must be assayed at the same time.

Each tube can only be used once.

###### 1. Calibrator curve:

Pipette 25 µL of each calibrator into the corresponding tubes.

###### 2. Unknowns and control sera:

Pipette 25 µL of each sample or control sera into the corresponding tubes.

###### 3. Add 500 µL of $^{125}\text{I}$ – 11 DOC tracer to each tube.

###### 4. Vortex, cover and incubate 2.5 hours at 37 °C ± 2 °C.

###### 5. Carefully aspirate or decant the solution of all tubes (except total counts tubes).

###### 6. Add 2 mL of washing solution to each tube. Aspirate or decant carefully (except total counts tubes).

###### 7. Add 2 mL of washing solution to each tube. Aspirate or decant carefully (except total counts tubes).

###### 8. Count the radioactivity fixed in each tube for at least 60 seconds.

##### 5.3 Data processing

Determine the mean count rate for each set of duplicate tubes.

Calculate the ratio B/B0 as follows:

$$\text{B/B0 \%} = [\text{CSL or Sample cpm} / \overline{\text{B0 (Std 0) cpm}}] \times 100$$

Draw the calibrator curve by plotting the ratio B/B0 % (linear scale) obtained for each calibrator versus its respective concentration expressed in ng/mL.

11-DOC concentrations in samples may be read directly from the calibrator curve.

If a computer is used to calculate the results, the data can be fitted to the appropriate equation, 4-parameter logistic function curve fitting is the preferred method

#### 5.4 Example of a typical assay:

	Content (ng/mL)	cpm 1 <sup>st</sup> duplicate	cpm 2 <sup>nd</sup> duplicate	Mean count rate	B/Bo (%)	11-DOC (ng/mL)
Total counts	-	63015	62510	62762	-	-
CAL 0	0	31865	31445	31655	100	
CAL 1	0.3	27429	28004	27711	87.6	
CAL 2	1.5	20079	19803	19941	62.9	
CAL 3	5	13246	13152	13199	41.7	
CAL 4	20	6369	6509	6439	20.3	
CAL 5	65	3464	3480	3472	11	
C <sub>1</sub> Low	-	22977	22569	22773	72	1
C <sub>2</sub> High	-	13711	13699	13705	43.3	4.5
Sample	-	4172	4242	4214	13.3	45.2

**Example of a typical assay, do not use for calculations**

## 6 PERFORMANCE CHARACTERISTICS

### 6.1 Specificity

Steroid	% Cross-reactivity
11-desoxycortisol	100
17 $\alpha$ Hydroxyprogesterone	5.6
Desoxycorticosterone	0.46
Progesterone	0.59
Cortisol	0.09
Estradiol-17 $\beta$	0.03
Androstenedione, Corticosterone, Dexamethasone, Cortisone, Testosterone, DHEA-S, Aldosterone	N.D.

### 6.2 Minimum detectable concentration of 11-DOC

The minimum detectable concentration has been assayed at 0.11 ng/mL and corresponds to the concentration given by two standard deviations below the mean cpm of 20 replicate determinations of the zero calibrators.

### 6.3 Recovery test

When sera of known 11-DOC contents have their 11-DOC supplemented by addition of 11-DOC in equal volumes (1/1), a satisfactory correlation between theoretical and assayed 11-DOC is obtained.

Added 11-DOC (ng/mL) (1:1 in serum sample)	1.5	5	20	65
Theoretical (ng/mL)	0.85	2.6	10.1	32.6
Assayed 11-DOC (ng/mL)	0.7	2.5	10.2	30.2
% recovery	82.4	96.2	101	92.6

### 6.4 Dilution test

The dilution test indicates that there is immunological identity between the 11-DOC present in the sample and the 11-DOC used to calibrate the test.

Dilution Factor	1	1/2	1/4	1/8	1/16	1/32	1/64
Assayed 11-DOC (ng/mL)	40.1	18.4	9.8	4.7	2.5	1.1	0.6
Expected 11-DOC (ng/mL)	-	20.05	10	5	2.5	1.25	0.63
% Recovery	-	92	98	94	100	88	96

## 6.5 Reproducibility

	Within assay variation 5 replicates		Between assay variation 8 replicates	
	Mean (ng/mL)	% CV	Mean (ng/mL)	% CV
Pool 1	1.1	6.2	1.95	8.7
Pool 2	3.7	5.2	5.48	11.5
Pool 3	28.3	7.7	36.85	15.1

## 7 LIMITATION OF THE PROCEDURE

- The results obtained from this or any other diagnostic kit should be used and interpreted only in the context of an overall clinical picture.
- Do not use lipemic, haemolyzed, icteric or turbid specimens.

## 8 EXPECTED VALUES

These values are given only for guidance; each laboratory should establish its own normal range of values.

The concentration range (2.5 to 97.5% percentile) obtained on 714 patients without stimulation:  
0.43 to 7.56 ng/mL (1.24 to 21.82 nmol/L).

According to literature, ICE (UK) and J.endocrinol.Invest.2003; 26:25-30, after stimulation with Metopirone®:  
> 200 nmol/L (> 69.3 ng/mL).

## 9 WARNING AND PRECAUTION

### For in vitro diagnostic use only

#### Safety

This kit contains  $^{125}\text{I}$  (half-life: 60 days), emitting ionizing X (28 keV) and  $\gamma$  (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with safety regulations.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

For more information, see Safety Data Sheet (SDS).

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