

DxDATA™

Automated Chemiluminescent
Immunoassay Analyzing System

NEW



FEATURES

Small-size: POCT

Accurate: Magnetic particle chemiluminescent immunoassay method

Convenient: Single cartridge, low maintenance, no liquid path, easy to operate, easy to maintain

TEST MENU

TBI (Traumatic Brain Injury)

GFAP
PGP 9.5

Inflammation

IL-6
PCT
SAA
CRP

Covid-19 Products:

SARS-CoV-2 IgM
SARS-CoV-2 IgG
SARS-CoV-2 NAb

Cardiovascular & Cerebrovascular Disease

CKMB
cTnI
D-Dimer
hsCRP
Myoglobin
NT-ProBNP

Thyroid Function

TSH
T3*
fT3*
T4*
fT4*
Anti-TG
Anti-TPO

SPECIFICATIONS

Size: 610 X 410 X 520 mm

Weight: 35 kg

Sample throughput: 8 channels

Sample type: Serum, Plasma, Whole blood

Detection time: 17-29 minutes

LIS/HIS: Unidirectional

Diabetes

C-Peptide*
HbA1c
Insulin

Gonadal Hormone

AMH
Estradiol*
β-hCG
FSH*
LH*
Progesterone*
Prolactin*
Testosterone

Maternal & Child Care

Folic Acid
Vitamin B12
Ferritin

Bone Metabolism

PTH
25-OH Vitamin D

Tumor marker

AFP
CEA
CA-125
CA 19-9